Fotoplast® S IO

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Produktbeschreibung

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
Directions for use

Product-Description

Fotoplast® S IO is a light-curing one component resin for the production of thin-walled shells for ITE-hearing aids.

For curing the material a light source with a light spectrum between 320 – 400 nm is necessary. The polymerisation times depend on the transparency of the negative forms, the light intensity of the polymerisation unit, and on the desired wall thickness of the shell.

For the shell forming in Dreve-light units (Polylux 1000N; Polylux 2000 and Polylux Carousel III) a polymerisation time of about 20 seconds (Polylux 3000, about 40 seconds) is necessary to achieve a wall thickness of 0.7 mm. A longer polymerisation time will increase the thickness.

After pouring out the material a later polymerisation of 2 minutes (Polylux 3000, 4 minutes) is necessary. The development of an inhibition layer inside the shells can be avoided either by filling the shell with Glysol or by polymerisation under protection gas atmosphere.

Processing

- Pour Fotoplast® S IO in a prepared form for the processing of an ITE shell.
- Cover the form with a cover plate.
- Put the form into a light polymerisation unit (Polylux) and start the polymerisation process.
- After the polymerisation time take the form out of the light-polymerisation unit.
- Remove the cover plate and pour out the liquid material out of the shell. The material can be reused.
- For avoiding the inhibition layer fill the shell with Glysol.
- For later polymerisation put the form with the shell (without cover plate) into the light-polymerisation unit.
- After finishing the polymerisation pour out the Glysol, take the shell out of the form and go on with further processing.

Advice: Shake well before use and let it rest until the material is bubble-free.

Storage: Keep bottle tightly closed and protect from light!

Fotoplast® S IO clear transparent:

Hazard statements

H315 Causes skin irritation.
H319 Causes serious eye irritation.
H317 May cause an allergic skin reaction.
H335 May cause respiratory irritation.

Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves / eye protection / face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: 1,4 Butanediol dimethacrylate, hydroxyl-propylmethacrylate, monoester with propane-1,2-diol, Acrylic resin, Diphenyl(2,4,6-trimethylbenzoyl)-phosphin oxide, tripropylene glycol diacrylate

Fotoplast® S IO colored transparent:

Hazard statements

H413 May cause long lasting harmful effects to aquatic life.

Precautionary statements

P273 Avoid release to the environment.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Additional information:

Contains Diphenyl(2,4,6-trimethylbenzoyl)phosphin oxide. May produce an allergic reaction.

Fotoplast® S IO opaque:

Hazard statements

H315 Causes skin irritation.
H319 Causes serious eye irritation.
H317 May cause an allergic skin reaction.
H335 May cause respiratory irritation.

Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves / eye protection / face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P405 Store locked up.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Contains: ethoxylated bisphenol A dimethacrylate, 2-hydroxyethyl methacrylate, hydroxypropylmethacrylate, monoester with propane-1,2-diol, Diphenyl(2,4,6-trimethylbenzoyl)phosphin oxide

Lot number / Best before date:

The lot number and the best before date are indicated on each Fotoplast® S IO packaging. In case of claims please always indicate the lot number of the product. Do not use the product after expiry of the best before date.

Dreve Otoplastik GmbH is not liable for any damages caused by improper application of the material.

2017-01-19 / Rev. 1
Dreve Otoplastik GmbH
Max-Planck-Str. 31 • 59423 Unna / Germany
Tel.: + 49 2303 8807-0 • Fax: + 49 2303 82909
E-Mail: otoplastik@dreve.de • Internet: www.dreve.com

Danger

CE
0044