

EPA Registration #

- 67619-24 (Liquid)
- 67619-25 (Wipes)

Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant



Preparation/Revision Date
April 2016

Overview

Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant is a bactericidal, virucidal and fungicidal disinfectant for use on hard, non-porous surfaces. The formula has EPA-registered efficacy against over 40 pathogens, with 30 second – 1 minute kill times on some of the most common bacteria and viruses found in healthcare settings. It is available in both ready-to-use liquid and pre-moistened wipes for everyday cleaning and disinfecting. The spray can also be used to sanitize soft surfaces.

Ingredients

Active Ingredient: Hydrogen Peroxide.....	1.4%
Other Ingredients.....	98.6%
Total (not including weight of wipe).....	100.0%

Efficacy

	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Liquids	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes
Disinfection of Hard, Non-Porous Surfaces		
Bacteria contact times ¹	30 sec-1 minute	30 sec-1 minute
Mycobacterium (TB) contact time	4 minutes	5 minutes
Virus contact times ²	30 sec-1 minute	30 sec- 1 minute
Norovirus contact time	1 minute	3 minutes
Antibiotic-resistant organism kill claims	13	13
Total kill claims	43	42
Soft Surface Sanitization		
Bacteria Kill Time	30 seconds	N/A

1. Excluding TB, *Staphylococcus epidermidis*.

2. Excluding Norovirus.

Data based on comparison of EPA federal masters as of 2/2014. Organisms selected based on top 14 HAI-associated pathogens reported to the National Healthcare Safety Network from 2009–2010. 2/1/14.

EPA-Registered Disinfecting Efficacy Against:

KILL TIMES
Liquid Wipes

		Liquid	Wipes
Bacteria	<i>Acinetobacter baumannii</i>	30 sec	1 min
	<i>Campylobacter jejuni</i>	30 sec	30 sec
	<i>Enterobacter aerogenes</i>	30 sec	30 sec
	Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)	30 sec	30 sec
	Multidrug-resistant <i>Enterococcus faecium</i>	30 sec	30 sec
	ESBL-producing <i>Escherichia coli</i>	30 sec	30 sec
	<i>Escherichia coli</i>	30 sec	30 sec
	<i>Escherichia coli</i> O157:H7	30 sec	30 sec
	<i>Escherichia coli</i> (NDM-1)	30 sec	30 sec
	<i>Klebsiella oxytoca</i>	30 sec	30 sec
	Carbapenem-resistant <i>Klebsiella pneumoniae</i>	30 sec	30 sec
	<i>Klebsiella pneumoniae</i> (NDM-1)	30 sec	30 sec
	ESBL-producing <i>Klebsiella pneumoniae</i>	30 sec	30 sec
	Multidrug-resistant <i>Klebsiella pneumoniae</i>	30 sec	30 sec
	<i>Pseudomonas aeruginosa</i>	30 sec	30 sec
	<i>Salmonella enterica</i>	30 sec	30 sec
	<i>Staphylococcus aureus</i>	1 min	1 min
	Linezolid-resistant <i>Staphylococcus aureus</i> (LRSA)	1 min	1 min
	Methicillin-resistant <i>Staphylococcus aureus</i>	1 min	1 min
	Methicillin-resistant <i>Staphylococcus aureus</i> (CA-MRSA USA300)	1 min	1 min
<i>Stenotrophomonas maltophilia</i>	30 sec	30 sec	
Penicillin-resistant <i>Streptococcus pneumoniae</i>	30 sec	30 sec	
<i>Streptococcus pyogenes</i>	30 sec	30 sec	
Viruses	Avian Influenza A Virus (H5N1)	30 sec	1 min
	Cytomegalovirus	30 sec	30 sec
	Hepatitis B Virus (HBV)*	30 sec	30 sec
	Hepatitis C Virus (HCV)	30 sec	30 sec
	Herpes Simplex Virus Type 1 (HSV-1)	30 sec	30 sec
	Herpes Simplex Virus Type 2 (HSV-2)	30 sec	30 sec
	Human Immunodeficiency Virus Type 1 (HIV-1)	30 sec	30 sec
	Human Coronavirus	30 sec	30 sec
	Influenza A Virus (H1N1)	30 sec	30 sec
	Influenza A Virus/Hong Kong	30 sec	30 sec
	Influenza B Virus	30 sec	30 sec
	Norovirus*	1 min	3 min
	Respiratory Syncytial Virus (RSV)	30 sec	30 sec
	Rhinovirus	1 min	1 min
	Rotavirus	1 min	1 min
TB FUNGI	<i>Candida albicans</i>	3 min	5 min
	<i>Trichophyton mentagrophytes</i> (Athlete's Foot Fungus)	3 min	5 min
	<i>Mycobacterium bovis</i> (TB)	4 min	5 min

* Obtained via a surrogate using Feline calicivirus.

Toxicity

HMIS

Health Hazard 1 Flammability 0 Physical Hazard 0 Personal Protection A

EPA Statement

This chemical is a pesticide product registered by the Environmental Protection Agency and is subject to certain labeling requirements under federal pesticide law. These requirements differ from the classification criteria and hazard information required for safety data sheets and for workplace labels of non-pesticide chemicals. Following is the hazard information as required on the pesticide label:

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet.

Directions for Use

Wipes



Hard, non-porous surfaces: Remove gross soil, if visible.* WIPE surface with towel until completely wet. Gloves should be worn. To DISINFECT, let stand for the required contact time as listed on the product back label.

Pull-Top Liquid



Hard, non-porous surfaces: Remove gross soil, if visible.* Squeeze the bottle to distribute the liquid evenly on the surface. To DISINFECT, let stand for the required contact time as listed on the product back label. WIPE with a clean, damp cloth. Allow to air dry.

Spray



Hard, non-porous surfaces: Remove gross soil, if visible.* Spray 6-8" from surface until surface is completely wet. To DISINFECT, let stand for the required contact time as listed on the product back label. WIPE with a clean, damp cloth. Allow to air dry.

Soft surface sanitization: Hold bottle 6-8" from surface. Using a smooth zigzag motion, spray mist across fabric until wet. Do not saturate. Allow fabric to air dry. Bacteria will be killed in 30 seconds.

*This is a one-step product. Precleaning is only required if blood, vomit or fecal matter or other soil is visible.



Visit us: www.cloroxhealthcare.com

Call: 1-800-234-7700

Email: healthcare@clorox.com

© 2016 Clorox Professional Products Company, 1221 Broadway, Oakland, CA 94612