

**PDI**

# SANI-CLOTH® AF<sup>3</sup>

## GERMICIDAL DISPOSABLE WIPE

Bactericidal, Tuberculocidal, Virucidal\*\* in 3 minutes  
Fragrance Free

REORDER NO. P13872

**Bacteria**

■ *Bordetella bronchiseptica* ■ *Bordetella pertussis* ■ *Burkholderia cepacia* ■ *Campylobacter jejuni* ■ *Enterobacter aerogenes* ■ *Escherichia coli* O157:H7 ■ *Escherichia coli* ■ *Klebsiella pneumoniae* ■ *Listeria monocytogenes* ■ *Proteus vulgaris* ■ *Pseudomonas aeruginosa* ■ *Salmonella enterica* ■ *Serratia marcescens* ■ *Shigella dysenteriae* ■ *Staphylococcus aureus* ■ *Streptococcus pyogenes* ■ *Vibrio cholera* ■ *Yersinia enterocolitica*

**Multi-Drug Resistant Bacteria**

■ *Acinetobacter baumannii*, Multi-Drug Resistant ■ Community Acquired Methicillin Resistant *Staphylococcus aureus* (CA-MRSA) NARSA NRS123 (Genotype USA 400) ■ Community Acquired Methicillin Resistant *Staphylococcus aureus* (CA-MRSA) NARSA NRS384 (Genotype USA 300) ■ ESBL Resistant *Escherichia coli* ■ *Escherichia coli* - NDM-1 Positive CDC 1001728 ■ ESBL Resistant *Klebsiella pneumoniae* ■ *Klebsiella pneumoniae* - Carbapenem Resistant ■ *Klebsiella pneumoniae* - NDM-1 Positive CDC 1000527 ■ Methicillin Resistant *Staphylococcus aureus* (MRSA) ■ *Streptococcus pneumoniae* - Penicillin Resistant ■ Vancomycin Resistant *Staphylococcus aureus* (VRSA) NARSA VRS1 ■ Vancomycin Resistant *Enterococcus faecalis* (VRE)

**Viruses\*\***

■ \*\*Avian Influenza A H5N1 virus ■ \*\*Cytomegalovirus, Strain AD-169 ■ \*\*Herpes simplex virus, type 2, Strain G ■ \*\*Human Coronavirus, Strain 229E ■ \*\*Influenza A virus/Hong Kong Strain ■ \*\*Influenza B virus, Strain B/Hong Kong /5/72 ■ \*\*Respiratory syncytial virus (RSV), Strain Long ■ \*\*Rotavirus Strain WA ■ Kills Pandemic 2009 H1N1 Influenza A Virus

**Bloodborne Pathogens**

■ \*\*Hepatitis B virus (HBV) - Duck HBV<sup>†</sup> ■ \*\*Hepatitis C virus (Human)(HCV) - Bovine Diarrhea Virus<sup>†</sup> ■ \*\*HIV-1 (AIDS virus)<sup>†</sup> Strain HTLV-IIIg

**Pathogenic Fungi**

■ *Candida albicans*

**TB**

■ *Mycobacterium bovis* - BCG (TB)

**ACTIVE INGREDIENTS:**

n-Alkyl (68% C<sub>12</sub>, 32% C<sub>14</sub>) dimethyl ethylbenzyl ammonium chlorides ..... 0.14%

n-Alkyl (60% C<sub>14</sub>, 30% C<sub>16</sub>, 5% C<sub>12</sub>, 5% C<sub>18</sub>) dimethyl benzyl ammonium chlorides ..... 0.14%

**OTHER INGREDIENTS** ..... 99.72%

**TOTAL** ..... 100.00%

Does not include the weight of the wipe.

**Net Contents**

**160** Large Wipes 6 x 6.75 in (15.2 x 17.1 cm)  
Net Wt. 1.0 lb. 11.4 oz. (777 g)

KEEP OUT OF REACH OF CHILDREN

**CAUTION** See back panel for additional precautionary statements

80PT3800

\*KILLS HIV-1 (AIDS VIRUS), HEPATITIS B VIRUS (HBV) AND HEPATITIS C VIRUS (HCV) ON PRE-CLEANED ENVIRONMENTAL SURFACES, OBJECTS AND PRE-CLEANED EXTERNAL SURFACES OF ULTRASOUND TRANSDUCERS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in three (3) minutes at room temperature (68°-77°F) in healthcare or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV) and Hepatitis C Virus (HCV).

**Use sites:** Hospital and Healthcare Settings: Ambulances, Ambulatory Surgical Centers (ASC), Dental Offices, Dialysis Clinics, Healthcare Settings, Home Health Care, Hospitals, Laboratory, Nursing homes, Physician's offices and Rehabilitation.

**Critical Care Areas:** CCU, Emergency Rooms, E.R., ICU, Neonatal Intensive Care Units (NICU), Operating Rooms, Pediatric Intensive Care Units (PICU) and Surgery.

**Hospital, Healthcare, and Critical Care Use Sites:** May be used on hard nonporous surfaces of: bed railings; blood glucose meters; cabinets; carts; chairs; counters; dental unit instrument trays; exam tables; gurneys; infant incubators; infant warmers; isolettes; IV poles; operating room tables and lights; physical therapy (PT) equipment surfaces; stethoscopes; stretchers; tables; telephones; toilet seats; and ultrasound transducers and probes.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection.

**To Dispense Wipes:** Remove lid and discard inner seal from container. Find center of wipe roll, remove first wipe for use, twist corner of next wipe into a point and thread through the hole in the container lid. Pull through about one inch. Replace lid. Dispense remaining wipes as necessary by pulling out at an angle. When not in use keep center cap of lid closed to prevent moisture loss.

**DIRECTIONS FOR USE**

It is a violation of Federal law to use this product in a manner inconsistent with its labeling.

**CLEANING AND DEODORIZING:** Wipe surface with towelette until clean. Let air dry. For disinfecting hard non-porous surfaces.

**TO DISINFECT:** Use a wipe to remove heavy soil. Unfold a clean wipe and thoroughly wet surface. Treated surface must remain visibly wet for a full three (3) minutes. Use additional wipe(s) if needed to assure continuous 3 minute wet contact time.

**ULTRASOUND TRANSDUCERS:** For use on external surfaces of Ultrasound Transducers only

To disinfect: Pre-clean ultrasound gel from external transducer to be disinfected by using a clean wash cloth. Then thoroughly wet the probe with a towelette. Treated surface must remain visibly wet for a full three (3) minutes. Use additional wipe(s) if needed to assure continuous 3 minute wet contact time. Let air dry.

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HEPATITIS B VIRUS AND HEPATITIS C VIRUS OF SURFACES/OBJECTS SOILED WITH BLOOD/BODY FLUIDS**

**Personal protection:** When using this product, wear disposable protective gloves, protective gowns, face masks, or eye coverings as appropriate when handling HIV-1 (AIDS Virus), HBV or HCV infected blood or body fluids.

**Cleaning procedure:** All blood and other body fluids must be thoroughly cleaned from surfaces and objects before disinfection by the germicidal cloth. Open, and unfold first germicidal cloth to remove heavy soil.

**Disposal of infectious materials:** Dispose of used towelette in accordance with local regulations for infectious waste disposal.

**Contact time:** Use second germicidal cloth to thoroughly wet surface. Allow to remain wet three (3) minutes, let air dry. Although efficacy at a 30 second contact time has been shown to be adequate against HIV-1 (AIDS VIRUS), this time is not sufficient for the other non-HIV organisms listed on this label. Therefore, a 3 minute wet contact time must be used for all listed non-HIV organisms.

**PRECAUTIONARY STATEMENTS:****Hazards to humans and domestic animals**

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes and clothing. Wash hands thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using toilet.

**FIRST AID**

**If in eye:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing eye. Call a poison control center or doctor for treatment advice. Have the product container or label with you when calling a poison control center or doctor or going for treatment.

**STORAGE AND DISPOSAL**

Do not contaminate water, food, or feed by storage and disposal.

**Storage:** Store in dry location. When not in use keep center cap of lid closed to prevent moisture loss. **Towelette**

**Disposal:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet. **Container Disposal:** Nonrefillable container. Do not reuse or refill this container. Triple rinse and offer for recycling. If recycling is not available, put in trash collection.

EPA REG. NO. 9480-9

EPA EST. NO. A= 9480-NY-1, B= 72956-WI-1

C= 72956-AR-1

D=8251-OH-1

Alpha character will precede batch

code on product

MADE IN USA

NOT A SKIN WIPE

FOR USE ON HARD NONPOROUS SURFACES ONLY



(01)10310819001066

8AP13801

Manufactured for:  
Professional Disposables International, Inc.  
Two Nice-Pak Park • Orangeburg, NY 10962-1376 • USA  
1-800-999-6423