Audiometer

# QUICK USER'S HANDBOOK

# **Release history**

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this page is intentionally blank

## Table of contents

1.Introduction and intended use	63
2.Description	64
2.1 Precaution	
2.2 Warranty and maintenance	
2.3 Unpacking and checking the R17A	
2.4 Accessories	
25 Connecting the accessories	71
2.6 Power Sources	
2.7 Installation	72
2.8 Battery indicator	
2.9 Symbols present on device	73
3. Components and connectors	74
3.1 Power on and off the device	
3.2 PIN Code request	75
4. Performing tests	76
4.1 Pure Tone Audiometry (AC)	
4.1.1 Ambient noise monitoring	
4.1.2 External microphone calibration	
4.2 Speech Audiometry	
4.3 Autothreshold Audiometry	
4.4 Operating position of the device	
4.5 RE-Play Audiometry	
4.6 Weber	
4.7 Otoscopy	
5. Setup settings	
5.1 General Setup	
5.2 Audiometry Setup	
5.3 Audiometry Calibration	
6. Patient database management	
7. Print management	
8. Data transmission to PC	
9. Main Labels	
10. Messages and warnings	
12. Cafe and maintenance	
12. Safety regulations	
12.1 Electrical Safety	
12.2 Device nanaling and operation	
12.3 Medsurement Sarety	
13. Technical specifications	
14. Electromagnetic compatibility	

## 1.Introduction and intended use

RESONANCE thanks you for having purchased one of its medical devices.

This manual contains information and warnings aimed at ensuring the correct and safe use of the R17A Audiometer available in the following versions:

Description	Product code
R17A AC Screening	4300104000
R17A BC Basic Diagnostic	4300104010
	4300104050
R17A HDA Basic Diagnostic	4300104020

Carefully store this manual away in a safe place for each future consultation.

The R17A Audiometer is designed and manufactured to meet or exceed all quality and safety requirements, and has been certified with CE-symbol according to European Medical De-vices Directive. All the functions are software-controlled, upgrading the software at a later date will be simple and free of charge.

The R117A Audiometer is intended to be used by audiologists, hearing healthcare professionals, or trained technicians.

The audiometer is a medical device used to assess different forms of deafness and diseases affecting the auditory system. It generates stimulus and masking signals with the possibility to vary the pressure level within a wide and dynamic value range. This test requires the active patient's participation, who must provide an answer to the various input signals, presented using headphones or other transducers.

The user must therefore have an adequate knowledge of audiological diagnostic equipment, expertise on the examination methods used and knowledge of English.

This device is intended to be used in medical environments with low levels of background noise. The operator must exercise caution in interpreting results where background noise has not been sufficiently isolated.

It is recommended that the device be operated within an ambient temperature range of 15 -35 degree Celsius.

This operating manual is not intended as a training manual for audiometry and tympanometry. The reader should consult standard audiology texts for the theory and application of the screening tests provided by this device.

## 2.Description

### 2.1 Precaution

**WARNING** indicates a hazardous situation which, if not avoided, could result in death or serious injury.

**NOTICE** is used to address practices not related to personal injury.

This instrument has been built and tested in compliance with product standards and left the production factory in perfect technical safety conditions.

- Prior to testing, the performance of the unit should be checked by the operator on their own ear. This is to ensure that the unit produces results expected results.
- The device transducers and parts in contact with the patient should be cleaned following every use. Appropriate antiseptic solution should be used in compliance to the maintenance and cleaning section of this manual.
- The device must be used only with the transducers provided by the manufacturer. Use of alternative transducers, or their replacement, will require re-calibration.
- The device is intended to be powered from Li-ion internal rechargeable battery or from a medical CE power supply. Use of power supplies other than those provided by the manufacturer is not permitted. When used, position so that the power plug can be easily disconnected.
- Although the instrument fulfils the relevant EMC requirements, precautions should be taken to avoid unnecessary exposure to strong electro-magnetic fields, e.g. from mobile phones, microwaves, radiotherapy devices, etc.; these may induce malfunctions in their performance. If the device is used adjacent to other equipment it must be observed that no mutual disturbance appears.
- R17A does not have specific ingress protection. Care should be taken to avoid liquid ingress. Immersion of the device may result in serious damage.

Any possible modification and/or alteration to this device is prohibited, without prior authorization thereto by the Manufacturer.

60601-1.

If the R17A is connected to other devices with their own power source, the system conforms to IEC 60601-1 must be checked under the responsibility of the installer; in the event that isolation devices are used to achieve such conformity, these must comply with IEC 60601-1. Any person who connects external equipment to the device connectors has formed a system and is therefore responsible for the system to comply with the requirements of IEC

If in doubt, contact your RESONANCE technician or local distributor for help.

Device not to be used in areas in which flammable anesthetics or flammable agents for disinfection or skin cleaning are applied.

Electric and electronic waste may contain hazardous substances and therefore must be disposed of separately. Such products will be marked with the crossed-out wheeled bin shown. The cooperation of the user is important in order to ensure a high level of reuse and recycling of electric and electronic waste. Failure to recycle such waste products in an appropriate way may endanger the environment and consequently the health of human beings.

R 17A User's Manual – Page 64

## 2.2 Warranty and maintenance

**RESONANCE** warrants that:

- The R17A is free from defects in material and workmanship under normal use and service for a period of 36 months from the date of device's installation by authorized personnel to the first purchaser.
- Accessories are free from defects in material and workmanship under normal use and service for a period of 12 months from the date of device's installation by authorized personnel to the first purchaser.

Furthermore, this warranty shall not apply to, and M.R.S. shall not be responsible for, any loss arising in connection with the purchase or use of any RESONANCE device that has been:

- repaired by anyone other than an authorized RESONANCE service representative;
- altered in any way so as, in RESONANCE judgment, to affect its stability or reliability;
- subject to misuse or negligence or accident, or which has had the serial number altered, effaced or removed; or improperly maintained; or used in any manner other than in accordance with the instructions furnished by RESONANCE.
- This warranty is in lieu of all other warranties, express or implied, and of all other obligations or liabilities of M.R.S., and M.R.S. does not give or grant, directly or indirectly, the authority to any representative or other person to assume on behalf of M.R.S. any other liability in connection with the sale of RESONANCE products.
- The R17A may be repaired only by your dealer or by a service company authorized by RESONANCE, recommended by your dealer.
- Only original RESONANCE spare or consumable parts must be used for service or repair the R17A.

In order to ensure that your device works correctly, the R17A Audiometer should be checked and calibrated at least once a year. We recommend to send the R17A to your RESO-NANCE dealer for this service, including the probe and all of the accessories.

The verification of instrument calibration must be conducted by qualified service personnel authorized by the manufacturer via appropriate instrumentation. This must provide the level control and the frequency of all signals emitted by the transducers supplied with the instrument.

Should transducers require replacement, use only parts provided by the manufacturer. Replacement of transducers will require re-calibration of the instrument.

## 2.3 Unpacking and checking the R17A

On checking for packaging and content damage immediately after the device's reception; should any sign of damage or tampering be noted, please notify the carrier immediately.

If the device has been damaged during transportation, the device should be checked for any electrical or mechanical defects.

Should the device be dropped or fall, it may sustain hidden damage that may result in hazardous operation.

In case of any doubt, please do not use the device and either contact your Resonance Dealer, or e-mail the Resonance Customer Support Service directly at <a href="mailto:support@resonance-audiology.com">support@resonance-audiology.com</a>

Keep all original packaging for future use.

The identification label on the rear plate holds the serial number. This should be checked and written down for later service claims.

### 2.4 Accessories

List of accessories DD45 Headset HAD 280 Headset (R17A HAD only) Bone-vibrator B71W/BC-1 alternatively (R17A BC and R17A HAD only) Patient response push button (R17A BC and R17A HAD only) SD-Card with speech material included

**Optional available** IP30 Insert Phones

IP30 Insert Phones External USB keyboard Wi-Fi thermal printer

Otoscope

Free field amplifier and loudspeakers Kit Tips MS clinical type

Medical CE approved power supply Cradle Multilingual CE Operation Manual MDS Software USB pen drive including: DHR documents, MDS software installer, Multilingual CE Operation Manual Carrying bg

#### Available tests Audiometry

Otoscopy (optional test)

\*Autotheshold automated audiometry \*Speech-Audiometry \*Re-Play Audiometry Weber (not available for R17A (AC)

\* = optional test, depending on device model

### R1/A AC Accessories

Image	Description	Standard /Optional	Plug
	DD45 Headset*	٠	ACR (4) ACL
90	Insert Phones IP30 *	o	ACR (4) ACL
9	Patient response push botton*	0	PAT.RESP.(5)
	Medical CE approved power supply: Model: MENB1010A0503 (xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord	٠	5 VDC (1)
ØX	As alternative: Medical CE approved power supply: Model: ME10A0503 (xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord		5 VDC (1)
•••	Charging Cradle •		/
(astarca)	USB pen drive •		For PC use
- searce of	Carrying bag	٠	/
	Wi-Fi Thermal printer	0	/
	External USB Keyboard	0	<b>€ द •</b> (2) ○
	RE-VO Video Otoscope	0	• <b>C</b> • (2)

	Free Field passive loud speakers	0	/
i sag	Free Field amplifier	0	ACR (4) ACL (6)
<ul> <li>= standard</li> </ul>	<ul> <li>o = optional</li> <li>- = not available</li> </ul>	e * = applied par	:

The applied parts are represented by the transducers ad pushbutton which are placed in contact with the patient. They are indicated by an asterisk (\*).

### R17A BC Accessories

Image	ge Description S		Plug	
DD45 Headset*		٠	ACR (4) ACL	
Insert Phones IP30 *		0	ACR (4) ACL	
$\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{\mathbf{$	Bone vibration B71W/BC-1*	•	BC (3)	
9	Patient response push botton*	•	PAT.RESP.(5)	
	Medical CE approved power supply:         Model: MENB1010A0503 (xyz)         (xyz) = B01 changeable direct plug-in         (xyz) = Q01 with power cord         (xyz) = N01 with power cord		5 VDC (1)	
ØX	As alternative: Medical CE approved power supply: Model: ME10A0503 (xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord	٠	5 VDC (1)	
	Charging Cradle	•	/	
( riscource a)	USB pen drive	•	For PC use	

R 17A User's Manual – Page 68

	Carrying bag	٠	/
e la compañía de la compa	Wi-Fi Thermal printer	O	/
	External USB Keyboard	O	<b>•</b> ← <b>(</b> 2) ○
	RE-VO Video Otoscope	0	(2)
	Free Field passive loud speakers	O	/
Ling (P	Free Field amplifier	O	ACR (4) ACL (6)
• = standard	• = optional - = not available	* = applied part	

The applied parts are represented by the transducers ad pushbutton which are placed in contact with the patient. They are indicated by an asterisk (\*).

### R17A HAD Accessories

Image	Description	Standard /Optional	Plug
HDA 280 Headset*		٠	ACR (4) ACL
Insert Phones IP30 *		o	ACR (4) ACL
$\bigcirc$	Bone vibration B71W/BC-1*		BC (3)
9	Patient response push botton*	٠	PAT.RESP.(5)
	Medical CE approved power supply: Model: MENB1010A0503 (xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord	•	5 VDC (1)

ØX	As alternative: Medical CE approved power supply: Model: ME10A0503 (xyz) (xyz) = B01 changeable direct plug-in (xyz) = Q01 with power cord (xyz) = N01 with power cord	•	5 VDC (1)
	Charging Cradle •		/
(riscource t)	USB pen drive	٠	For PC use
	Carrying bag	٠	/
2	Wi-Fi Thermal printer	0	/
	External USB Keyboard	0	● <b>← ᢏ_</b> → (2) ○
	RE-VO Video Otoscope	0	(2)
	Free Field passive loud speakers	0	/
	Free Field amplifier	0	ACR (4) ACL (6)
• = standard	<pre>o = optional -= not available</pre>	* = applied part	

The applied parts are represented by the transducers ad pushbutton which are placed in contact with the patient. They are indicated by an asterisk (\*).

## 25 Connecting the accessories

Reference	Symbol	Function
1	5 VDC	Socket for external power supply
2	•	USB Host port for connection of pen drive (for software update) or external keyboard, Wi-Fi dongle, otoscope
3	BC/Probe	Bone vibrator plug.
4	ACR/Contra	Right transducer plug, right side output for Free Field amplifier connection
5	PAT.RESP.	Patient response push bottom mini usb plug. This connector can be used for wired data transfer
6	ACL	Left transducer plug left side output for Free Field amplifier connection

### 2.6 Power Sources

The device can operate with its own internal rechargeable battery (Li-ion 3.6 V / 3 Ah), or connected to the medical power supply, or placed on the charging cradle.

The average recharge time of an empty battery is less than 4 hours, both with direct feeder link or with its charging cradle. Two hours are enough to reach 80% charge.

To increase the battery duration, the display brightness is reduced to 20% after 60 seconds of inactivity (no action on the touch screen).

After 2 minutes of inactivity starting from the reduction of display brightness, the device switches to sleep mode. Simply touch the screen to reactivate the device. After 8 minutes of inactivity, and if not connected to the battery charger, the instrument switches off automatically. Press the ON / OFF button to restart the device.

## 2.7 Installation

Ensure all transducers are connected prior to starting the device.

Turn the device on and check the battery status: at first switching on, the battery indicator may provide an approximate value. A full charge cycle is suggested prior to its use.

## 2.8 Battery indicator

The battery charge status is showed on the display with a dedicated icon. When less than 10%, the battery indicator flashes to underline that recharging is necessary.

At low charge levels, the device displays a Low Battery warning message. Some tests may not be available following the low-battery warning.

Although the device is able to operate even when the battery is running very low, it is not recommended to use it with less than 5% battery autonomy and it is recommended to connect it to the medical power supply or to its charging cradle.
## 2.9 Symbols present on device

	5 VDC Direct current (Medical CE approved power supply)
<b>(</b>	Refer to instruction
$\triangle$	Refer to manual for transducer and power supply connection
	Class II device (Symbol refers to external medical power supply with double insulation)
*	BF type applied parts
X	Electric and electronic waste may contain hazardous substances and therefore has to be collected separately
	Together with name and address indicates the manufacturer of the instrument
<b>CE</b> 0051	Conformity marks to directive 93/42/CEE – notification body 0051 (IMQ)
2015	Year of production of the device
0	USB Host port
PAT. RESP.	Patient response push button plug

# 3. Components and connectors

	Description
1	Bone vibrator plug
2	Right transducer plug
3	Left transducer plug
4	7" TFT Resistive touch screen display
5	Integrated microphone: used for communication between operator and patient
6	Main switch
7	Led: shows that battery is charging
8	Socket for external medical power supply
9	Patient response push button mini usb plug.
	This connector can be used for wired data transfer
10	USB Host port for connection of pen drive (for software) or external keyboard. Wi-Fi dongle
	and otoscope

### 3.1 Power on and off the device

To turn on the device, press briefly the Main switch (6) and wait for the Pin code request on the display. To switch off the device, press the Main switch (6) for about 2 "and then release the button.

Before powering on the device, ensure that all the accessories are already plugged-in.

### 3.2 PIN Code request

When switching the device on a PIN code is requested to access both the device functionalities and the data it contains. The default PIN code set by the manufacturer is 1234. After the PIN code has been typed in the starting screen is presented. The PIN code can be changed or disabled through the dedicated function available in the setup menu. For additional details see chapter 5.1

In case wrong codes are typed in, after 3 failures the possibility of new attempts is prevented and unlocking is only possible with the PUK code provided along with the device. Should the PUK code be lost, it can be retrieved by communicating to the manufacturer the device serial and product numbers to be found on the label on the lower side of the unit.

## 4. Performing tests

### 4.1 Pure Tone Audiometry (AC)

Pure Tone Audiometry represents the basic test of auditory investigation aiming at evaluating the auditory function and assessing the hearing threshold. The term hearing threshold indicates the minimum level perceived by the test subject.

The test is performed by administering pure tones using the included earphones (Air Conduction – AC). Threshold estimation can be done either for one ear and then for the other, for all tones, or on both sides for each tone.

The test can estimate the AC loss but cannot distinguish between abnormality in the conductive mechanism and sensorineural mechanism.

Before perform the test, eliminate any obstructions which will interfere with the placement of the earphone cushions on the ear (e.g. hair, eyeglasses, rings).

Place the right headphone (identified by red color on the transducer) on the right side and the left (identified by blue color on the transducer) on the left.

Ensure that the cushions are clean to avoid any potential infection's transmission.

Adjust the headband of the headphones sitting on the two metallic bars, to achieve the correct pressure on the ear and the correct alignment of transducer to ear canal.

During BC threshold determination, place the bone vibrator on the mastoid area where the tones are best heard by the patient.

Instruct the patient, before performing the test, to signal pressing the pushbutton (or to raise his/her hand) when he/she hears a sound, and to keep the pushbutton depressed (or his/her hand up) as long as he/she hears it.

Threshold determination: The test begins at 1000 Hz on the patient's better ear.

Present pure tones at 1000, 2000, 4000, 8000 Hz then 500, 250 and 125 Hz\*.

Start presenting the tone at a level loud enough to ensure audibility. The patient will thus be able to identify the type of signal to be expected exactly and give clear responses.

Should levels close to threshold be used when starting the test, the patient may fail to respond as he/she is unsure of what he/she should hear. Should the level prove insufficient to elicit clear responses, increase the level until the patient responds reliably.

To assess threshold level, proceed as follows:

- When the stimulus is clearly perceived by the patient, reduce progressively its level in 10 dB steps.
- When reaching a level, at which the patient fails to respond, further decrease the level by 10 dB and, interrupting the tone manually, begin increasing its level in 5 dB steps, until the patient responds again.
- Decrease this level once more until the patient fails to respond.

The true pure tone test threshold corresponds to the level evoking patient's response in 50% of presentations.

In tinnitus cases, use of pulsed tones (rather than continuous) is preferable.

Avoid any noisy environments because background noise can produce false test results, especially at lower frequencies.

### Performing tone audiometric tests

Turn the device on and, from the main page, press the TONE AUDIOMETRY button (see image on the left) to get to the test. The display will show the parameters set for both channels and the active function button

### **Stimulus configuration**

Press the **STEADY** button to change the type of signal to be emitted: continuous, pulse (2 Hz rate), warble. At every click on the button, the icon changes identifying the type of signal.

### Output configuration

Gently press directly on the audiogram to change the side of the stimulus: right/left. Masking will be automatically sent to the opposite ear.

By shortly pressing the icon for output selection (see image on the left), it is possible to select AC (phones) for the threshold estimation, AC UCL (phones) for uncomfortable level or BC (bone vibrator). By keeping it pressed for a few seconds you can select the free field output. You can go back to the previous selection by pressing the icon again for a few seconds.

The same output connectors are used for both headsets and free field. Before starting the test make sure the appropriate transducer has been connected.

### **Frequency settings**

Press the **ARROW** keys to increase or decrease the frequency of the signals emitted (tone for channel CH1 and narrow band noise for channel CH2). The frequency set is shown between the two arrow keys.

### Intensity settings

Press the ARROW keys to increase or decrease the intensity of the CH1 signal (stimulus) and CH2 (masking). The level reached is indicated between the two arrow keys.

### Signal presentation and storage of threshold levels

Hearing threshold levels can be determined by presenting test signals to the patient with earphones (air conduction – AC). Once the frequency to be analyzed has been set and the intensity of the stimulus defined, press SIGNAL (see image on the left) to generate CH1 tone. The button turns yellow when pressed, to indicate the sending of the signal. The masking on the opposite ear is always active (except when in "Off" condition). The functionality of this button can be modified in the set up.

Determined the threshold at a given frequency, press the button shown to save the point on the audiogram. The button pressed changes its colour to green to confirm the registration. In cases when the output level has been gained without patient's response, it is possible to store the point with the symbol "Under" pressing and holding the button for 1 second.

It is possible to delete any point of the audiogram by moving the cursor on it and pressing and holding the button for 3 seconds. A warning message will ask for confirmation to proceed with deleting the point.

The audiogram area included between the threshold (AC) and the uncomfortable level (UCL) curves is highlighted in colors, to simplify the identification of the patient's hearing capability.

### Communicate to the patient

**TALK OVER** button allows the user to talk to the patient via the transducer. A microphone is integrated in the device.

#### SETUP

Press the setup icon to access the setup page where you can choose useful settings for test execution.

After a point is stored on the audiogram, the frequency automatically switches to the next one; with this key it is possible to define which level to set. 3 selections are possible: keep the previous frequency level, re-set the level to 40dB or reset to the value saved as default.

Pressing this button, you can save the parameters previously set (frequency, stimulus) as default. A dialog window will ask for confirmation before saving the parameters

These keys allow you to define the functionality of the **SIGNAL** key for signal presentation. The first key is used to establish the operating mode, like if the signal is generated when the key is pressed or if it is always present and interrupted when key is pressed. Use the arrow keys to extend the action once the key is released for a period of 0 to 2.5 seconds.

### Exit and Save

Pressing the **EXIT** button, a dialog box will pop up and ask if you wish to save and/or print the test out on the thermal printer.

Press **YES** or **NO.** If answer is YES, R17A will show patient database page to select the patient to be associated to the test performed, confirmation is required pressing this button. Please, refer to section 6 of this user manual for further information about the patient database management.

### Information shown on the display

Indicates the name of the patient undergoing the examination, as previously selected. Touch the icon to get direct access to patient database.

Indicates battery residual charge and current time.

When the patient presses the response pushbutton a green line appears on the display to highlight the answer.

### **Keyboard operation**

The R17A allows for the connection of an external USB keyboard to perform the tone audiometry testing without using the touch screen. The keys and key combinations to be used are the following:

- Stimulus selection: press key I to change stimulus
- Output selection: press **R** for the right side or **L** for the left side (you can also switch sides by pressing **TAB**), press O for the output.
- Frequency selection: left and right arrow keys
- Stimulus intensity: **up** and **down arrow** keys
- Masking level: **shift + up** and **down** arrow keys. As an alternative you can also activate and deactivate the masking function by pressing **M**.
- Stimulus presentation: spacebar. When briefly pressed it activates or deactivates the continuous tone presentation, when kept pressed it behaves like the SIGNAL button on the screen. Press ENTER to save the value.
- Communication with the patient: press **T** to activate and deactivate this function

### 4.1.1 Ambient noise monitoring

Environmental noise monitoring aims to measure in real time the noise present in the environment in which the examination is being carried out, estimating its intensity at different frequencies in order to assess whether the minimum conditions exist to guarantee the accuracy of the test. It is possible to establish minimum noise requirements for each octave band and if these are exceeded, an indication on the screen will inform that there is an unsuitable environment locking the possibility to store the threshold level. Ambient noise is measured using a special calibrated microphone connected to the audiometer with a USB cable, allowing it to be positioned as close as possible to the patient. Environmental noise monitoring is an option of the tonal audiometry test that can be activated by license, see chap. 5.

The environmental noise monitor is displayed in the graph of the non-active side; changing sides will reverse the display. In vertical position the display of ambient noise is always deactivated.

For each octave the acceptable limit in gray color is shown in the background, while the value measured in real time is indicated by green or red bars depending on whether the intensity is lower or higher than the set threshold.

In the local setup menu, there is also the possibility to select different setups, or different curves with different threshold levels to better adapt to the actual test conditions, for example in the case of using a silent booth.

R 17A User's Manual – Page 80

This button permits to enable the ambient noise measure. It is possible to enable or disable the noise measure at any time.

### Threshold levels Setup for different frequency bands

Pressing the button shown it is possible to change the threshold levels for different frequencies. The following setting are available:

**Phones:** use of the phones supplied with the device without any noise reduction system **ADC:** use of the DD45 phones with Audiocups noise reduction enclosures.

1: use of a silent booth with a noise reduction of about 30dB.

2: use of a silent booth with a noise reduction of about 40dB.

Setups 1 and 2 have to be considered as indicative; please refer to section 13 for default values.

It is possible to modify the threshold values for Setups 1 and 2 based on the characteristics of the silent cabin used. The values are contained in 2 text files named Room1\_ and Room\_2 on the USB pen drive supplied with the device. It is sufficient to modify them with a common text editor taking care to maintain the formatting. Once the pen drive has been inserted into the USB of the audiometer, the button shown will appear, pressing it will permanently import the values into the audiometer.

### 4.1.2 External microphone calibration

The first time the environmental noise monitoring function is enabled, the sensitivity and the serial number of the microphone must be set.

A pen drive with calibration data is supplied with the microphone, simply turn on the device and insert the pen drive to import the calibration data.

This operation must be performed only the first time the microphone is connected, or when the license is activated; the microphone is subject to annual calibration, such as headphones or other transducers, do ref. to the chap. 5.6 for more information.

### 4.2 Speech Audiometry

Speech audiometric test has become a fundamental tool in hearing-loss assessment, it also provides information regarding discomfort or tolerance to speech stimuli and information on word recognition abilities. Information gained by speech audiometry can help determine proper gain and maximum output of hearing aids and other amplifying devices for patients with significant hearing losses.

Speech tests can be performed using registered speech material, available in several languages.

Information on patient's speech hearing and understanding can be gained tracing a "Articulation Curve"; after a short exercise at appropriate intensity, a curve that measures the understanding at various levels of intensity can be obtained. Up to 10 values per each curve can be stored.

Speech audiogram: The intensity or loudness of the sound is measured in decibels (dB) on the abscissa, on the ordinate the percentage of discrimination.

**NOTE:** The Articulation Curve varies according to phonemes content and meaning of the speech materiel used.

It is possible to conduct tests using different masking signals as Speech Noise (SN) or White Noise (WN) to define speech-detection threshold. SN focuses its spectrum of energy in those frequencies that are characteristic of the voice, whereas White Noise uses the entire spectrum of frequencies with constant energy.

### **Test execution**

Press the **SPEECH** icon on the main page to access the test. The display will show the parameters set for both channels and the active function buttons.

### Signal settings

Press the **SETUP** icon to access the setup page.

Press the **ARROWS** to change the language of the recorded speech material. The phonemes list will be showed on the display. By touching this area you can choose the preferred list used at the start of the test; the subsequent lists will be automatically selected during the execution of the exam to avoid the use of the same phonemes.

Press this icon to select the masking type sent to the patient trough channel 2. Masking options are: White Noise (WN), Speech Noise (SN) or it is possible to send the same phonemes on masking channel as well (SD card).

By shortly pressing the icon for the output selection (see image on the left) the output can be selected: AC (phones) or BC (bone vibrator) By keeping it pressed for a few seconds you can select the free field output.

You can go back to the previous selection by pressing the icon again for a few seconds.

The same output connectors are used for both headsets and free field.

Before starting the test make sure the appropriate transducer has been connected.

Press this icon to set the time lapse between the reproduction of phonemes from the same list. The interval can be of 2, 4, 6 seconds or determined manually by the operator; in this case the reproduction of the single phoneme will be activated by pressing the **PLAY** icon.

Press the **SETUP** icon to exit the setup page.

### **Output selection**

Gently press directly on the audiogram to change the side of the stimulus: right/left. Masking will be automatically sent to the opposite ear.

### **Intensity settings**

Press the **ARROW** keys to increase or decrease the intensity of the CH1 signal (stimulus) and CH2 (masking). The level reached is indicated between the two arrow keys.

R 17A User's Manual – Page 83

### Signal presentation and storage of threshold levels

Once the speech material and its intensity has been set, press the PLAY icon (see image on the left) to send the signal to the patient. The masking on the opposite ear is always active (except when in "Off" condition).

Press the **REPEAT**, **PAUSA** and **STOP** icon during the test to repeat the last phoneme, pause the test or stop the list reproduction.

Press the **CONFIRM** icon during the test execution to confirm the right patient's response and increase the count of positive responses, don't press if the response is not correct.

These indicators show the numbers of correct answers given by the patient, number of total words submitted and the score as a percentage.

At the end of the phonemic list presentation, press the button shown to save the point on the audiogram. The button pressed changes its colour to green to confirm the registration.

### Communicate to the patient

**TALK OVER** button allows the user to talk to the patient via the transducer. A microphone is integrated in the device.

### Exit and Save

Pressing the **EXIT** button, a dialog box will pop up and ask if you wish to save and/or print the test out on the thermal printer.

Press **YES** or **NO**. If answer is YES, R17A will show patient database page to select the patient to be associated to the test performed, confirmation is required pressing this button. Please, refer to section 6 of this user manual for further information about the patient database management.

### Information shown on the display

Indicates the name of the patient undergoing the examination, as previously selected. Touch the icon to get direct access to patient database.

Indicates battery residual charge and current time.

### 4.3 Autothreshold Audiometry

The automatic threshold test, in addition to manual testing, allows the automatic detection of the hearing threshold by varying the level of the stimulus based on patient's response. R17A incorporates this patient controlled automatic threshold test, complying with EN ISO 8253-1, in a 10 db decrease and 5 dB increase procedure, without masking. This test is a procedure used to determine pure tone thresholds automatically, based on patient's response at each frequency, through pressing on patient's button, connected to the device via nimi-USB socket (9).

How to conduct the test: the patient needs to be instructed on how to use patient's button. The patient will be asked to keep the button pressed when he hears the stimulus and release it when he no longer hears the stimulus.

### **Test execution**

Press the **TONE AUTO** icon (see adjacent image) on main menu to enter the Auto-threshold Audiometry test.

### Stimulus configuration

Press the **STEADY** button to change the type of signal to be emitted: continuous, pulse (2 Hz rate), warble. At every click on the button, the icon changes identifying the type of signal.

Press the **START** icon to start the test, automatically the device will generate a 40 dB HL 1 KHz stimulus. The frequencies will change according this sequence 1, 2, 4, 8, 1, 0.5 e 0.25 kHz; only the frequencies enabled in the setup page will be used.

Test starts with a pure tone at 1000 Hz (40 dB HL) presented for 3 seconds with the purpose to train the patient. Intensity will be increased in 5 dB steps and decrement of 10 dB steps.

When the patient releases the pushbutton (stimulus no longer heard), the automatic threshold detecting procedure starts by presenting pure tones with a duration of 1.5 seconds and random interruption time. When the minimum level (at which the patient is no more able to hear the sound) has been reached, the level is increased and the first valid response is stored in the audiogram as the threshold.

During the test, intensity, frequency and patient response are shown on display

R 17A User's Manual – Page 85

Press the **STOP** icon to stop the test.

Press the **PAUSE** icon to temporarily stop the test.

### Communicate to the patient

**TALK OVER** button allows the user to talk to the patient via the transducer. A microphone is integrated in the device.

### **Exit and Save**

Pressing the **EXIT** button, a dialog box will pop up and ask if you wish to save and/or print the test out on the thermal printer.

Press **YES** or **NO**. If answer is YES, R17A will show patient database page to select the patient to be associated to the test performed, confirmation is required pressing this button. Please, refer to section 6 of this user manual for further information about the patient database management.

### Information shown on the display

Indicates the name of the patient undergoing the examination, as previously selected. Touch the icon to get direct access to patient database.

Indicates battery residual charge and current time.

When the patient presses the respons pushbutton a green line appears on the display to highlight the answer.

### 4.4 Operating position of the device

The way in which tone and speech audiometry tests appear on the screen depends on the orientation of the device: if the device is horizontal (operating position when docked in its support), the screen will appear like the ones describe until now in this manual.

By turning the device by 90° (display is now vertical), after having entered the test conduction page, the screen display will appear like the following:

Functions remain unvaried for tone audiometric and autothreshold tests, in which the commands present in the vertical visualisation are the same as those seen in the horizontal visualisation. The only variable, for these two tests, is that in the vertical visualization, the indications relative to the two sides, left and right, appear on a single graph.

For speech audiometry the way the data is presented changes: in the vertical visualisation, data are in a table. Column AC (or BC) indicates the levels of speech, column M shows masking levels and column S shows the patient's positive answers as a percentage.

Keeping the vertical visualisation in the speech audiometric test, the data in the table refer to the transducer selected in the test setup page (AC or BC). The selection of the side to work on (right or left) is done by touching one of the tables with the data and the active side is indicated by the colour of the bar above the tables.

R 17A User's Manual – Page 87

### 4.5 RE-Play Audiometry

RESONANCE-PLAY is a pure tone interactive audiometry conceived to actively involve the child during the test. The idea is that the child is no longer the patient but a bandmaster, he's not the one tested but the one making the auditions for the elements of a band. A short "fairy tale" has been conceived to create the right atmosphere and prepare the child to the test.

As a bandmaster the child is asked to do an audition to the animals playing different instruments and select the ones working to be added to the band.

By touching the screen, an animal that plays an instrument will be displayed and simultaneously an acoustic signal is sent. If the child hears the sound, the instrument is working and can be added to the toybox by clicking on the relevant image, if the child doesn't hear the sound, means that the instrument is not working and must be trashed.

The child is involved in the test and his attention is focused on the stimuli recognition. At the end of the test a music will be played as a final award. Please, refer to the "Great Bandmaster" leaflet for more information about the test execution.

### **Test execution**

Press the  $\ensuremath{\text{RE-Play}}$  icon on main menu to enter the test. Setup settings are displayed on screen.

A familiarization procedure to train the child should be performed before executing the test.

### **Training selection**

Press the **MODE** icon to choose the training mode.

3 training mode are available:

**Mode 1:** the familiarization is performed with a preset stimuli and frequencies sequence. The device provides a 4 stimuli pattern: a clearly audible level (50 dB HL), a not audible level (-10 dB HL), one at 30 dB HL and finally a not audible level (-10 dB HL). Frequencies are: 1000 Hz, 2000 Hz and 500 Hz. **Mode 2:** the device works as in the real test without saving the result. The child familiarization is performed on the full range of stimuli level.

R 17A User's Manual – Page 88

**Mode 3:** the device emits the stimuli as in "Mode 1" but 20 dB HL higher (70 dB, 50 dB). This training mode can be adopted to test patient with hearing loss.

Start

Press the START icon to start the training.

### Stimulus configuration

Press the **STEADY** button to change the type of signal to be emitted: continuous, pulse (2 Hz rate), warble. At every click on the button, the icon changes identifying the type of signal.

### **Output configuration**

Press the **R** icon to change the side. Options available are: Right (R), Left (L) or both (R/L).

By shortly pressing the icon for output selection (see image on the left), it is possible to select AC (phones) for the threshold estimation, AC UCL (phones) for uncomfortable level or BC (bone vibrator). By keeping it pressed for a few seconds you can select the free field ou-tput. You can go back to the previous selection by pressing the icon again for a few seconds.

The same output connectors are used for both headsets and free field. Before starting the test make sure the appropriate transducer has been connected.

Frequencies settings

Select the frequencies to test by pressing the relevant box. Active frequencies are highlighted with green color, non-active ones are red. 1000 Hz can't be disabled.

Press the **ARROW** keys to increase or decrease the minimum and maximum level for the test. The level is indicated between the two

### **Intensity settings**

Mode settings

### Press the **ARROW** keys to change the mode settings:

arrow keys.

**AUTO** RE-Play is performed in auto-threshold mode by changing the level according to patient response; the minimum heard value will be stored as threshold. 5/10/20 RE-Play is performed in threshold mode by increasing the level of a predefined step from min to max level. DC (double confirm) RE-Play is performed in auto-threshold mode starting at 40 dB in a 10 dB decrease and 5 dB increase procedure. A double confirm at the same level is necessary to store the threshold.

### Frequency and level information

Press the **EYE** icon to enable or disable the display of level (dB) and frequency (Hz) information during the test. This information helps the doctor following the test progress, but in some cases, it could be understood by the patient as well, altering the test result.

### **Preset setup settings**

Preset	
1	
2	
3	

It's possible to store up to 3 different setup settings. Once selected the button be-comes green. It changes to yellow if some parameters are changed, pressing for 5 seconds the new settings are stored as new protocol. The factory presets are:

1)0.5, 1, 2, 4 Khz; 20 dB HL Min; 20 dB HL Max; Auto mode.

2)0.5, 1, 2, 4 Khz; 0 dB HL Min; 80 dB HL Max; DC mode.

3)0.5, 1, 2, 4 Khz; 0 dB HL Min; 80 dB HL Max; Auto mode.

Press the **START** icon to start the test.

### **Test execution**

In the displayed grid, different tiles represent different frequencies; the numbered tiles are related to the enabled frequencies to test, the tiles with an animal image refer to the disabled ones. During the test all the tiles will be filled with animal images in order to complete the puzzle.

R 17A User's Manual – Page 90

Touch the screen to start the test.

Touch the animal playing an instrument (shown in the center of the display) to emit the stimulus. An animation will be displayed to entertain the child as a stimulus reinforcement.

At the stimulus presentation, two options are available on screen:

The patient has to touch the toy-box if he has heard the stimulus, otherwise the trash bin if the stimulus hasn't been heard.

The test proceeds by changing the stimulus level according to the selected setup settings.

In case of error, the doctor can undo the last response and go back to the previous condition by pressing the patient pushbutton (connected to the device).

After the first frequency has been tested, the grid displays a new completed tile (animal image), touching the display, the test will proceed with the next frequency. At the end, an animated image di-splaying all the animals will be shown as reward and a melody will be played.

Press **EXIT** button, a dialog window will ask if you wish to save the test. Please, refer to the other audiometry test sections for more information about results saving.

On the **SETUP** page, a button to display the audiogram is now visible. Press this button to show the test results.

Below the audiogram graph, an indicator of the test reliability will be displayed. It is based on response time and on the response to fake stimuli randomly presented:

Reliable

Partially Reliable

Not Reliable

### Exit and Save

Press **EXIT** button, a dialog window will ask if you wish to save the test. Please, refer to the other audiometry test sections for more information about results saving.

### 4.6 Weber

Weber is a screening test generally used to check for unilateral hearing loss. The test is carried out by placing the bone vibrator on the patient's forehead, sending the stimulus to a comfortable level and asking the patient from which ear he perceives it. The test can be performed at different frequencies.

After placing the bone vibrator on the patient's forehead, select the desired stimulus and set a level that is comfortable for the patient. Send the signal by pressing the SIGNAL button asking the patient which ear hears the sound from. Press the corresponding button according to the answer: right ear, both (center), left ear or not heard.

It is possible to perform the test at different frequencies and save the results.

### 4.7 Otoscopy

The R17A Audiometer can provide a very efficient tool to analyze the external ear conditions, when coupled with the RE-VO otoscope for the image acquisition of the external ear-canal and of the tympanic membrane. Like all the other audiological tests, images can be stored in the local database and later on transferred into the MDS Software.

### Connection

The otoscope is equipped with a dedicated cable: all you have to do in order to connect it to the device is inserting the mini-USB connector in the otoscope and the USB connector in the R17A USB Host; it automatically switches on when an examination is started, and you do not need to press any buttons to start.

### **Speculum insertion**

The otoscope is supplied with a set of multiple size disposable specula making it suitable for both young and grown-up patients.

Place the speculum on the otoscope tip and press it until it is fully inserted. In case the image appearing on the screen is not centered, softly move the speculum tip around, until the correct position is achieved.

Specula are single-use and need to be replaced each time after using. Only specula supplied by MRS can be used. In case any other part of the device comes into contact with the patient, just disinfect. The unit is made of biocompatible materials.

### **Test performance**

In order to start the test, switch the device on and press the icon illustrated on the side, on the screen.

When the test is started all LEDs light up. Brightness is automatically adjusted to obtain an optimal image quality however; you can also manually change it at your will by using the + and – buttons positioned on the handle. By keeping the – button pressed for a few seconds you can switch from 2 LEDs mode to 4 LEDs mode and vice versa

Hold the otoscope with its focus ring facing downwards so that the image appears straight (i.e. not rotated). The focus ring allows for both focus adjustments and image blowing-up.

Softly insert the speculum into the ear canal and adjust the focus using the focus ring until you obtain a sharp image.

For capturing and saving the image into the gallery, press the screen on the image or icon illustrated on the side. You can save up to 9 images for both the right ear and the left ear. In case you capture more images, the last 9 ones will be saved overwriting the previous ones.

To switch side, press the **R** or **L** icons.

To review the captured images, press the Gallery icon.

The gallery shows all saved images and in order to see them in extended view you just need to touch them. Once a full-screen image is obtained you can further blow it up by double -clicking on the screen. The zoomed image will become larger than the screen, but you can still see it all by keeping your finger firmly on the screen while rotating it in the different directions.

To go back to the full-screen image, double-click on the screen again.

R 17A User's Manual – Page 94

By pressing the Gallery icon you go back to all saved images in small size.

Press the Otoscope icon to leave the Gallery and go back to the otoscope main screen.

In order to delete an image keep your finger pressed on the small image you wish to delete for a few seconds.

### Test saving and exit

Pressing the **EXIT** button, a dialog box will pop up and ask if you wish to save and/or print the test out on the thermal printer.

Press the icon corresponding to your wish: if you confirm you want to save the test, the R15C will automatically address you to the patient database where the patient being tested can be selected and confirmed. For additional information on the patient database, please refer to section #6 in this manual.

## 5. Setup settings

Start Start Tore Audiometry Speech Audiometry France Fran

After switching on the device, set the parameters that affect the test operation.

Press the **SETUP** icon (see adjacent image) to enter the setup page, there the settings are grouped as follow: GENERAL; AUDIOMETRY

### 5.1 General Setup

Press one of the tabs to access the respective Setup page:

### Bluetooth

Press the **BLUETOOTH** icon (see adjacent image) to activate database data transmission to PC, if set to OFF transmission is disabled.

R 17A User's Manual – Page 96

#### Autosave

Press the **AUTOSAVE** icon (see adjacent image) to enable the autosave mode. If enabled, once selected a patient, every performed test will be automatically stored into the data-base after pressing the **EXIT** icon; if disabled a pop up window will be displayed asking for confirmation.

#### Language

Through the **LANGUAGE** button (see adjacent image) you can select the interface language. By pressing this icon a menu with all the available languages will be opened.

#### **Rotation block**

Some of the tests provide for their visualization either horizontally or vertically: the visualization is automatically set by rotating the instrument. This option allows to block the rotation keeping the display horizontal.

Wi-Fi

By pressing the **WIFI** icon (see adjacent image) a window allowing the connection to a Wi-Fi network will. This icon is active only when a Wi-Fi USB dongle is inserted in the USB host port. This option can be used to check the connection with the wireless printer.

#### **Power management**

In order to improve the battery operating autonomy, you can set 3 different power modes: energy saving, balanced mode and performing mode. Based on the selected mode the display backlight and the unit automatic switching off will change. When the unit is connected to the power supply and the performing mode is selected, the display remains constantly backlit.

### **PIN code change**

By pressing this icon, you can change the PIN code by typing in any 4-digit numeric code or disable the PIN code request at power on. Please note that in order to change or disable the PIN you need to type in the one currently in use first.

#### Patient data mandatory fields

Pressing this key gives access to a menu where it is possible to define the minimum fields necessary for the creation of a new patient. It is possible to set only a numeric code or activate the fields like name, surname, etc.

### Print-out

By pressing the **PRINT** icon (see adjacent image) you can set the print options for the optional portable printer.

You can select the printing-out of the following fields:

- Unit serial number and logo
- Patient's data
- Calibration expiry date
- Lines for comments after the heading
- Lines for comments after every test

I		09:44
Portable Printer		
Status: Disconnected		- l
Print patient data		
Print next calibraton date	Report Logo	$\sim$
Print lines for notes after each test		
		D.

### **Report header customization**

You can customize the report header you print by importing a graphic file that contains your logo or information about your company. The graphic file must be in BMP format, black and white (1 bit x pixel) and have a maximum size of 600x400 pixels. Create a folder named PrinterLogo in the root of the pen drive, copy the file or files to be used inside and insert it in the USB HOST port of the device. It is necessary to copy the file to the device memory to be able to use it for printing.

Pressing the LOGO button (see adjacent image) will open the page for selecting the file

Using the **ARROW** keys, it is possible to choose the graphic file from those present both in the device and in the pen drive

If the selected file is on the pen drive, the + key will turn orange, pressing it will copy it to the device memory.

Once the desired file has been selected, press the **CONFIRM** key (see adjacent image) to set it as the print header.

The **DELETE** button will allow you to delete the selected file.

If the printer is connected, using this key it is possible to print the test page with its logo.

### Licenses

Press the **LICENSES** icon (see adjacent image) to check the active test licenses. The following is displayed:

To add a new license, press the **ADD NEW** button and enter the license code:

		100% 4 10:54
	ENTER YOUR KEY	
		Add new
1	2 3 4 5 6 7 8	90 🖾
QW	VERTYUIC	DPĞÜ
<b>→</b> I	A S D F G H J	K L Ş İ
z x	K C V B N M Ö Ç	· • •
ÀÜÊ		. , /

The **ADD NEW** button should be pressed again to confirm.

The device recognizes the validity of the license and will display a confirmation message stating the test that has been added.

An error is displayed if the license is not valid.

If there are no active licenses, the tests can be used in demonstration mode up to 10 times.

The number in brackets shows the number or demonstrations remaining.

### Date and time

Press the DATE icon (see adjacent image) to set date and time.

Use the arrows to adjust the date and time.

It is possible to change the date format and related separator characters using the appropriate field using the left-right arrows

Press this icon to confirm.

### About

Pressing the **ABOUT** icon displays the devices' serial number, firmware and software versions. In addition, when entering this page, the device executes an internal test of the main hardware functionalities and shows the results.

61	7 <b>a</b>	
Firmware Release Software Release	200 200	
Product number: Bluetooth address:	4200104010 00:00:39:02:03:27	resonance 🔊
Serial number:	III7AProduct11	
Coliberting environ	837A 88/01/2018	by M.R.S. a.r.L
	Margaret .	Via Cesore Batteti, 134
Battery charge	30%	24025 Gerange - BU
DSP status:	OK .	
SD CARD atotus:	05	www.resonance-audiology.com
CODEC status	ok.	upport≣misonance-audiology.com +39.035.712091
Reliability test:	ok :	

R 17A User's Manual – Page 101

### 5.2 Audiometry Setup

Press the AUDIOMETRY tab to enter the respective setup page:

### Frequencies

Press the **FREQUENCIES** icon (see adjacent image) to enable the frequencies to be used during the Pure Tone audiometry test.

Enabled be	quencies			
s 🔲	.en 🛄	.vai 🛄	ra 🛄	NARE 🔲
2	,A930	***	5-1841 🔯	we 🔲
	****	soge 🔯		
		0		

Select the frequencies and press the **CONFIRM** icon (see adjacent image) to save these settings.

The TRACK icon allows to activate or deactivate tracking between two channels: ON varies the level of the signal in channel 1varying to the same extent the level of channel 2

### Steps dB

Press the STEPS dB icon to set the increase of the level of the signal of the two channels (1 dB, 2 dB, 5 dB and 7 dB)

### Display

Press the **DISPLAY** icon to choose the audiogram background during tone audiometry: visualisation of the maximum level at different frequencies, visualisation of the hearing loss level, using grey bands of different intensities (according to the BIAP 1997 – Bureau International di Audiophonologie classification); visualisation of vocal speech area (Speech Banana).

### Standard

Press the **STANDARD** icon to set the device to perform the tests according to the European norm IEC 60645-1 or ANSI S3.6

### Headphones

Press the icon to choose the headphones type: DD45, TDH39, HDA280, IP30

This function allows switching between different contralateral transducers. Ensure that any transducer to be connected is calibrated to the device. To change the settings, keep the button pressed for some seconds.

### Symbols

Press the **SYMBOLS** icon to choose the symbol for the marker on the audiogram (Tone or speech) that identifies the right-side bone vibrator. (not active on R17A series)

### Stimulus level keys functionality

The button allows you to change the functionality of the arrow keys to change the level of the stimulus: up arrow to increase and down to decrease or vice versa.

### Calibration

Press the **CALIBRATION** icon (see adjacent image) to enter the audiometry calibration page. Only authorized personnel are allowed to perform calibration.

### 5.3 Audiometry Calibration

Entering the setup audiometry menu and pressing the **CALIBRATION** icon (see adjacent image), you will be able to access to the calibration page.

Only authorized personnel are allowed to perform a new calibration.

The device calibration must be performed by trained technical staff only. Information concerning how to perform a calibration are available on the Technical Manual

## 6. Patient database management

The device has an internal memory which can store patient data and test results. The database is accessed by pressing the **PATIENT** button on the home screen.

In the database, all data stored can be displayed.

### Database

The left column lists patients, the center column shows test sessions associated to the patient, and the right-hand column displays demographics.

Using the **ARROW** keys to navigate between patients, results for the highlighted patient will be displayed

Patients list can be sorted by name or company name (if inserted in the patient data). Press the icon shown on the left to change the sorting option.

### **Patient data**

The right side of the screen displays the selected patient's personal data. The lower part contains some icons that activate the specific functions described below.

### Create new patient.

Use the + button to add a new patient to the database. When executing this command, a virtual alphanumeric keyboard will appear which allows data to be entered.

To display the accented characters, simply press and hold the key for a few seconds: a popup will appear with the alternatives (e.g. hold down the A key to enter À, Á, Â, Â, ...) Using the options key, you can select different keyboard layouts with the corresponding character sets (e.g. Qwerty, Arabic, Cyrillic ...)

By shortly pressing the **DELETE** icon you will delete the highlighted test. Once all tests have been deleted, the patient can be removed from the database by pressing the DELETE icon again. As an alternative you can delete a patient and all the related tests by keeping the **DELETE** key for a couple of seconds and then confirming the action.

Multiple test elimination: if you wish to delete multiple tests belonging to the same patient, select the patient, keep the first test you wish to delete pressed for a few seconds and until it turns orange, then add the other tests you wish to delete by touching them shortly. When all the selected tests will turn orange, you will be able to delete them all by pressing **DELETE**. The action can be aborted by pressing the Exit button.

Patient data editing. Pressing the EDIT icon allows data editing for the selected patient.

Test results display. Press the VIEW icon (see adjacent image) to display a specific test result for the selected patient.

**Patient activation**. Press the icon shown on the left to activate the selected patient. The selected patient's name will be displayed on the top left of the display once you enter one of the pages that allow the tests execution. Tests that have been performed will be saved to the associated patient.

## 7. Print management

The test results can be printed through the PC Suite Resonance MDS, or directly from the unit to the optional portable wireless printer.

To print out via the MDS, Transfer the data from the device to the PC Suite Resonance MDS and then follow the instructions related to the "Session report" MDS function.

If you wish to print directly from the unit, switch the printer on and insert the Wi-Fi dongle in the R17A USB port. Please note that the two devices require 20-30 seconds to connect.

Once unit and printer are connected, enter the patient database menu, select the test you wish to print out and press the **PRINT** icon you will see on the display. It is possible to print report during the text exit by pressing the print button on the pup up menu.

The wireless printer uses the 2.4 GHz Wi-Fi technology. Using it in areas where too many Wi-Fi networks, or other sources of electromagnetic interference such as mobile phones or other similar devices can be found, can provoke malfunctions.

## 8. Data transmission to PC

Data transfer from the device to PC is possible via Bluetooth connection.

To transfer test results stored within the R17A and subsequently view, edit, print or store them in a NOAH database, the PC suite Resonance MDS must be installed on to a computer. The installation file named **MRS Resonance Setup.exe** is supplied on the USB pen-drive which includes this operating manual. Refer to the installation & operating instructions provided with the PC suite Resonance MDS for further details.

### How to associate R17A to PC via Bluetooth

- 1. Activate the Bluetooth connection on the device by pressing Bluetooth button present in the General Setup page (see section 5.1 of this manual).
- 2. Make sure the Bluetooth data transmission is active on your PC.
- 3. Open **Windows settings** on PC.
- 4. Select the **Devices** icon.
- 5. Open **Bluetooth** devices management feature.
- 6. Wait until your PC detects the device. This may take several minutes. When finished, an icon showing the serial number of the device to be associated will be displayed.
- 7. Check that the serial number displayed corresponds to the one on the device identification label. If it matches, press the icon on screen to complete the association.
- 8. The unit will generate a **PASSCODE** which you will see displayed on the screen.
- 9. The **PASSCODE** will be displayed on the PC screen, too.
- 10. If the two codes are the same, press the confirm button on the unit first and then on the PC.
- 11. When this process is finalized, the unit is associated.

As an alternative, data can be transferred via cable through the unit mini USB port. If you choose this option be aware that by connecting a Personal Computer to a Medical Device you create a system which, if used in a patients' area need to comply with all safety requirements.

In case the R17A is connected to other devices equipped with their own external power supplies, the compliance of the whole system to the IEC 60601-1 standard needs to be verified under the technician's connecting the units responsibility. In case insulation devices are used, they have to comply with the IEC 60601-1 standard. Each person connecting an R17A to other external units creates a system and is held responsible for the compliance of the system itself to the IEC 60601-1 standard.

## 9. Main Labels

### Device labels

4 34 O

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RIGHT SIDE

FRONTAL VIEW

LEFT SIDE

BOTTOM HOUSING
## 10. Messages and warnings

Save the test in the database?

Is it possible to perform the print out. The button is active only if the printer is active and connected to the Wi-Fi adapter.

Save current settings in setup?

Update also the calibration date?

Calibration expired. Recalibration required. Continue anyway?

Communication to patient; level setting

Save this setting as default?

R 17A User's Manual – Page 109

Delete the patient from database?

Attention, optional test license is not correct

Optional test license is correct.

Attention, battery is running low. Connect the device to its medical power supply to recharge and continue the test.

Delete the selected point?

The free memory for saving exams is running out. It is possible to free up memory by transferring the data to the MDS suite and deleting them from the device

## 11. Care and maintenance

Switch off the device and disconnect the power plug before cleaning!

Always disconnect the mains conductor during the cleaning process, and be careful that no fluid enters inside of the device or the accessories; no alcohol or spirits should be used.

The device requires cleaning the housing with a damp cloth (do not use aggressive cleaning agents); avoid liquid ingress.

During testing operations, headset pads and probe tips come into direct contact with the patient. Although said components are biocompatible, they are not sterile and they must therefore be subject to special care:

After each test, it is necessary to disinfect the headset cushion before use on a new patient. This includes physically cleaning using of a recognized disinfectant. Individual manufacturer's instruction should be followed for use of this disinfecting agent to provide an appropriated level of cleanliness.

The insertion headsets tips are disposable: use only once and substitute them before use with a new patient.

Calibration and any possible repair, installation and/or update of the device must be strictly carried out only by specialized technical personnel, duly authorized by the R17A manufacturer and in full compliance with the terms and conditions specified in the Technical Manual.

Standard maintenance and device calibration operations are to be regularly performed at one-year intervals. Standard maintenance operations cover the inspection of the transducer cable conditions, of the mains power supply cable, of the earphone pads as well as an inspection check to ensure the device's outer section and/or cabinet are in good working order. All device calibration operations must be strictly performed only by qualified technical personnel, duly authorized by the manufacturer and implementing the appropriate apparatus for checking the levels and frequencies of all the output signals by the transducers supplied with the device. Once device calibration testing using the special purpose apparatus is complete, it is necessary to run a device function test.

Should transducers require replacement, use only parts provided by Resonance. Upon replacement of transducers, the entire system will require re-calibration

The device is provided with a lithium-ion battery type 1SP NCR-18650A. Said battery is fitted into the device and possible replacement thereof must strictly be performed only by qualified technical personnel. In case of replacement, use only the same battery type.

## 12. Safety regulations

## 12.1 Electrical Safety

The device applied parts comply to type BF of IEC 60601-1, in any case, the device is not intended to be used in environments dealing with explosive material or equipment.

## 12.2 Device handling and operation

The device should be checked once daily, as per operations described on section 11. The device is designed to be handled and operated by trained personnel such as Audiologist, ENT doctors or personnel with similar qualifications ONLY.

## 12.3 Measurement Safety

As a manufacturer, we recommend to perform an annual inspection and calibration in order to ensure safety and quality of the measurement.

Annual calibration and maintenance should only be carried out by RESONANCE authorized personnel. M.R.S. will not be liable for any failure to comply with the specified inspection date, according to the Medical product law. The use of uncalibrated and uninspected device is strictly prohibited

## 13. Technical specifications

Model	R17A AC Screening	R17A BC Basic Diagnostic	R17A HDA Basic Diagnostic
Code	4300104300	4300104310	4300104305 4300104315 4300104325
Test			
Pure Tone Audiometry	•	•	•
Speech Audiometry	•	0	•
Autothreshold	•	0	•
Re-Play	0	0	0
Weber	•	•	•
Accessories			
DD45 headset*	•	•	-
HDA 280 Headset*	-	-	•
IP30 Insert earphones*	0	0	0
Bone vibrator B71W/BC-1*	-	•	•
Internal SD card with phonemes	•	•	•
USB pen-drive with User Manual	•	•	•
Medical CE approved power supply	•	•	•
Cradle	•	•	•
Intrenal Bluetooth module 2.1/4.0	•	•	•
MDS software	•	•	•
Carrying bag	•	•	•
External USB keyboard	0	0	0
USB Cable	0	0	0
Wireless printer	0	0	0
RE-VO Otoscope	0	0	0
Free field amplifier	0	0	0
Passive free field loudspeakers	0	0	0

•= standard  $\circ$  = optional - = not available \* = applied parts

Applied parts are constituted by the transducers that convey the test signals. Said transducers are applied so that they go into direct contact with the patient. In the foregoing list, they are highlighted in bold script and are counter-marked with an asterisk (\*).

Refernce standard:	IEC 60645-1; ANSI S3.6 Type 4
Safety:	IEC 60601-1 (2006) + A1 (2012) Class II Type BF
EMC:	IEC 60601-1-2 (2014)
Medical CE Marks	Class IIa The Manufacturer's Quality Management System complies to the Annex II of Medical Device Directive 93/42/EEC. The notified body IMQ, Identification N. 0051, has approved the Quality System.

Pure Tone Audiometry			
Speech audiometry with speech material stored in the device internal memory			
Autothreshold audiometry with automatic threshold search (model depending)			
Outputs	Air conduction/ right free field ACR, left ACL, BC		
Inputs	Patient pushbutton (model depending)		
Tipe of signals Pure tone, warble, Speech material on internal memory			
Type of signal presentation Continuous, pulsed (2Hz rate)			
Masking Narrow band noise, speech noise, white noise			

	DD 4	45	HD	A 280	Insert p	phones	Free	e Field	B71W BC-1
Freq. Hz	Pure Tone (dB HL)	NBN (dB EM)	Pure Tone (dB HL)	NBN (dB HL)	Tone (dB HL)	NBN (dB EM)	Tone (dB HL)	NBN (dB EM)	Tone (dB HL)
125	65	40	70	55	80	55	70	60	-
250	85	60	90	70	90	70	80	70	45
500	105	75	110	85	100	80	90	80	65
750	110	80	110	90	110	85	90	80	70
1000	110	80	110	90	110	90	95	80	70
1500	110	80	110	90	110	90	95	80	70
2000	110	80	110	90	110	90	95	80	70
3000	110	80	110	85	110	90	95	80	70
4000	110	80	110	85	100	90	95	80	70
6000	105	80	100	75	95	80	80	70	50
8000	90	80	85	75	80	75	70	70	-
10000	-	-	80	70	-	-	-	-	-
12500	-	-	80	70	-	-	-	-	-
WN	100 (dE	SPL)	100 (	dB SPL)	110 (d	B SPL)	75 (c	B SPL)	-
SN	80	)		90	8	0		70	-
SD Card	90	)		90	8	0		75	45

Ambient noise monitor – MPANL						
Frequencies	DD45	IP30	HDA280	ADC, (dB)		
Hz	(dB)	(dB)	(dB)	DD45 only	ROOM1	ROOM2
125	58	78	60	+6	+21	+25
250	37	55	42	+8	+22	+30
500	18	37	26	+17	+22	+33
1000	23	36	23	+15	+28	+37
2000	30	36	19	+13	+29	+40
4000	36	41	34	+13	+36	+47
8000	33	47	34	+11	+40	+50

Note; The limits indicated are valid for tests from 500Hz to 8KHz, the values for the frequency bands of 125 and 250Hz are reduced in the case of tests with tones at frequencies of 125 and 250Hz.

Technical descriptions more detailed and complete as wiring diagrams, parts list, instructions for calibration and other materials necessary for the repair and maintenance of the instrument, are contained in the service manual which will be provided, upon request, technical personnel authorized by the manufacturer to carry out such operations

Dimensions	W x D x H: 180 x 135 x 41 mm	Weight: 520 g net
Environmental conditions		
Transport and storage:	Temperature: from -20 to *50°C Humidity: from 20 to 90% without condensation Pressure: from 500 hPa to 1060 hPa Temperature: from 20 to 90% without condensation Pressure: from 700 hPa to 1060 hPa	
Power:		
Batteries:	Internal rechargeable Lithium-lon battery 3,6V- 3,1 Ah Battery status indication	
Battery duration:	Normal use, up to 8 h Continuous use (test performance): 3 h	
Battery charge time: External power supply:	Less than 4 h; 80% in 2 h Medical power supply input 100/240 VAC 50-60 Hz- output +5 VDC/2A	
Construction	Self-extinguishing plastic ABS cabinet	

## 14. Electromagnetic compatibility

Guidance and manufacturer's declaration- Electromagnetic emissions

The R17A is suitable for use in the specified electromagnetic environment. The purchaser or user of the device should assure that is used in an electromagnetic environment as described below:

Emissions test	Compliance	Electromagnetic environment
RF emissions CISPR11	Group 1	This R17A uses RF energy only for its internal function. Therefore, the RF emission is very low and not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR11	Class B	This R17A is suitable for use in all establishments, including domestic
Hermonic emissions IEC/EN 61000-3-2	Complies	establishments and those directly connected to the public low-voltage power supply network
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	that supplies buildings used for domestic purposes.

Guidance and manufacturer's declaration- Electromagnetic immunity

This R17A is suitable for use in the specified electromagnetic environment. The purchaser or user of the device should assure that it is used in an electromagnetic environment as described below:

Immunity test	IEC 60601-1-2 Test level	Compliance level	Electromagnetic Environment
Electrostatic discharge (ESD) IEC/EN 61000-4-2	8 kV contact 2/4/8/15 kV air	IEC 60601-1-2 Test level	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Radiated electromagnetic field IEC 61000-4-3	3 V7m 80 MHz to 2.7 GHz	IEC 60601-1-2 Test level	Portable and mobile FR communications equipment should be used no closer to any part of the R15C including cables. Minimum distance 30 cm.
Electrical fast transient/burst IEC/EN 61000-4-4	2kV for power supply lines 1 kV for input/output lines >3m	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	0.5/1 kV differential mode 0.5/1/2 kV common mode	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	3V 150 kHz to 80 MHz	IEC 60601-1-2 Test level	Portable and mobile RF communications equipment should be used no closer to any part of the R17A including cables. Minimum distance 30 cm.

Voltage dips short interruptions and voltage variations on power supply imput lines IEC 61000-4-11	0% Un. For 0.5 cycle 0% Un. For 1 cycle 70% Un. For 25 cycles 0% Un. For 5 s	IEC 60601-1-2 Test level	Mains power quality should be that for a typical commercial or hospital environment. If the user of the R17A requires continued operation during power mains interruptions, it is recommended that the R17A be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	30 A/m	IEC 60601-1-2 Test level	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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