
Audiometer

QUICK USER'S HANDBOOK

R279R379

CE
0051

resonance 

Release history

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this page intentionally left blank

Table of contents

Operating options and Scope of use	32
Safety warnings and precautions	33
Controls and components - front section	34
Connectors and components - rear section	36
Installation	37
Parameter Setups	38
Performance of Audiometry Tests	38
Loading paper in the thermal printer	43
Print management and .pdf report	44
Main labels	44
Symbols present on the device	45
Messages that may show up	46
Maintenance and cleaning	46
Environmental conditions	48
Main technical specifications	48
Models	50
Accessories, test and applied parts	51
Electromagnetic emission	52

IMPORTANT NOTICES

Resonance® thanks you for having purchased one of its medical devices!

Resonance® hereby advises that in order to ensure optimized device performance and full customer satisfaction, prior to proceeding with device installation and use you are kindly requested to carefully read and understand all the notices and warnings given in this User's Handbook, as they provide important directions on device safety, care and maintenance.

Carefully store this manual away in a safe place for each future consultation.

This device has been manufactured and factory tested, fully compliant to product standards and regulations and has been delivered from its manufacturing premises in perfect technical safety conditions. Once the device has been removed from its packaging, please check to ensure it has not been damaged and that no parts are missing.

In case of any doubt, please do not use the device and either contact your Resonance Dealer, or e-mail the Resonance Customer Support Service directly at support@resonance-audiology.com

If the device drops or falls accidentally and/or is subject to violent impacts, it may also undergo hidden and/or invisible damages, incurring in hazardous device operation.

Operating options and Scope of use

The R27A/R37A measuring unit provides for clinical audiometric testing.

The extensive battery of clinical audiometry tests includes:

Pure Tone and Speech Audiometry, Auto threshold testing, ABLB, DLI, MLB, SISI, Tone decay testing, Békésy automatic audiometry, Multi-frequency audiometry, Gap Detection Thresholds, Stenger test performance, QuickSIN, DLF and High Frequency audiometry up to 20000 Hz (only for R37A-HF model).

Use of the unit is reserved to audiologists, medical and technical hearing aid and conservation professionals and audiology specialists.

The unit is meant for use in clinical environments and at low ambient noise levels. For some types of tests (audiometry), very low ambient noise levels are required.

During patient testing processes, operators are required to exercise maximum care and attention. An operating ambient temperature ranging between 15 and 35 °C (59 - 95 °F) is recommended.

Safety warnings and precautions

Any possible modification and/or alteration to this device is prohibited, without prior authorization thereto by the Manufacturer.

In the event that the R27A/R37A is connected up to other devices having individual, external power supply sources, it is necessary to test compliance of the overall system to the EN 60601-1-1 Standard provisions. Should the implementation of any insulation or grounding devices be required to meet conformity, said devices must also conform to Standard EN 60601-1.

During test performance, always check the unit's stimulus frequency intensity setting prior to conveyance to the patient. This is to ensure application of acceptable patient stimulation intensity levels.

Prior to starting a new patient, proceed with cleaning the sections of the transducers that come into direct contact with the patient (e.g. earphone pads). All cleaning operations need to be performed using an appropriate antiseptic solution, compliant to the instructions provided directly by the antiseptic product Manufacturer. Please also consult the maintenance and cleaning section in this manual. The R27A/R37A device is not suitable for use with flammable anesthetics. The R27A/R37A is furthermore not provided with protection against the ingress of liquids: the possible penetration of drops and/or liquid spray substances may result in serious damages occurring to the unit.


The unit must always be used with the unit's own transducers, provided on-supply by the Manufacturer. Please be advised that the unit has been calibrated for use with the transducers that are supplied. If they are replaced, recalibration of the unit is necessary. In the event that the transducers are connected to the device through extension cables, as is required in the event of connections to soundproof cabins, it is recommended to proceed with unit recalibration.

Although the unit is compliant to the EMC (electromagnetic compatibility) regulations currently in force, prevention of unit exposure to other electromagnetic sources such as cell phones is highly recommended. Should the unit be in the vicinity of or placed alongside other devices, prior to using the unit please ensure that there are no interferences between the two systems.

Controls and components

front section


1. 7", color graphic TFT display module.
2. Built-in thermal printer.
3. Integrated microphone module for patient talk over (automatically excluded upon hookup of the external microphone into the appropriate jack connector (26) located on the back panel).
4. Channel 1 (dB) and 2 (dB) adjustment knobs: the channel 1 adjustment knob enables channel 1 dB level adjustments and selection of the items on the menu pages. The channel 2 adjustment knob enables channel 2 dB level adjustments and selection of the menu page items. Both adjustment knobs enable sound level controls via the volume control menu page (CONTROLS key 14).
5. F1.....F6 function keys: enable execution of the function displayed on the bottom

section of the screen frame, just above the corresponding keys. Said functions will vary according to the type of test being run or the menu that is displayed.

6. TEST: when in testing mode, this key enables return to the main selection page.
7. ENTER: key for accessing the tests from the main menu.
8. RIGHT – LEFT inputs for channels 1 and 2: these keys enable selection of the ear being tested for both channels. The selected stimulus signal is thus sent to the transducer corresponding to the selected ear.
9. OUTPUT – INPUT for channels 1 and 2: enables individual selection of the type of stimulus and the transducer used for testing.
10. ARROW KEYS: used for the various selections provided on the menu pages. For multifrequency audiometric testing they are used to set and enter the frequency step and to change the base frequency setting.
11. NORM/REV for channels 1 and 2: these keys enable stimulus transmission to the transducers. Normal mode (output signal present when pressing SIGNAL key 12). Reverse mode (output signal not present, when pressing SIGNAL key 12).
12. SIGNAL for channels 1 and 2 (touch keys): in audiometry testing, with the channel in normal mode, these keys enable transmission of the respective channel signals to the transducers for the length of time they are pressed. In Reverse mode, these keys deactivate the transducer output signal for the length of time they are pressed.
13. FREQUENCY: these keys enable selection of the signal frequency sent to the transducers.
14. CONTROLS: this key provides access to the external signal adjustment and calibration page, for the monitor-headset, for the patient-to-operator mike (e.g. for talk back) and for the operator-to-patient mike (e.g. for talk over).
15. T. OVER: this provides direct access to the talk over microphone volume adjustment controls (for operator to patient communication).
16. START/STOP: this key enables start/stop of some of the test options (for Autothreshold, SISI, Tone Decay, Bekesy and Speech Audiometry in clinical audiometry testing).
17. STORE: this key enables storage and display of the values measured during audiogram testing.
18. PRINT: if pressed from the main selection page, this key enables printing of all the tests performed (via thermal printer, or by creating a pdf report), otherwise printout of a single test is possible if the key is pressed directly from the test page (with printout via the built-in thermal printer).

Connectors and components

rear section

19. PHONES (R – L): n.2 6.35 mm stereo jack outputs for hookup of the DD45, TDH39, HDA280 or Insert earphones
20. HF PHONES: two 6.35 mm stereo jack outputs for hookup of the high frequency headphones (optionals) for R37A-HF model only. For models R27A and R37A manufactured since 2014 (detectable by the serial number of the audiometer), the two outputs are reserved to the insert headphones. For audiometers manufactured before 2014, the outputs are unused and closed with appropriate caps
21. INSERT: 6.35 mm stereo jack output for hookup of the insertion transducer for sound masking
22. BONE VIBRATOR: 6.35 mm stereo jack output for hookup of the bone conductor unit
23. FREE FIELD (R – L): RCA outputs for amplifier hookup, required for performance of tests operating in free field mode

24. CD TAPE (EXT1 – EXT2): RCA inputs for connection of external source signals
25. TALK BACK: 3.5 mm stereo jack input for hookup of the patient microphone
26. T. OVER: 3.5 mm stereo jack input for hookup of the operator microphone
27. MONITOR: 3.5 mm stereo jack output for hookup of the monitor-headset
28. PAT. RESP: 6.35 mm stereo jack input for connection of the patient response pushbutton.
29. PAT. RESP. 2: Input 6,35mm stereo jack. For audiometers manufactured since 2014 (detectable by the serial number of the audiometer), the input is reserved to the second patient response button (optional). For audiometers manufactured before 2014, the input is unused and closed with appropriate cap.
30. SD CARD: Slot for the external SD card; supplied with the device, providing a broad range of multilanguage speech recordings used for speech audiometry testing
31.  : USB host port for hook-up of a USB keyboard and/or an external USB data storage pen drive
32. : USB slave port for PC communication and data transfer
33. ON-OFF switch: on-off switch with fuse holders and mains power supply cable connection
34. Operator goose-neck microphone connection: area provided for attaching an optional operator microphone

Installation

Prior to switching the device on, please proceed with connecting the transducers required for testing.

The following is required for performance of the main audiometry tests:
Headset (DD45, TDH39, HDA280 or Insert), Bone conductor, Patient response pushbutton.

Connect the power supply cable to the device, checking to ensure that the cable is turned correctly to connect to the power mains, then press on the rear switch (33) to power the device on.

According to the EN60601-1 Standard provisions, the device falls under class 1. It is therefore necessary that the power grid, which it is connected, is provided with a grounding connection.

Parameter Setups

Once the device is on, proceed with the principal settings involved with test performance.

Press the F1 key (5) to access the setup menu. Via the adjustment knobs (4) or the down arrow, select the PARAMETER SETUP function and then press ENTER (7). Use the side arrow keys (10) to select the various options and/or to change selected options. Use the down arrow key to move to the next setup function.

On SETUP page, available options are grouped in separate sections.

Device general parameters

Language: changes the end-user interface language.

Date and time: changes the system date and time that is currently displayed.

Audiometry parameters

AC headphone: selects the type of transducers to be connected to the unit (DD45, TDH39, HDA280 or Insert).

BC right side and BC masked: selects the type of symbols required for the bone conductor graph mapping.

BC marks joined: if active, the points on the bone conduction audiogram will be connected, otherwise the points will not be connected.

Autom. Change Freq.: if active, by selecting a point on the audiogram with the STORE key (17), the frequency of the test will change automatically. If not selected, to change frequency the tester will have to manually change after pressing STORE (17).

Test Data: if active, Tonal and Speech Audiometry results are printed both in graph and numeric format.

Enable Frequencies Audiometry: activates or deactivates the available frequencies only for pure tone audiometry tests.

The setting entries can be confirmed by pressing F5 (5). Time and date settings can be stored by pressing F3 (5). Press F1 to return to the main Test Selection page.

Performance of Audiometry Tests

Once the device has been powered up, audiometry tests can be selected via the display by acting on the CH2 adjustment knob.

Pressing F3 soft key, additional Audiometry Tests are displayed.

Test selection is performed by operating the CH1 (4) adjustment knob or via up and down arrow keys (10).

To access the test page press ENTER (7) key. It is also possible, in Test Menu, to select a "priority" test by pressing on softkey F2. Please, follow the displayed instructions to complete the setting.

Note: In the description of the different tests, the following guidelines, while referring in particular to the "Pure Tone Audiometry", are generally applicable to all the audiometry tests.

Pure Tone Audiometry

Selection of the type of signals required and relative transmission to the patient

The selection of the type of signal and the transducers used, is performed for both channels by using the INPUT and OUTPUT (9) keys.

With the F1 (MODE) key, it is possible to activate the pulsed-tone mode. The pulsed-tone pulsing rate can be changed with the F2 (PULSE RATE).

To choose which ear the signals are to be sent to is entered via the RIGHT and LEFT keys of the two channels (8). Press again the selected side button to enter R+L mode, this will provide transmission of the channel signal simultaneously to both ears (R+L). To change the pure tone and NBN frequency use the FREQUENCY

(13) key. The intensity of the signals transmitted to the patient are adjusted for both channels via the adjustment knobs (4).

The patient signal transmission mode is entered by the NORM/REV (11) and SIGNAL (12) keys for both channels: in the NORMAL mode, transmission occurs by pressing on the SIGNAL keys. In the REVERSE mode, signal transmission is on and can be interrupted by pressing the SIGNAL keys.

For high-intensity pure-tones (20 dBHL below maximum level), the channel goes into EXTENDED mode and the tone can be transmitted to the patient only for a duration of one second. In the EXTENDED mode, the NORM/ REV (11) key is no longer active.

Signal intensity and channel synchronization

To set the LOCKED and TRACKING function: PARAMETERS (5) softkey enables access. Pressing the softkeys will activate the respective functions.

If LOCKED is activated, it enables control of the CH2 stimulus emission from CH1 (the CH2 NORM/REV and SIGNAL keys are deactivated).

If TRACKING is activated, it enables control of the CH2 intensity levels from CH1: by acting on the CH1 dB-adjustment knob (4), the CH2 signal emission intensity is adjusted by the same value. The CH2 adjustment knob (4) is active.

Storing values on an audiogram

To store an audiogram plot occurs by pressing on STORE (17). The value-plot cannot be stored on the audiogram if the R+L function is activated (i.e. if signal transmission is provided simultaneously to both ears).

Deleting curves and generated plots

To delete plots recorded on an audiogram, select the frequency and the level where the plot is generated then press the DELETE PLOT (5) softkey that will show on the display (1).

To delete curves previously generated on an audiogram, press the TEST MANAGEMENT (5) softkey.

Another two softkeys will be displayed on screen, i.e. DELETE CURVE and DELETE ALL (5). To delete one single curve, select output (9) and the corresponding side (8), then press the CANCEL CURVE (5) softkey.

To cancel all the curves generated and displayed (1), press the CANCEL ALL (5) softkey.

Test report storage

All the curves generated by the different tests are stored and remain visible on display until the unit is powered OFF. Powering the device OFF will result in cancellation of the tests performed. Therefore, prior to powering the unit off, either download the acquired test reports onto a PC or print out the test reports via the built-in thermal printer (2) or through a .pdf report creation.

Storage of test parameter settings

Press the TEST MANAGEMENT softkey to access the PARAMETER SETTINGS (5) softkey. Pressing this key will store the input, output, frequency side and levels being used for the test that is currently underway. The same parameter settings will be automatically maintained each time the same test is accessed.

Storage shall be subject to acceptance by the operator, to be confirmed via the ENTER key (7).

Patient-Operator Communication

The operator can communicate with the patient either via a built-in microphone (3) or via a headset integrating the operator microphone, which is connected up to the appropriate connector located on the unit's rear panel (26). Once the headset microphone is connected via the rear panel, the built-in microphone is excluded. The patient will be able to hear the operator via the transducers used for test performance (headphones, bone conductor, speakers). For operator-to-patient communication, press on the T. OVER (15) key, for level adjustment, use the CH1 or CH2 (4) adjustment knobs.

The patient speaks to the operator over the patient microphone that is connected to the unit via the rear panel (25). The operator hears the patient either via the built-in monitor speakers or via monitor-headset, which is connected to the appropriate connector located on the unit's rear panel (27). Once the monitor-headset is connected via the rear panel, the built-in speaker is excluded.

The operator microphone, the patient microphone and the monitor-headset volume levels can be adjusted by selecting the CONTROLS (14) key which enables access to the controls page. Selection and adjustment is performed via the arrow keys (10) or the adjustment knobs (4) for both channels.

Monitor Function

The operator can listen to the signals transmitted to the patient either via the monitor-speakers built into the unit or via a monitor-headset connected to the unit's rear panel. The monitor is activated or deactivated via the controls page, by selecting the MONITOR softkey (F4). The control page is accessed by pressing the CONTROLS (14) key.

Speech Audiometry

Speech material threshold regulation

The test can be performed using speech material from the following sources:

A CD-player connected to the unit (EXT1 and EXT2 inputs); an operator microphone connected to the appropriate connector on the unit's rear panel, SD card storing speech material supplied with the unit and connected into the appropriate slot on the unit's rear panel (30).

For adjustment of the external signal levels (EXT1 and EXT2) coming either from the CD-player or from the operator microphone: Connect the external signal source to the CD TAPE connectors on the rear panel (24). From the test page, switch both channels (11) to ON.

Via the INPUT (9) keys, select EXT1 and EXT2.

Press the CONTROLS (14) key to access the controls page.

Proceed with playing the CD track with the calibration signal.

Using the adjustment knobs (4) for both channels, adjust the level of the two EXTERNAL GAIN bars until the two bar graphs alongside the two adjustment bars, read 0 dB.

Press the SAVE (5) softkey to store the level then press GO BACK (5), to return to the test page.

Level adjustment of the speech material provided on the SD card does not require any further adjustments by the operator. By the keys INPUT (9) and SIDE (8) of the two channels, it is possible to send the masking noise and the speech material in the same ear in order to perform a Speech in Noise test.

Test conduction and percentage answers

Using the operator microphone or an external signal source (via MIKE or the EXT1 and 2 inputs), the operator first needs to set the number of words to be transmitted to the patient, by acting on left and right (10) arrow keys. Once the phoneme has been conveyed to the patient, the correct reply shall be transmitted to the program by pressing on the up (10) arrow key. The percentage is detected automatically. The relative audiogram plotting is generated by pressing the STORE (17) key.

With integrated SD card speech material (SD CARD input), the available speech material play lists and speech material genres are shown on display. Using the arrow keys (10) it is possible to select the speech material genre (e.g. two-syllable words for Italian), the language and the list. To start the test, press the PLAY (5) softkey. During the test, correct replies can be entered either using the arrow keys or by pressing on the ENTER (5) softkey. If this last option is used, an indicator will also be displayed alongside the corresponding word in the list being run. While the test is being run, it is also possible:

- to repeat transmission of a phoneme that has just been transmitted (via the REPEAT softkey)
- to go into test pause mode (via the PAUSE softkey),
- to change the time interval between phonemes (via the TIME+ and TIME- softkey)
- to interrupt the test by pressing on the STOP (5) softkey.

The phonemes conveyed to the patient shall be shown on display individually. It is also possible for the operator to listen to them either via the monitor-headset or via the built-in monitor-speakers. Pressing F2 softkey is possible to store the SRT value for the active output.

Auto-threshold Test (Hugson-Westlake Auto-threshold)

The automatic threshold test is started either by pressing on the START/STOP (16) key or by pressing on the START (5) softkey. To interrupt the autothreshold test being run, either press on the START/STOP (16) key or press on the STOP (5) softkey. Threshold readings for each frequency shall depend on the response given by the patient via the patient response pushbutton that is connected via the connector (28) on the unit's back panel. The patient will have to keep the button pressed when he hears the stimulus and release it when he no longer hears the stimulus.

Test ABLB (Alternated Binaural Loudness Balance or Fowler Test)

The test is carried out only using AC headsets alternating the same acoustic stimulus (only in pure tones) between the two ears (right and left), with the possibility of varying independently the level of the stimulus

sent to the two ears. It is possible to change the presentation rate of the stimulus by pressing F2. It is possible to carry out the test for four different stimulation frequencies for each ear, each one will occupy a "Frame" on the exam window. To change the working frame after changing the frequency, press F1.

Test MLB (Alternated Monoaural Loudness Balance or Reger Test)

The test is carried out only using AC headsets sending to a single ear (right or left) alternately two acoustic stimuli (only in pure tones) at two different frequencies, with the possibility of varying independently the level of the two stimuli at a different frequency. To vary the tone of the reference frequency use FREQUENCY (13) keys; to vary the tone of the variable frequency use the left and right arrowhead buttons (10). Frame and stimuli presentation rate management is done as seen before for the ABLB test.

Multifrequency Audiometry

This test is useful if the audiogram shows an hearing loss at a specific frequency. In this case it is possible to investigate this specific frequency with intervals ranging from 1 Hz to 100 Hz. Concerning the commands, the test is similar to the tonal audiometry test seen before (INPUT, OUTPUT, SIDE, etc...). To change the basal frequency with normalized steps, also used in tonal audiometry, use the FREQUENCY keys (13). To vary the frequency around the basal frequency in preset increasing steps, use the arrowhead buttons left/right (10). The frequency increase steps can be set using the up/down buttons (10).

Test DLI (Difference Limen for Intensity) or Luscher Test

DLI allows evaluating the subject's ability to detect small short intensity increments. The subject receives a pure tone at a fixed intensity and frequency modulated with successive increases of intensity varying from 0.1 dB to 5 dB with a presentation frequency varying from 0.5 Hz to 2 Hz. To vary the intensity of the increments press F1. To vary the presentation frequency of the increments press F2.

SISI test (Short Increment Sensitivity Index or Jerger Test)

A test similar to DLI in which the subject receives a pure tone at a fixed intensity and frequency, modulated with successive increases of intensity variable from 0.25 dB to 5 dB with a presentation frequency varying from 0.2 Hz to 0.5 Hz or random with fixed durations of 200 msec. The subject receives 20 stimulus increases for each frequency and is required to give an answer when he/she can identify an increase. To start the test in automatic mode, press START/STOP (16), verify that the increment value is set on 1 db or less. If this value is higher only manual mode will be available for patient training purpose. A percentage of the correct given answers will be available at the end of the sequence (SISI Score). To vary the increment size press F1; to vary the frequency rate press F2.

Tone Decay test

This test evaluates the ability of the subject to perceive a stimulus for a determined length of time. The subject receives an acoustic stimulus at a certain intensity and frequency for 60 sec. The subject must give a response until he/she can hear the signal. It is possible to increase the intensity of the signal, with 5 dB steps, from the moment the subject cannot hear it anymore. To increase the intensity of the stimulus during the test use the channel 1 knob (4).

Stenger

With this test it is possible to send the same type of stimulus to both ears simultaneously. This test is used to evaluate the masking effect of a signal sent to an ear when the opposite ear receives the same signal but at a different intensity. In addition to the Pure Tones, it is possible to use the same input signals present in speech audiometry (speech material from external players, from the operator's microphone, from the SD card present in the device). For this test it is not possible to visualize audiograms that were previously stored.

Bekesy Automatic Audiometry Test

Test used to diagnose problems relative to the acoustic nerve.

The test is conducted in automatic mode with the subject's collaboration. The subject has to give a

response only when he/she hears the output stimulus. As a stimulus, pure tones at different frequencies sent in continuous, pulsed (presentation frequency 2 Hz), lot (presentation frequency 1 sec mode, ON=200 msec, OFF=800 sec) are used. The increment of the output stimulus will be continuous and with variable speed (1 dB/sec, 2.5 dB/sec, 5 dB/sec). The frequency variation can be fixed (frequency steps in octaves), or continuous (frequency sweep). It will also be possible to select the presentation sense of the stimulus: from high to low frequency or vice versa. The overall exam time can be between 30 to 60 sec. for the fixed presentation mode and 60 sec. for the sweep mode.

Once inside the test page, use the following commands for the different settings:

For stimulus presentation mode (continuous, pulsed, lot) use F1.

To increment speed of the stimulus press F4>F5 then F5 for variations.

For the frequency variation (fixed or sweep) press F4>F5 then F4 to choose.

For the presentation sense (low>high or high>low) press F4 then F3 to choose.

For test duration (30sec or 60 sec) press F4>F5 then F6 to choose.

Test GAP (Gap Detection Threshold)

Gap detection testing can be done several ways, but the premise is to introduce two stimuli (tones, clicks, broad band noise) of some duration with a variable gap between them.

The subject is to identify when the temporal gap when the two stimuli are heard as one sound or two. The threshold at which the gap is sufficiently wide to be heard as two sounds is the "gap detection threshold." Research on speech perception indicates that persons should be able to hear transitions of formant frequencies of speech in the 0 to 20 millisecond range, so that the assumption is that gap detection thresholds greater than 20 milliseconds are abnormal and a temporal processing disorder may be present. The greater the gap detection threshold, e.g. up to 300 msec, the greater the probability that a person will have difficulty with speech discrimination.

To change the GAP use F1 softkey or left and right (10) arrow keys.

To change the Rate use F2 softkey.

To perform the test set CH1 stimulus (WN or SN) at a comfortable level and press START/STOP (16). Patient should press the patient push button if he can recognize the GAP. At the end of the test a score will be automatically calculated based on heard stimuli percentage; the value can be stored by pressing STORE (17) key.

Test DLF (Difference Limen for Frequency)

The minimum detectable change, or difference limen, of frequency of pure tone (DLF) is a fundamental limit of frequency resolution and is, hence, an important criterion to examine general theories of hearing. There have been two common ways of measuring frequency discrimination. One involves the presentation of two successive steady tones with slightly different frequencies. The subject is asked to judge whether the first or the second had the higher frequency. The order of the tones is varied randomly from trial to trial, and the DLF is usually taken as that frequency separation between the pulses for which the subject achieves a certain percentage of correct responses.

A second measure uses tones which are frequency modulated (FM) at a low rate (For example: 5Hz modulation frequency with a 500-Hz carrier). Usually, two successive tones are presented, one modulated and the other unmodulated. The amount (or depth) of modulation necessary for detection is determined.

To set the stimuli frequency variation use F1 softkey.

To set stimuli presentation rate use F2 softkey.

To perform the test set CH1 stimulus at a comfortable level and increase the DLF percentage until the patient can perceive the difference between the two frequencies. The value can be stored by pressing STORE (17) key.

Quick SIN Test

Quick SIN is a speech-in-noise test used to quickly evaluate patient capability to understand speech in noisy environment. This capability can not be determined through the standard speech audiometry

tests. Quick SIN should be performed on adults in addition to standard audiometric tests. It is available only in English language and requires an optional licence to be activated.

The test is performed by submitting 6 sentences, each of them composed of 5 keywords, at a comfortable level, while the background noise is increased. A score is calculated based on recognized keywords, the total score assesses the patient understanding capability. Speech material is provided on SD-Card, for an easy test execution.

Loading paper in the thermal printer


Lever printer locked in print position (down).

1. Turn on the Resonance device.
2. Open the printer lid. With the printer lever locked in print position (down), put the paper into the roller as shown in Picture 1, and then press the F6 key to advance the paper.

3. Once the paper is through to the other side of the roller, unlock the printer lever and manually position the paper so that the edges of the two sides are aligned as shown in Picture 2. Once you have the alignment, lock again the printer lever in print position (down).

4. Lower a little the printer lid, Insert the paper into the proper slot on the lid and close it.

Printer lever unlocked (up)

Print management and .pdf report

After performing the test, a summary can be printed via the built-in thermal printer, by generating a report in pdf sent directly to a USB stick connected to the instrument, or by downloading to the software manager. An external USB keyboard to edit data directly from the device is available as optional.

Print Summary by thermal printer

When on the page of any of the test, press the TEST button (6) to move in one of two pages SELECTION TEST. Press PRINT key (18) to access the management page printing.

Using the up and down arrow keys to highlight texts you do not want to print and deselect them by pressing ENTER (7). Always using the up and down arrow keys move through the fields (editable via the external USB keyboard). Once you have highlighted the scope for the type of release (pdf or thermal), using the left and right arrow keys select "Print Thermal" and press the PRINT key (18). The test will be printed in sequence the thermal printer.

Print summary reports via pdf

Repeat the above operations above until the selection of the type of printing. Select the item "Print PDF" and press the PRINT key (18). A pdf file will be saved in the root of the USB stick and appointed with the first three characters of the name followed by the patient's current date and time (if no patient name was added to the file will be named with "unknown" followed by the current date and time).

Custom logo and header


A custom logo can appear on the report pdf, from the root of the USB stick, there must be a folder named "images" within an image file with the desired logo. The image file must have a .jpg extension, be named **Company_logo.jpg** and be 200 x 50 pixels in size.

In addition to the logo, so that it can appear on the report pdf a custom text with user data of the instrument, in the same images folder, located on the USB stick, there must be a file text called **Company_data.ini** showing the desired text (4 lines of text). Once the USB stick is connected to the instrument the image and the text will be imported directly into the .pdf template.


Main labels

Product identification label

(located on the underside)

Label located in the power supply section

Symbols present on the device

	ON (mains power supply)
	OFF (mains power supply)
	Alternating current
	Please consult the operating instructions
	Application of type B components
	Please consult the manual for instructions on which transducers to select for connection.
	This symbol placed on the device indicates that when it is scrapped, the device must be subject to "separate municipal collection". The User is therefore required to deliver the scrapped device (either directly or indirectly) to the separate waste collection canterers set up for the purpose by the local or municipal waste management organizations, or to give the device back to the dealer against purchase of a new, equivalent type of device.
	RoHS mark ensuring compliance to RoHS Directive
	Together with the relative name and address, this symbol identifies the device Manufacturer.
	CE mark ensuring compliance to the MDD 93/42/EEC Directive - Notified certification body no. 0051 (IMQ).
	Year of production of the device
	USB Host port
	USB Slave port
	Serial number

Messages that may show up

Message on display	When it appears	Message type	Meaning	Troubleshooting
DO YOU CONFIRM ?	This query appears after entry of plot cancellation or parameter change requests	System message	It means that a plot cancellation or a parameter change request has been entered	Press either the YES or the NO when the query appears
CHANGE CALIBRATION DATE ?	This message appears at the end of the device calibration procedure, once the "SAVE" softkey has been pressed and values have been changed.	System message	It provides the option of updating periodic calibration dates	Press either the YES or the NO when the query appears
ENTER PASSWORD	This message appears when attempting first-time access to the device calibration procedure, once the unit has been switched on	System message	It instructs you to enter your access password to unlock the device calibration procedure	Enter your password via the softkeys
WRONG PASSWORD	This message appears when attempting to access the device calibration procedure using a wrong password	System message	It means that you have entered the wrong password for access to the Calibration menus	Press the Enter key, then re-enter the password correctly
PARAMETERS CHANGED	This message is displayed upon startup of the device in the event that the system detects possible errors in one or more than one parameter relating to the test. The incorrect parameters will be displayed simultaneously with the message, on the same page	System message	It means that one or more than one of the files relating to the test parameters is either corrupt or missing	Press the Enter key. The corrupt files will be replaced with the system default files and it will then be possible to proceed using the device. If you are instructed to do so, it will be necessary to run the device calibration procedure

Maintenance and cleaning

- Clean the outer shell of the device using a clean, damp cloth (do not use aggressive detergent solutions). The penetration of liquids must be strictly avoided.
- During testing operations, headset pads and probe tips come into direct contact with the patient. Although said components are biocompatible, they are not sterile and they must therefore be subject to special care:
 - After each test session it is necessary to disinfect the headset pads before using them to test a new patient. To do so, use only commonly marketed non-allergenic disinfectant solutions and follow the directions supplied by the solution manufacturer.
 - The transducers and insertion headsets tips (IP30 and Insert) are disposable: use only once and substitute them before use with a new patient.
- Calibration and any possible repair, installation and/or update of the device must be strictly carried out only by specialized technical personnel, duly authorized by the manufacturer and in full compliance with the terms and conditions specified in the Technical Manual.

- Standard maintenance and device calibration operations are to be regularly performed at one-year intervals. Standard maintenance operations cover the inspection of the transducer cable conditions, of the mains power supply cable, of the earphone pads as well as an inspection check to ensure the device's outer section and/or cabinet are in good working order. All device calibration operations must be strictly performed only by qualified technical personnel, duly authorized by the manufacturer and implementing the appropriate apparatus for checking the levels and frequencies of all the output signals by the transducers supplied on-issue with the device. Once device calibration testing using the special purpose apparatus is complete, it is necessary to run a device function test.
- The R27A/R37A device can be used combined with a soundproof or silent cabin as per the ambient condition requirements suitable for audiometric testing. Prior to connecting the device, check to ensure that the sockets are compatible with the device features.
- Device calibration is warranted by the manufacturer only if the transducers implemented are directly connected to the device. Should transducers be connected to the device through the inter-positioning or insertion of extension cables and/or interconnection sockets, as may be the case in the event of connections to soundproof cabins, it is necessary to first proceed with unit recalibration prior to using the device.
- In the event of replacement of transducers and/or other device components, exclusive use of Resonance Original spares supplied directly by the manufacturer is highly recommended. Whenever transducers are replaced, recalibration of the device by qualified technical personnel is necessary.
- The device is provided with a lithium battery type CR2032 which enables constant updating of the date and time data displayed by the device. Said battery is welded onto the device PCB and possible replacement thereof must strictly be performed only by qualified technical personnel. Failure to observe this instruction may result in damages during battery replacement, thus preventing correct display of the date and time data by the device.
- Blown fuses must be replaced only with the same type of fuses, as specified on the device itself.

Maintenance and cleaning schedule

Daily	Weekly
<ul style="list-style-type: none"> • Presence of signals in the transducers provided • Cushions cleanness 	<ul style="list-style-type: none"> • Good condition of all accessories cables and related jack-plugs
Monthly	Annually
<ul style="list-style-type: none"> • General device condition • Cushions condition 	<ul style="list-style-type: none"> • Perform calibration and routine maintenance check

Environmental conditions

Storage and handling	temperature range	from -20° C to +50° C / from -4°F to +122°F
	humidity range	from 20% to 80%, without condensation
	pressure range	from 500 hPa to 1060 hPa
Operating conditions	temperature range	from +15° C to +35° C / from +59°F to +95°F
	humidity range	from 30% to 90%, without condensation
	pressure range	from 700 hPa to 1060 hPa

Main technical specifications

Dimensions

- L x W x H: 370x290x180 mm Weight: 3.5 kg net

Power Supply

- 110-240 V AC 50/60 Hz 40 VA Fuses: 2 x T 1 A L 250 V

Reference Standards

- MDD 93/42/EEC and update 2007/47/CEE: Class IIa (as referred to in Annex IX, rule 10 of said MDD 93/42 EEC)
- Safety: to IEC 60601-1:2005 + AMD1:2012 Class 1 Type B
- EMC: to EN 60601-1-2 (2012)
- Audiometer: to EN 60645-1; EN 60645-2 and ANSI S3.6 Type 1A
- Manufacturer quality system: to EN ISO 9001:2015 and UNI EN ISO 13485:2016 (Notified certification body: IMQ)

Tests

- Pure tone Audiometry and Speech Audiometry (speech material recorded on SD card)
- ABLB; MLB; DLI; SJSI
- Tone Decay, Multifrequency, Stenger
- Automatic threshold (Autothreshold) Audiometry, Békésy
- GAP, DLF, Quick SIN (optional)

Outputs

- ACR, ACL, BC, INSERT, Free Field, monitor phone

Inputs

- CD tape (Ext1 and Ext2), mike talk over

Type of signals

- pure tone, warble, mike, external, speech material recorded on SD card

Masking

- Narrow band noise, speech noise, white noise, external signal.

Max frequencies and max levels

- Tone frequency accuracy: ±1%
- dB level variation: 1, 2, 5 dB

freq. (Hz)	DD45-TDH39		HDA280		HDA300		INSERT	BC B71 (dBuN)		Free Field	
	Tones (dBHL)	Noise bands (dBEM)	Tones (dBHL)	Noise bands (dBEM)	Tones	Noise bands (dBEM)	Noise bands (dBEM)	Tones	Noise bands	Tones (dBHL)	Noise bands (dBEM)
125	80	55	80	55	----	----	----	----	----	70	70
250	100	75	100	75	----	----	75	45	40	85	85
500	120	90	120	90	----	----	95	70	60	95	95
750	120	95	120	95	----	----	95	75	65	95	95
1000	120	95	120	95	----	----	95	80	70	95	95
1500	120	95	120	95	----	----	95	80	75	95	95
2000	120	95	120	95	----	----	95	80	75	95	95
3000	120	95	120	95	----	----	95	80	70	95	95
4000	120	95	120	95	----	----	95	80	70	95	95
6000	115	90	110	90	----	----	80	60	55	95	95
8000	100	70	100	70	90	80	----	----	----	85	85
9000	----	----	----	----	90	80	----	----	----	90	80
10000	----	----	90	70	90	80	----	----	----	90	80
11200	----	----	----	----	90	80	----	----	----	90	80
12500	----	----	80	70	80	70	----	----	----	90	80
14000	----	----	----	----	70	60	----	----	----	85	75
16000	----	----	----	----	50	40	----	----	----	55	45
18000	----	----	----	----	40	30	----	----	----	----	----
20000	----	----	----	----	30	20	----	----	----	----	----
WN	110 (dBspl)		110 (dBspl)		----		110 (dBspl)	100		95 (dBspl)	
SN	90		90		----		90	60		95	
EXT	105		105		----		----	60		90	

DLI and SISI Tests

Increment levels SISI: 0.25, 0.5, 0.75, 1, 1.5, 2, 3, 4, 5 dB

SISI increment recurrence rates: 0.2, 0.5 Hz, random

Increment levels DLI: 0.1, 0.2, 0.3, 0.4, 0.5, 0.6, 0.7, 0.8, 0.9, 1, 1.5, 2, 3, 4, 5 dB

DLI increment recurrence rates: 0.5, 1, 2 Hz

Tone Decay Test

Test duration: 60 sec.

Multifrequency Audiometry Test

Frequency variation step: 1, 10, 20, 50, 100 Hz.

Bekesy Test

Frequencies: from 125Hz to 8 KHz in one octave step

Stimuli presentation: steady, pulsed (2 Hz), LOT (rate 1Hz – 200 msec ON / 800 msec OFF)

Direction: low to high frequency / high to low frequency

Intensity changes: 1, 2.5, 5 dB/sec.

Mode: Sweep (60 sec. duration), Fixed (30 or 60 sec. duration)

GAP Test

GAP duration: 1, 2, 5, 7.5, 10, 15, 20, 35, 50, 75, 100, 200, 500 ms

Repetition rate: 0.2, 0.5 Hz, random

DLF Test

Frequency variation: 0, 0.25%, 0.5%, 1%, 2%, 4%, 6%, 8%, 10%, 15%, 20%, 30%

Repetition rate: 0.2, 0.5, 1Hz

Technical descriptions more detailed and complete as wiring diagrams, parts list, instructions for calibration and other materials necessary for the repair and maintenance of the instrument, are contained in the service manual which will be provided, upon request, technical personnel authorized by the manufacturer to carry out such operations.

Models

SERIES	MODEL	CODE
R27A	R27A TDH Diagnostic	4300102500
	R27A DD45 Diagnostic	4300102505
	R27A DD45 Diagnostic	4300102510
	R27A ADC Diagnostic	4300102520
	R27A ADC Diagnostic	4300102525
	R27A HDA Diagnostic	4300102540
	R27A HDA Diagnostic	4300102550
R37A	R37A TDH Clinical	4300102600
	R37A DD45 Clinical	4300102605
	R37A ADC Clinical	4300102620
	R37A HDA Clinical	4300102640
	R37A HF Clinical	4300102660
	R37A HF Clinical	4300102665

Accessories, test and applied parts

● = standard ○ = optional — = not available * = applied parts	FAMILY		
	R27A	R37A	R37A-HF
ACCESSORIES			
Adult easy adjust headset DD45 or TDH39 *	●	●	●
Insert earphones *	—	—	○
HDA300 HF headset	—	—	●
HDA 280 headset (depending on ordered model) *	●	●	●
ADC audiocups (depending on ordered model) *	●	●	●
Bone Vibrator B71W *	●	●	●
Insert for masking *	○	●	●
Operator Headset (headphone + micro) *	○	●	●
Patient Microphone *	○ 4300102510	●	●
Goose-neck operator microphone	○	○	○
Patient push-button *	●	●	●
Thermal printer	● - 4300102510	●	●
Thermal paper roll	● - 4300102510	●	●
SD Card with Speech material	● - 4300102510	●	●
Power supply cable	●	●	●
Device dust cover	●	●	●
Multilingual Quick user's handbook	●	●	●
USB Pen Drive	●	●	●
Resonance® MDS software (demo version)	●	●	●
DHR documents regarding calibration data and CE compliance	●	●	●
External USB HUB Keyboard	○	○	○
Carrying bag	○	○	○
Resonance® MDS software License	○	○	○
TEST			
TONE AUDIOMETRY	●	●	●
SPEECH AUDIOMETRY	● - 4300102510	●	●
ABLB	● - 4300102510	●	●
AUTOTHRESHOLD	●	●	●
SISI	●	●	●
MULTIFREQUENCY	○	●	●
BEKESY	—	●	●
MLB	—	●	●
DLI	—	●	●
STENGER	●	●	●
TONE DECAY	—	●	●
GAP	●	●	●
DLF	●	●	●
QUICK SIN	○ - 4300102510	○	○

Applied parts are constituted by the transducers that convey the test signals. Said transducers are applied so that they go into direct contact with the patient. In the foregoing list, they are highlighted in bold script and are counter-marked with an asterisk (*).

Following is a display of the transducers that constitute the applied parts of the device, with a description of the name and number of the relative connectors located on the unit's rear section for connection of the individual transducers.

Transducers for audiometric tests

Bone vibrator
"Bone vibrator" (22)

Insert for masking
"Insert" (21)

HDA 280 headset
"AC Phones R-L" (19)

Adult easy adjust headset
"AC Phones R-L" (19)

Insert earphones
"AC Phones R-L" (19)

ADC audiocups
"AC Phones R-L" (19)

Patient Microphone
"Talk Back" (25)

Patient push-button
"Pat. Resp" (28)

Electromagnetic emission

Guidance and manufacturer's declaration - Electromagnetic EMISSION

The R27A/R37A audiometer is intended for use in the electromagnetic environment specified below. The customer or the user of the R27A/R37A should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment - guidance
RF Emissions CISPR11	Group 1	The R27A/R37A uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR11	Class B	The R27A/R37A is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies building used for domestic purposes.
Harmonic Emissions IEC/EN 61000-3-2	Complies	
Voltage fluctuations/flicker emission IEC 61000-3-3	Complies	

Producer's guide and declaration - Electromagnetic IMMUNITY

The R27A/R37A is suitable for use in the specified electromagnetic environment. The purchaser or user of the The R27A/R37A should assure that it is used in an electromagnetic environment as described below:

Immunity test	IEC 60601-1-2 Test level	Compliance Level	Electromagnetic Environment
Electrostatic discharge (ESD) IEC 61000-4-2	8 kV contact 2/4/8/15 kV air	IEC 60601-1-2 Test level	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated electromagnetic field IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	IEC 60601-1-2 Test level	Portable and mobile RF communications equipment should be used no closer to any part of the R27A/R37A including cables. Recommended distance: $d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum power output of the transmitter in watt (W) based on the transmitter's producer and d is the recommended distance in meters (m).
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines > 3 m	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	0.5/1 kV differential mode 0.5/1/2 kV common mode	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V 150 kHz to 80 MHz	IEC 60601-1-2 Test level	Portable and mobile RF communications equipment should be used no closer to any part of the R27A/ R37A including cables. Recommended distance $d = 1,2\sqrt{P}$

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% Un for 0.5 cycle 0 % Un for 1 cycle 70 % Un for 25 cycles 0 % Un for 5 s	IEC 60601-1-2 Test level	Mains power quality should be that of a typical commercial or hospital environment. If the user of the R27A/R37A requires continued operation during power mains interruptions, it is recommended that the R27A/R37A be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	IEC 60601-1-2 Test level	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: At 80 MHz and 800 MHz, the higher frequency range applies.
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the R27A/R37A is used exceeds the applicable RF compliance level above, the R27A/R37A should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the R27A/R37A.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distance between portable and mobile
RF communications equipment and the R27A/R37A**

The R27A/R37A audiometer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the R27A/R37A can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the R27A/R37A as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 KHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.2
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

this page intentionally left blank

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