



TRUMPET

REAL EAR MEASUREMENT SYSTEM AUDIOMETER

USER MANUAL



Read this manual carefully before using the instrument. Pay particular attention to the instructions given in Chapter 1 ("Safety: warnings and information") and in Chapter 3 ("Installation").



Internal inspection and servicing of the instrument must be entrusted entirely to technicians approved by Inventis S.r.l.

Copyright: Inventis S.r.l. is the owner of the copyright on this manual. It is forbidden to copy, reproduce or alter the manual, in its entirety or in part, without the express written authorization of Inventis S.r.l..

Inventis ® is a registered trademark owned by Inventis S.r.l..

(€ 0123



Inventis S.r.l. Corso Stati Uniti, 1/3 35127 Padova, Italia Tel.: 049.8962844 Fax: 049.8966343 www.inventis.it info@inventis.it

Revision:	
Date:	

22 2022.05.26

Summary

Summa	<i>ryv</i>
CHAPT	TER 1 Safety: warnings and information
1.1	User Manual1
1.2	Operator responsibilities1
1.3	Intended use of the device2
1.4	Indication for use and end users of the device2
1.5	Precautions
1.6	Disposal7
1.7	Conformity7
1.8	FCC/ISED use restrictions and warnings7
1.9	Symbols on labels10
CHAPT	TER 2 Introduction
2.1	Description of the product11
2.2	Included and optional parts13
2.3	Preparation for an REM test22
CHAPT	TER 3 Installation
3.1	Connections
CHAPT	TER 4 Maintenance 43
4.1	Calibration and equalization (REM) 43
4.2	Periodic checks
4.3	Cleaning and maintenance of transducers

4.4	Replaceable parts	46
4.5	Repairs and technical assistance	46
4.6	Enable new licenses	46
4.7	Firmware upgrade	47
4.8	Calibration of microphones by code	48
APPE	NDIX A: Technical Specifications	49
REM	specifications	49
Audio	meter specifications	51
Other	specifications	57
APPE 	NDIX B: Specifications of speaker for free field te	sts 63
APPE	NDIX C: Troubleshooting	65
APPE	ENDIX D: Electromagnetic emissions	71

CHAPTER 1

Safety: warnings and information

1.1 USER MANUAL

Be sure to read this manual through completely, so that all of the features offered by the Trumpet instrument can be used to their full potential. In particular, be sure to read this chapter in its entirety, as it contains information and warnings that are of fundamental importance in ensuring safe and correct use of the device.

The safety warning symbol illustrated below is used in this manual to draw the attention of the reader to information of particular importance in matters of safety, and to guard against incorrect use.



1.2 OPERATOR RESPONSIBILITIES

The Trumpet device guarantees consistent and dependable performance only when used in accordance with the instructions and procedures described in this manual.

In the event of the device malfunctioning and needing to undergo repair or maintenance work, it must be disconnected from the electrical power supply, and not used again until after this work has been completed. When defective or failed parts must be renewed, accept only original replacement parts supplied by Inventis s.r.l. Any repairs must be entrusted exclusively to Inventis or to service technicians approved by Inventis.

No parts of the device must be modified or replaced without authorization from Inventis.

The user of the device assumes full responsibility for any malfunction resulting from improper use or operation, likewise from maintenance or repair work performed by third parties other than Inventis s.r.l. or Service Centres approved by Inventis. Inventis s.r.l. and approved Service Centres will answer for the performance and reliability of the equipment only if:

- 1. adjustments, modifications or repairs have been entrusted to persons authorized by Inventis;
- 2. the electrical system and earthing of the installation are in compliance with the specifications of standards for electro-medical appliances.

1.3 INTENDED USE OF THE DEVICE

Trumpet is a real ear measurement system. A real ear measurement system allows the determination of the intensity of sounds reaching the tympanic membrane of the patient. To this end it uses a couple of microphones, one located at the level of the patient's earlobe (the reference microphone), and one inserted in the patient ear canal through a silicon tube (the probe microphone). The operator, using the system, delivers to the patient sound stimuli different in intensity and in frequency content, and measures the intensity of the sound received by the two microphones. A real ear measurement system is generally used to quantify the gain provided by a hearing aid, when this is worn by the patient.

Trumpet can be also an audiometer. An audiometer is a device that helps the operator in defining the patient's auditory sensitivity by generating and delivering to the patient sound stimuli of different types and intensities for diagnostic purposes.

1.4 INDICATION FOR USE AND END USERS OF THE DEVICE

Trumpet is intended for use by healthcare ENT professionals in hospitals, ENT clinics and audiology offices or by hearing care professionals in conducting hearing evaluations, in assisting in diagnosis of possible otologic disorders, and in performing the verification of fitting of a hearing aid.

There is no patient population restriction in the use of the device; always be sure to perform an otoscopy before using the device. These tests must be conducted in a quiet environment to avoid artifacts and to ensure that errors are not committed when determining the hearing threshold.

1.5 PRECAUTIONS

To ensure correct and safe use of the device, the following precautions must be observed.

1.5.1 Installation and general precautions



Make certain that the required ambient conditions are met (during transport, storage and operation).



The Trumpet will not be protected if exposed during use to flammable anaesthetic gases or similar products. Risk of explosion.



Avoid installing and using the device near sources of strong electromagnetic fields: these could interfere with the operation of the equipment.



Use only original parts supplied by Inventis s.r.l., unless specifically indicated otherwise.



Use only the power adapter supplied by Inventis s.r.l.



The Trumpet is a medical device: if connected to a computer (or any external device) located within the "patient area" (as defined in EN 60601-1-1), this likewise must be a medical device, or protected by an isolating transformer, in order to ensure that the combination of computer (external device) + Trumpet is in compliance with EN standard 60601-1-1.



The Trumpet can be used in conjunction with a soundproof booth to conduct tests under optimum acoustic conditions. Before connecting the instrument to a soundproof booth, check that the sockets are compatible with the specifications prescribed for each connector. The Trumpet has been tested under the following ambient conditions:



Temperature between 18 °C and 28 °C Relative humidity between 20% and 80% Atmospheric pressure between 81.3 kPa and 106.3 kPa.



The Trumpet is designed for use after being positioned on a stable work surface (such as a desk) or a suitable supporting structure (such as a stand).



The Trumpet must be installed and operated taking account of the information regarding electromagnetic compatibility (EMC) provided in Appendix D of this manual



The proximity of portable and mobile appliances used for RF communications can affect the operational efficiency of the Trumpet box. Refer to the information regarding electromagnetic compatibility (EMC) provided in Appendix D of this manual.



The wireless REM probe is an electromagnetic field emitter which could interfere with implantable medical devices (pacemakers, ICDs, etc.). Pay particular attention to the use of the wireless REM probe on patients wearing these devices. If it is suspected that the probe is causing any kind of interference, stop the operation immediately.

1.5.2 Calibration



The calibration of the audiometer and the REM probes is dependable only for the transducers supplied with the instrument. If a transducer is replaced, the instrument must be recalibrated.



The calibration of the audiometer and the REM probes is dependable for transducers supplied with the product, if connected directly to the instrument, without any interposition of extension leads and without the passage from connectors to panel (as habitually occurs in soundproof booth installations). If the transducers are not connected directly to the Trumpet, a new calibration procedure will be required before the instrument is used.



Calibration of the output in free field conditions can be considered dependable as long as the acoustic properties of the surroundings are not altered significantly. If this condition cannot be maintained, or the instrument is moved, a fresh calibration procedure will be required.



When conducting audiometry tests, the presence of a transducer that has not be calibrated will be flagged by the software. It will not be possible to present any stimulus to the patient using transducers that have not been calibrated.



Take note of the calibration interval indicated for the audiometer and the REM probes. Use of the instrument beyond the calibration interval expiry date can lead to unreliable diagnoses.

1.5.3 Hygiene and infections



The earpieces of ER-3C or IP30 earphones are disposable; do not use the same earpiece for different patients. Dispose of earpieces after use.



Disinfect the cushions of headphones between one patient and the next, otherwise there is the risk that infections could be passed on.



The silicone tubes connected to the microphone of REM probes are disposable; <u>do not use the same tubes for different patients</u>, as there is the risk that infections could be passed on. Dispose of tubes after use.



Probe tubes can be reused for the <u>same</u> patient: clean using a cloth moistened lightly with alcohol before inserting into the ear. Check also that there are no obstructions internally of the tube: if the bore is not completely clear, use a new tube.



Use only biocompatible silicone tubes supplied by Inventis srl.



An ear inspection must always be performed before inserting the probe tube, to make certain that the auditory canal is free of any obstruction caused by accumulations of wax, and that there is no inflammation or possible infection.

1.5.4 Avoiding discomfort for the patient

The Trumpet is capable of generating pure tones at an intensity potentially harmful to the patient (up to 90-100 dB SPL). Take particular care over the intensity of the tone before presenting it: studies have shown that exposure to a sound reaching the auditory canal at an intensity greater than or equal to 85 dB SPL for more than 7 hours can result in damage to hearing, but if these levels are amplified by a hearing aid, the same kind of damage can occur in as little as 30 seconds.

To avoid any such risks when conducting an REM test, make certain that the thresholds of the hearing aid have been set correctly, if necessary after running the MPO (Maximum Power Output) test using the relative function included in the Maestro software application (refer to *Maestro – Fitting and hearing aid test functionalities – User Manual*). At all events, the software includes an option for checking that UCL (UnComfortable loudness Level) values attributed to the patient are not exceeded for a given length of time, which interrupts the sound output if such a situation should occur (refer to *Maestro – Fitting and hearing aid test functionalities – User Manual*). If the user has an audiometer license, it is possible to determine the UCL thresholds for the patient (refer to *Maestro – Impedance audiometry functionalities - User manual*).

During the course of an audiometry test, the system will prompt for further confirmation before the intensity of the sound emitted can be increased to a level higher than 100 dB HL.

When conducting REM tests at high levels, used to verify the limits of the hearing aid (as in the case of MPO), the duration of the emission must not exceed 15 seconds.

It is important to remain watchful for any adverse reaction of the patient, so that action can be taken swiftly if clear signs of discomfort or pain are shown during the course of the test, either interrupting the signal at source or distancing the transducer from the test area.

1.5.5 Use



The microphones used for talk-over and talk-back functions shall be placed at a distance not greater than 15 cm from the mouth, to reduce the background noise.

1.6 DISPOSAL

Like any other electronic device, the Trumpet contains certain extremely hazardous substances such as cadmium or mercury, albeit in extremely small proportions. If released into the normal refuse collection system without suitable preliminary treatment, these substances can cause serious damage to the environment and to health. At the end of its service life, accordingly, each component of the device must go through a sorted collection process.

This means that the user should deliver (or despatch) waste items to the sorted collection centres set up by local authorities, or alternatively, hand them back to the dealer when purchasing a new appliance of the same or similar type.

Thanks to the sorted collection of waste items and the subsequent processing, recovery and disposal operations they undergo, appliances can be made from recycled materials, and any negative impact of improper waste management on the environment and on health can be suitably limited.

1.7 CONFORMITY

The Inventis Trumpet device is a class II device, in accordance with Annex IX of the medical devices directive 93/42/EEC as amended and supplemented by directive 2007/47/EC.

Inventis s.r.l. is a company certified to ISO 13485.

1.8 FCC/ISED USE RESTRICTIONS AND WARNINGS

With reference to:

FCC ID: 2AVOO-RE1RA IC: 25857-RE1RA HVINs: Trumpet REM & AUD Wireless, Trumpet REM, Trumpet AUD, Trumpet REM & AUD, Trumpet REM Wireless

and

FCC ID: 2AVOO-RE1RAWP IC: 25857-RE1RAWP HVIN: WPROBE Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Responsible party's contact located in the Canada:

Canadian Certification Consulting, Inc. 2210 Horizon Drive, Suite 17 West Kelowna, BC V1Z 3L4 Tel. No.: 1-250-575-1719 E-mail: info@can-cert.com

This device contains license-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's license exempt RSS(s). Operation is subject to the following two conditions: (1) This device may not cause interference. (2) This device must accept any interference, including interference that may cause undesired operation of the device.

L'emetteur/recepteur exempt de licence contenu dans le present appareil est conforme aux CNR d'Innovation, Sciences et Developpement economique

Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisee aux deux conditions suivantes: (1) L'appareil ne doit pas produire de brouillage; (2) L'appareil doit accepter tout brouillage radioelectrique subi, meme si le brouillage est susceptible d'en compromettre le fonctionnement.

ICES-003 Class B Notice -Avis NMB-003 Classe B:

This Class B digital device complies with Canadian ICES-003

Cet appareil numerique classe B est conforme à la norme Canadien NMB-003.

CAN ICES-3(B) /NMB-3(B)

1.9 SYMBOLS ON LABELS



Name and address of manufacturer.

Warning: the use of this instrument calls for certain precautions; to ensure safe use, consult the accompanying documentation.

Read the operator manual.



This symbol indicates that the product is subject to the requirements of Directive 2012/19/EU on waste electrical and electronic equipment (WEEE). In the event of this product being sold and/or scrapped, it must not be disposed of as ordinary household or industrial waste, but collected separately.



Device with applied parts, Type B (EN60601-1)



Mark indicating conformity with Council Directive 93/42/EEC concerning medical devices (as amended and expanded by Directive 2007/47/EC) – Class IIa device, notified body 0123 (TÜV SÜD Product Service GmbH).

IP (Ingress Protection) code: this device is protected against the penetration of objects larger than 12.5 mm; it is not protected against the penetration of liquids.

Rx only *Warning: under US federal law, the device can only be sold through or on the instruction of authorized healthcare professionals.*

MODEL Model of the device

REF Catalogue code

Serial number of the appliance. The serial number is made up of 13 alphanumeric characters indicating the model, the year of manufacture and an incremental number. In particular, the serial number comprises these segments:

- first 5 characters: Inventis product code
- characters 6 and 7: year of manufacture ("16" denotes 2016)
- characters 8-13: incremental number



SN

UDI code



CHAPTER 2

Introduction

2.1 DESCRIPTION OF THE PRODUCT

The Trumpet is an REM device allowing in-situ evaluation of all the parameters required to ensure optimum fitting of a hearing aid. It is characterized by a small and lightweight set of probes, and a powerful speaker mounted to the front panel of the instrument.

The Trumpet is also a diagnostic audiometer that can be used to perform pure tone and speech audiometry by air or bone conduction.

At least one of the two modules (REM or audiometer) is purchased together with the instrument, with the option of integrating both at a later date by purchasing the relative license.

The instrument is controlled by way of the Maestro application, which can be installed only on a Windows PC.

Measurements	Inputs	Outputs
Standard/advanced Real Ear Measurements	Prerecorded files PN, WN, ISTS, ICRA	Built-in speaker
LIVE measurements	Prerecorded files PN, WN, ISTS, ICRA, live speech	Built-in speaker

The following table summarizes the tests that are available with the REM module of the Trumpet:

During an REM test, the instrument can operate only in FF (Free Field) mode, and the power adapter must be connected (see chapter entitled 3 *Installation*).

The following table summarizes the tests that can be conducted in audiometer mode.

Test	Channel	Inputs	Outputs
PURE TONE	Ch. 1	Pure tone, warble, NBN	AC, BC, FF
AUDIOMETRY	Ch. 2	NBN, WN	AC, BC, FF
AUTO- THRESHOLD	Ch. 1	Pure tone	AC
SPEECH	Ch. 1	MIC, (USB1, USB2)	AC, BC, FF
AUDIOMETRY	Ch. 2	SN, WN, MIC (USB1, USB2)	AC, BC, FF
MASTER	Ch. 1	MIC, (USB1, USB2)	AC
HEARING AID	Ch. 2	MIC, (USB1, USB2)	AC
Optional tests			
QuickSIN	Ch. 1	Pure tone, warble	AC, FF
QuickSilv	Ch. 2	NBN	AC, FF

If the power adapter of the Trumpet is disconnected during an audiometry exam, only AC (Air Conduction) and BC (Bone Conduction) modes remain enabled, whereas FF (Free Field) mode is disabled. To enable free field mode, the power adapter must be reconnected (see Chapter 3 *Installation*). If the instrument is powered up only by way of the USB cable, the maximum limits of the enabled outputs will be reduced by 10 dB.

2.2 INCLUDED AND OPTIONAL PARTS

The Trumpet can be supplied with all or only certain of the parts described, depending on the module purchased. If a license is purchased at a later time, then the corresponding parts will also be supplied (those not included in the version of the instrument purchased initially).

REM version	<u>Wired</u>	<u>Wireless</u>
Trumpet unit	•	•
Wired REM probe	•	-
Wireless REM probe	-	•
Set of probe parts (includes probe tubes, insertion guides and spare O-rings)	•	•
USB cable	•	•
Monitor headphones with microphone	•	•
Medical grade power adapter 15V/2A	•	•
Inventis Software Suite	•	•
User manuals	•	•

The tables below list the parts included in the related modules:

Audiometer version

Trumpet unit

Supra-aural headphones (Telephonics TDH-39 or Radioear DD45)

Set of connectors for external speakers

Bone vibrator B71

Patient response button

Monitor headphones with microphone

Clip-on talk-back microphone

Inventis Software Suite

USB cable

Medical grade power adapter 15V/2A

User manuals

The following are optional parts:

Optional parts

Trumpet soft carrying case

Wall hanger for Trumpet

Trumpet stand

Optional parts – Only Trumpet AUD

Insert earphones Etymotic ER-3C

Insert earphones RadioEar IP30

Pair of SIARE Delta 60 passive loudspeakers

LSX additional loudspeaker

REM license for Trumpet AUD – Includes the REM probe

Wireless REM license for Trumpet AUD – Includes the wireless REM probe

Amplivox Audiocups noise excluding enclosures for TDH-39 / DD45 headphones

QuickSIN[™] test license

Wall hanger for Trumpet and AUD transducers

Trumpet AUD silent room kit

Optional parts – Only Trumpet REM

RECD complete kit

Wireless probe upgrade package

 $AUD\ license\ for\ Trumpet\ REM-Includes\ the\ transducer\ set$



Compatibility has been verified only with the transducers and parts supplied by Inventis S.r.l. Do not use other transducers or parts.

2.2.1 Probes (REM)

Two probes are located at the ends of the REM headset. Detailed description below:



- **B.** Probe tube
- C. O-ring for attachment to ear
- **D.** Reference microphone

F. Guide housing for insertion of tube **G.** Probe calibration guides



It is advisable, where possible, to use the REM probes inserted in the headset

Insertion of probe tube

The Trumpet comes with a set of probe tubes (REM license). To insert this part, position the blue part of the tube over the metal cylinder that provides the housing for the tube (A in the illustration) and push in gently until fully home.



Do not apply excessive force when inserting the tube, as this could damage the actual probe microphone.

2.2.2 Wireless probes (REM)

As an alternative to wired probes, wireless probes can be used for REM measurements:



A. Button with status LED

B. Bluetooth transmitter (battery compartment on the back)

To turn on the probes, simply press the transparent button in the center. The status LED will start flashing blue as it looks for the Trumpet, and then will turn orange when it establishes a connection. If the Trumpet is not found (it must be connected to the power supply and a running PC via USB and the Maestro REM module must be in use), the probes will continue to flash blue and then will turn off after one minute.



Do not remove the blue or red shell for any reason as you may irreparably damage the probe headset or wires. To charge the probe battery place them on the support as shown in the following image:



The probe LED will turn green to confirm that they are charging, and then flash orange and green if the probes are connected to the Trumpet or blue and green if they are not. If the LED does not turn green, try to slightly change the position of the probes so that they are aligned with the center of the instrument, or try to remove them from the charge position, wait five seconds, and then put them back in position making sure that they are exactly in the center.

If the probes are off they will be turned on automatically when placed in the charging position.



The Trumpet must be connected to the power supply and a working USB socket in order to recharge the probes. If you want to charge them with your PC turned off, connect the Trumpet to an externally powered USB hub.



ONLY use batteries supplied by Inventis S.r.l.



The use of unsuitable batteries or incorrect use (reverse polarity, damage to the outer casing, etc.) could lead to overheating or explosions.



It is strongly recommended to place the probes in the recharger when not using them in order to avoid reaching an excessively low charge and risk not being able to finish a test.



It is necessary to recharge the probes completely before using them for the first time after removing them from the packaging. Caution: if you skip this step the residual charge indicated may not be correct.

The battery must be inserted as shown below:



The "–" pole of the battery marked "VARTA" must be face down while the battery is being inserted.

If the probes were sold together with the Trumpet then they will already be paired with it, otherwise before using them you need to perform the pairing procedure:

- 1. Connect the Trumpet to the power supply and USB socket of the running PC.
- 2. Start Maestro and enter the REM module, then if necessary set up the wireless connection and enter the license.
- 3. Turn on the probes by pressing the transparent button in the center: the LED will flash blue;
- 4. Pressing the transparent button on the probes, place them in the charging position while continuing to hold the button, then wait 3 seconds and release it.
- 5. After 5-20 seconds the probe LED will flash orange and green, indicating that the pairing procedure was completed correctly and you can use them for the test.



Make sure that during the procedure there is no other Trumpet connected to the power supply and a USB socket nearby, otherwise the probes could pair with the wrong device.



It is necessary to re-pair the probes and your Trumpet even when they are shipped without the main device, to recalibrate the microphones.

If there are problems it is possible to switch off the probes manually by holding the transparent button for about 10 seconds. They also turn off whenever the Trumpet is not used by Maestro and they are not charging.

Color	Meaning
Off	Probes off
Flaghing blue	Probes searching for
	Trumpet, not charging
Flashing blue and	Probes searching for
green	Trumpet, charging
Flashing orange and	Probes connected to
green	Trumpet, charging
Fixed erange	Probes connected to
	Trumpet, not charging

The following table summarizes the possible states of the probe's LED:

When Maestro is closed the probes are disconnected from the Trumpet, so it is normal for them not to be able to connect (orange LED). If the probes do not connect after more than one minute when the Maestro REM module is open and the probe connection mode indicated is wireless (refer to the manual *Maestro – Hearing aid fitting and testing features – User Manual*), try switching them off and on again or pairing them again.

2.2.3 Equalization of active speaker (REM)

During the speaker equalization procedure, the headset with the REM probes must be placed in front of the speaker, at the same distance and in the same position as will be occupied by the patient's head when the test is actually conducted. For best results, it is advisable to fit the headset directly to the patient before running the equalization procedure. The distance of 1 m from the speaker is the optimum distance both for equalization purposes and for conducting the REM test. Make certain that there are no obstructions (e.g. hair) covering the microphones.

For more details, refer to the manual entitled *Maestro – Fitting and hearing aid test functionalities – User Manual.*

2.2.4 Probe tube calibration (REM)

The position for calibrating the probe microphones is shown below (refer to the manual *Maestro - Hearing aid fitting and testing features - User Manual* for more information on calibration). Insert the tip of the probe tube into the appropriate guides, as close as possible to the microphone of reference, and place the headset about 1 m from the active speaker, in front of it.



2.3 PREPARATION FOR AN REM TEST



The REM test must be conducted only by a person specialized in and authorized to perform audiology procedures. Use of the instrument by a person without the necessary expertise could result in damage to the hearing of the patient.



Change the probe for each new patient. Reusing the tube on different people could result in the spread of ear infections.



Do not push the tube too far into the patient's ear: there is the risk of perforating the eardrum. Refer to the ruler marked on the headset to determine the depth of insertion.



Always carry out an ear inspection before proceeding with the REM test, to make certain there are no obstructions in the ear that could affect the results of the examination.

2.3.1 Measurements in the ear

When preparing the patient for an REM test, observe the following procedure:

- 1. Position the patient at around 1 metre from and facing the speaker of the instrument, and in such a way that the ears are placed approximately at the same height as the instrument.
- 2. Position the black ring of the tube so that the placement corresponds to the required insertion distance (refer to the markings of the ruler on the headset as indicated in the image below). The recommended distances are:
 - a. Men: 27 mm
 - b. Women: 27 mm
 - c. Children: 20-25 mm



<u>NOTE:</u> the distances can vary depending on the patient. To avoid problems, in any event, start with an ear inspection: this will give some idea as to the dimensions of the auditory canal.

- 3. Position the probe headset on the patient, adjusting the loops (O-rings) so that the probe is positioned just under the ear lobe.
- 4. Insert the probe tube into the patient's auditory canal, to the point at which the black ring reaches the intertragic notch. If resistance is encountered, do not force the tube further.

To favour the correct insertion of the probe tube and ensure it does not snag on a wall of the auditory canal, there are insertion guides for this same purpose supplied with the instrument. Insert the guide in the relative hole (if necessary, refer to heading 2.2.1 *Probes (REM)* for further information) and position the tube as illustrated below:



In this way it will be possible to ensure the probe tube is kept straight when inserted into the auditory canal.

5. If required for the measurement (i.e. for REAR measurements), place the hearing aid in the patient's ear, making certain that the probe tube remains in the same position and is not occluded by the hearing aid.

The correct position of the headset and the tube is illustrated in the following image:





It is very important that the probe microphones are correctly calibrated every time a new tube is used. Incorrect or failed calibrations can influence the acquisition of curves and lead to false measurements.



It is advisable, where possible, to use the REM probes inserted in the headset

2.3.2 Measurements using 2cc coupler

Should it be impossible to rely on the cooperation of the patient for the entire duration of an REM test, measurements can be acquired by way of a coupler using the RECD unit, or a Drum instrument connected to the same PC as the Trumpet.

RECD unit

Before proceeding with the test, the speaker of the RECD unit must be equalized and the microphones calibrated. The operation requires that the coupler microphone and the reference microphone, suitably aligned one with another, should be positioned on the speaker of the RECD unit as illustrated in the image below.

For further information, refer to the heading entitled *Coupler measurements* in *Maestro – Fitting and hearing aid test functionalities – User manual.*



To enable correct connection of the coupler to the various types of hearing aid available on the market (e.g., BTE, CIC, ITE, ...), different types of adapters, illustrated in the image below, are available.



- Adapter (1) is designed for use with ITE, ITC and CIC hearing aids. Make certain that hearing aid is securely connected by using adhesive insulating putty in this way:



- Adapter (2) allows connection of BTE type hearing aids to coupler using an ear hook. Insert the end of the ear hook inside the provided transparent tube and connect its other end to the adapter.
- Adapter (3) is designed for use with RIC hearing aids. Here again, always use adhesive insulating putty to make certain of a secure connection in this way:



Once the hearing aid has been suitably connected to the coupler, the test can be conducted by positioning the device on the RECD unit (aligned centrally with the loudspeaker if possible) and ensuring that the reference microphone is placed as close as possible to the hearing aid microphone, but without actually touching it. Remember also that it is not advisable to use RECD in the case of open fit hearing aids.

An example of correct positioning is shown in the following illustration.



Before a coupler fitting procedure can be carried out using the RECD unit it is highly recommended to acquire an RECD curve, so that the target curve calculated by the fitting rule considers the differences between the coupler and the real ear.

The acquisition of an RECD curve involves two steps: measurement of a response in the ear and measurement of a response in the coupler. In particular, acquisition of the response measurement in the ear can be accomplished in two ways: using a disposable tip or using the earmold of the patient's own hearing aid. If the first option is selected, the probe tube must first be inserted in the auditory canal of the patient as described in Heading 2.3.1, following which the disposable tip, connected previously to the insert transducer, can be placed in the ear. This particular step must be carried out with special care, to ensure that the probe tube inserted previously is not displaced overmuch. The correct position of the tube and the transducer is illustrated in the following image:



Once the measurement in the ear has been obtained, the response in the coupler can be measured. Electing to use a disposable tip for acquisition of the measurement in the ear, the response in the coupler can be measured only with the adapter (1): this situation is identified as the HA1 configuration. Accordingly, connect the adapter (1) to the coupler, and the insert transducer to the tip. Finally, attach the tip to the adapter, if necessary with the aid of adhesive insulating putty to ensure the two parts are securely connected (the insulating putty need not be used if the tip covers the hole sufficiently, without leaving any gaps). The correct position is illustrated in the following image:



HA1-tip configuration

The response in the auditory canal can also be measured using the earmold of the patient's hearing aid, rather than a disposable tip. Likewise in this instance, the probe tube must first be inserted in the auditory canal of the patient as described in Heading 2.3.1, following which the mold, connected previously to the insert transducer, can be placed in the ear. Here again, special care must be taken during the operation to avoid displacing the probe tube inserted previously from its correct position. An example of the final set-up is illustrated in the following image:


Electing to use the patient's earmold for acquisition of the measurement in the ear, there are two distinct configurations possible for measuring the response in the coupler.

The first configuration, identified as **HA1**, requires the use of the adapter (1). As with the disposable tip, connect the adapter (1) to the coupler, and the earmold to the insert transducer. Finally, attach the earmold to the adapter, using adhesive insulating putty to ensure that the two parts are securely connected with no gaps left between them. The correct position is shown in the following image:



HA1-mold configuration

The second configuration, on the other hand, identified as **HA2**, requires the connection of the adapter (2) to the coupler. In this instance, the only action required is to connect the outlet of the adapter (2) to the insert transducer using the 25 mm silicone tube, as shown in the following figure:



HA2-mold configuration

When both measurements have been taken, in the ear and in the coupler, the difference between them allows automatic calculation of the RECD curve, which will be used during the course of the coupler mode fitting procedure. For further details on using the Maestro software for acquisition of the RECD curve, refer to the heading entitled *RECD* in the Operator manual *Maestro – Fitting and hearing aid test functionalities – User Manual*.

Drum

To enable measurements with the coupler using the Drum instrument, the same considerations apply as for the RECD unit, since the parts supplied with the product are the same (except for the 25 mm silicone tube, which is bent to a right angle rather than straight, but entirely similar in terms of use). The only action required is to connect the Drum to the same PC as the Trumpet.

For an overview of how the hearing aid is positioned correctly in the Drum instrument, refer to the Heading *Position the hearing aid in the test chamber* of the relative manual: *Drum – User Manual*.

Before performing an RECD measurement with the Drum, the software will check that the alignment of the microphones is still as required; if not, refer to the section on *Calibration and alignment of microphones* in the relative manual: Drum - User Manual.

CHAPTER 3

Installation

Whilst the installation of the Trumpet device is a relatively simple procedure, it should be entrusted to a person with the requisite skills. If the installation is not performed correctly, the system could be affected by safety problems when in use.

This chapter describes the procedure for installing the system.



Keep the packaging materials, in case the instrument should need to be sent to the dealer or to Inventis for any reason.

3.1 CONNECTIONS

3.1.1 Back Panel

Illustrated below are the connection sockets on the back panel of the instrument.



Plug the transducers and parts into the respective sockets as indicated in the following table:

Connector	Part
PAT. RESP	Patient response button
BC	Bone vibrator
AC L and AC R	Supra-aural headphones (TDH-39 or DD45) for Air Conduction: right (AC R) and left (AC L)
TALK BACK	Microphone used by patient
TALK OVER	Microphone used by operator to communicate with patient
SPK L (blue) & R (red)	External speaker for free field tests: right (R, red) and left (L, blue)
RECD/INSERT L & R	RECD transducers, right (R) and left (L) or Insert earphones (ER-3C, IP30) for Air Conduction: right (R) and left (L)
RECD BOX	External speaker for coupler measurements
CLIENT	Spare, for future uses
MON	Monitor headphones for operator
USB	USB cable for connection to personal computer
SPK (green)	Spare, for future uses
15V/2A ⊖-⊙-⊕	Power adapter

Connect the power adapter to the input and to the relative power cord, which should be plugged into an a.c. mains socket supplying the voltage indicated on the adapter label.



All connections must be made with the instrument switched off, or at least before plugging in the power cord and the USB cable.



Use only the power adapter supplied by Inventis s.r.l.



Connect the device only to medical equipment or to devices equipped with an isolating transformer.



Connect the device to an electrical power line carrying a voltage between 100 V and 240 V, with frequency between 50 Hz and 60 Hz.

3.1.2 Front Panel

A socket located at the bottom right hand corner of the Trumpet (the only connector on the front panel) is used to plug in the wired probes used for REM measurements.



Connector				F	Part			
PROBE	Set mea	of surer	probes nents	with	wires	used	for	REM

A Led positioned around the connector indicates the communication status between the Trumpet and the computer. Whenever the device is communicating with the PC controlling it, the Led will remain alight.

NOTE:

To perform the REM test, only the USB cable and power supply (and optionally the operator monitor headphones in the MON input) need to be connected to the rear panel, and the set of probes to the front panel (for wired probes).



3.1.3 Connection diagrams for REM mode



3.1.4 Connection diagrams for audiometry mode

AC/BC	
transducers	







3.1.5 Positioning the device

The Trumpet can be positioned:

- on a table,
- mounted on a stand,
- wall mounted using a standard VESA bracket (optional),
- wall-mounted using a dedicated bracket, which can also house the AUD transducers (optional)



When securing the Trumpet to a VESA wall mount, take care to use M4 screws that penetrate the threaded holes to a depth of no more than 6 mm.

CHAPTER 4

Maintenance



When cleaning operations of any kind are required, the instrument must first be disconnected and switched off by unplugging the power cord and the USB cable.



Internal inspection and servicing of the instrument must be entrusted entirely to technicians approved by Inventis S.r.l.



Transducers are manufactured utilizing ultra-fragile diaphragms that could be damaged in the event of impact. Handle with care during maintenance operations.

The Trumpet device do not require any special periodic maintenance other than calibration and normal cleaning, both of which are described in this chapter.

4.1 CALIBRATION AND EQUALIZATION (REM)



The calibration of the reference microphone must be validated by technicians approved by Inventis S.r.l. The operation should be carried out at least once every 12 months, and whenever the transducer is replaced.



Every time there is a change to the geometry of the surroundings in which the device is used, the operator must run the speaker equalization procedure (for further information, refer to "Maestro – Fitting and hearing aid test functionalities – User Manual").



Every time the probe tube is changed, consequently for each new patient, the probe microphone must be recalibrated (for further information, refer to "Maestro – Fitting and hearing aid test functionalities – User Manual").



If an RECD unit is being used, the calibration of the reference microphone must be validated by technicians approved by Inventis S.r.l. The operation should be carried out at least once every 12 months, and whenever the coupler is replaced.

4.2 PERIODIC CHECKS

To maximize the service life of the instrument, it is advisable that the following tests be carried out on a daily basis:



The tests must be conducted with the instrument positioned for normal use. For more information, refer to the Maestro software manuals.

- Before switching on the instrument, check that all the parts are properly connected and that the cables and/or connectors, as well as all other components, are neither broken nor damaged externally.
- Run the calibration procedure for the probe microphones with a new probe tube to make certain they are functioning correctly.
- Check subjectively that the outputs when using air conduction and bone conduction are the same for both channels and at all frequencies. Conduct this test at 10 or 15 dB HL, a level sufficient to verify that the stimulus signal has been presented. The person who carries out this test should have optimum hearing.
- Check each frequency at 60 dB HL in AC mode and at 40 dB HL in BC mode: there must be no distortion, noise or parasite signals.
- Check for correct operation of the patient response button.
- Check the voltage of the earphone headset and the bone vibrator
- Check communication with the patient.

4.3 CLEANING AND MAINTENANCE OF TRANSDUCERS



Do not use liquids or sprays to clean the Trumpet. When cleaning the device, use a soft cloth moistened with a mild detergent.

The headphone cushions and the bone vibrator are made of biocompatible material, but are not sterile. To avoid the possible spread of infection when

these items are to be used on a new patient, always wipe the surfaces with a hypoallergenic proprietary disinfectant, following the maker's instructions. For the DD45/TDH-39 cushions, it is recommended to use a denatured alcohol wipe or denatured alcohol with a microfiber cloth.

The earpieces of insert earphones are made of biocompatible material and must be used once only, then discarded in compliance with current waste disposal regulations.

Use a cloth moistened slightly with non-caustic medical grade disinfectant to clean the parts of the probe headset that come into contact with the patient (O-ring around the ear, back of probes against the face).



To avoid damaging the DD45/TDH39 headphones, do not push it against a flat straight surface as this can create vacuum and cause a damage to the transducer (suction cup effect).



Keep the device — and in particular the microphone components — well away from liquids of whatever nature.



Make certain that moisture cannot form inside the unit. Moisture forming internally of the unit can damage the instrument and expose the user or the patient to the risk of electric shock.



The earpieces of insert earphones are not sterile. The use of unsterilized earpieces can cause ear infections.



The probe tubes are disposable, and accordingly, must not be washed and reused but changed for each new patient.



No part of the Trumpet is suitable for sterilization in an autoclave, or for other thermal sterilization or disinfection procedures.

4.4 **REPLACEABLE PARTS**

The transducers and parts can be disconnected from the appliance. Should a fault develop in any one of these devices, the part must be switched off and isolated from the main unit, having first switched off the instrument and disconnected it from the mains power supply.



Only parts supplied by Inventis should be connected to the instrument.

4.5 REPAIRS AND TECHNICAL ASSISTANCE

Before contacting the service department, make certain that all the possible solutions in the "Troubleshooting" appendix C have been tried.

There are no requirements for the sterilization or disinfection of parts being returned to the manufacturer for repair and servicing.

Should the instrument need to be sent to the Inventis service department or returned to the dealer, it is important that the original packing should be used, and that all parts and transducers are enclosed.

Contact the dealer or the Inventis service department <u>before</u> shipping the product for servicing or repairs. In most instances it will be possible to remedy the problem without resorting to shipment.

4.6 ENABLE NEW LICENSES

When a new license key is purchased from Inventis, this can be used to enable new functionalities of the Trumpet system. The "Licenses Activation Tool" software can be requested from the Inventis Service department. On receiving the zip file, extract and run the RE1LicensesActivationTool.exe.

Be sure to have the license key available. Connect the Trumpet to the PC with the USB cable provided and plug the adapter into the mains power socket (for more information, consult the chapter 3 of the *Trumpet User Manual*). Close Maestro and proceed to run the "RE1LicensesActivationTool" software. The following window will appear:

(i) RE1 Licenses Activation To	ol Rev: 1.0.1.3535		_		×
File					
	Connected device S/N:	RE1RA16999999			
Active licenses on con	nected device				
✓ REM	AUD		QuickS	IN	
License Key]] -	Enable Lice	ense	

The licenses currently active on the device are shown in the central area of the window.

Enter the license key received from Inventis in the 5 fields (4 characters per field), exactly as written on the license document. Now click the "Enable License" button: if the license has been enabled, a message confirms the success of the operation and the relative check box in the central part of the screen will be ticked. If a message appears warning that the key is not valid, check that it has been entered correctly. Remember that for licenses to be enabled, the device must be connected to the computer by way of the USB cable provided.



If there is no doubt that the key has been entered correctly and the license still cannot be enabled, contact Inventis customer service

4.7 **FIRMWARE UPGRADE**

The procedure for updating the firmware of the device is described in the manual entitled *Maestro – Fitting and hearing aid test functionalities – User Manual*, and can be implemented using the "Firmware Upgrader" tool stored in the Maestro installation folder under *Service/Trumpet /FirmwareUpgrader*.

4.8 CALIBRATION OF MICROPHONES BY CODE

The procedure for calibrating the reference microphones of REM probes or the RECD coupler microphone by code is described in *Maestro – Fitting and hearing aid test functionalities – User Manual*, and can be carried out using the Maestro software.

APPENDIX A:

Technical Specifications

REM specifications

REM classification

Real Ear Measurement system, class 2a (MDD 93/42)

REM GENERAL SPECIFICATIONS		
Frequency range	125 Hz – 12 kHz	
Reference microphones measurement range	40 – 110 dB SPL	
Probe microphones measurement range	40 – 130 dB SPL	
Stimulus intensity at 1 m	50 dB SPL – 90 dB SPL, steps of 5 dB (+/- 3 dB)	
Equalization method	Modified Pressure Method with Concurrent Equalization (Stored Equalization for open fitting)	
Analysis	1/3 - 1/24 octave bands (selectable by the user)	
Probe tubes	Silicone: nominal diameter 1 mm, nominal length 80 mm	

SIGNALS AVAILABLE FOR REM	
Туре	
ICRA	•
ISTS	•
The rainbow Passage	•
Speech	•
Assorted ambient sounds	•
Freq. lowering-s	•

Freq. lowering-sh	•
Pink noise (PN)	•
White noise (WN).	•

REM MEASUREMENTS AVAILABLE		
REUR / REUG (single side and bilateral)	•	
REOR / REOG (single side and bilateral)	•	
REAR / REAG / REIG (single side and bilateral*)	•	
RECD	opt.	
MPO (single side and bilateral*)	•	
Advanced (single side and bilateral*)	•	
Live (single side and bilateral)	•	

*: bilateral available only for ITE measurements

AVAILABLE TARGET CURVE REQUIREMENTS	
NAL-NL1	•
NAL-NL2	•
DSL V5	•

OUTPUTS AND TRANSDUCERS AVAILABLE FOR REM		
On-site probe microphones (wired or wireless)	•	
Insert transducers for RECD	•	
Built-in speaker (free field)	•	
Coupler base (opt)	٠	

Audiometer specifications

Audiometer classification

Pure tone audiometer type 2 and speech audiometer class A (IEC 60645-1, ANSI S3.6)

SIGNALS AVAILABLE FOR AUDIOMETRY		
Type		
Tone	•	
Warble tone	•	
MIC input for speech audiometry (live speech)*	•	
Narrow band noise (NBN)	•	
White noise (WN).	•	
Speech noise (SN)	•	

*the microphone must be positioned in front of the operator's mouth.

SIGNAL CHARACTERISTICS		
Attenuator step	1, 2 and 5 dB	
	Continuous	
Presentation mode	Pulsed at a rate of 0.5 Hz, 1 Hz, 2 Hz, or operator-	
	settable rate	

OUTPUTS AND TRANSDUCERS AVAILABLE FOR AUDIOMETRY				
Output				
Air conduction (DD45 or TDH-39 headphones)	•			
Air conduction (ER-3C or IP30 insert earphones)	•			
Bone conduction (B-71 vibrator)	•			
Free field (built-in speaker, external passive speakers)	•			

AUDIOMETRY TESTS AVAILABLE				
Pure tone audiometry	•			
Autothreshold	•			
Speech audiometry	•			
Master Hearing Aid	•			
QuickSIN®	opt.			
Stenger	•			

	PURE TONES AND WARBLE									
AVAILABLE FREQUENCIES AND MAXIMUM										
	OUTPUTS									
Freq. (Hz)	AC TDH39 DD45 ¹ (dB HL)	AC ER-3C IP30 ¹ (dB HL)	BC ¹ (dB HL)	FF Int. ² (dB HL)	FF Ext. ² (dB HL)					
125	80	90	-	60	65					
250	100	105	45	75	80					
500	110	110	65	85	90					
750	115	115	70	85	90					
1.000	120	120	75	85	90					
1.500	120	120	80	85	90					
2.000	120	120	80	85	90					
3.000	120	120	75	85	90					
4.000	120	110	75	85	90					
6.000	105	100	55	85	90					
8.000	95	90	50	75	80					

SPEECH AUDIOMETRY MAXIMUM OUTPUTS							
AC AC AC ER-3C BC ¹ FF FF TDH-39 ¹ DD45 ¹ IP30 ¹ (dB HL) (dB HL							
100	100	100	60	70	80		

opt.

NBN, WN AND SN MASKING AVAILABLE FREQUENCIES AND MAXIMUM OUTPUTS							
Freq. (Hz)	AC TDH39 ¹ (dB EM)	AC DD45 ¹ (dB EM)	AC ER-3C IP30 ¹ (dB EM)	BC ¹ (dB FL)	FF Int. ² (dB HL)	FF Ext. ² (dB HL)	
125	60	60	70	-	55	60	
250	80	80	85	105	60	65	
500	95	95	95	110	75	80	
750	100	100	100	110	80	85	
1.000	105	105	105	110	80	85	
1.500	105	105	105	110	80	85	
2.000	105	105	105	110	80	85	
3.000	105	105	105	100	80	85	
4.000	105	105	100	100	80	85	
6.000	100	100	90	90	75	80	
8.000	90	90	80	80	65	70	
WN	90	90	80	45	70	75	
VV IN	dB EM	dB EM	dB EM	dB EM	dB HL	dB HL	
CN	90	90	80	45	70	75	
SIN	dB EM	dB EM	dB EM	dB EM	dB HL	dB HL	

⁽¹⁾ If the instrument is powered up only by way of the USB cable, the maximum limits will be reduced by 10 dB
 ⁽²⁾ With USB power supply only, the FF outputs are disabled

COMPATIBLE TRANSDUCERS						
Туре	Manufacturer	Model	Impedance			
Supra-aural headphones	Telephonics Corp.	TDH39	10 ohm (1 kHz)			
Supra-aural headphones	RadioEar Corp.	DD45	10 ohm (1 kHz)			
Insert earphones	Etymotic Research Inc.	ER-3C	10 ohm (1 kHz)			
Insert earphones	RadioEar Corp.	IP30	10 ohm (1 kHz)			
Bone vibrator	RadioEar Corp.	B71	10 ohm (1 kHz)			

AMBIENT NOISE ATTENUATION VALUES FOR TRANSDUCER						
Freq	TDH 39 / DD45 with MX41\AR or PN 51 cushions	ER-3C IP30				
[Hz]	[dB]	[dB]				
125	3	33.5				
250	5	34.5				
500	7	34.5				
750	-	-				
1000	15	35				
1500	-	-				
2000	26	33				
3000	-	-				
4000	32	39.5				
6000	-	_				
8000	24	43.5				

	FF	ISO 389-7 (ANSI S3.6)	[dB] re 20 µРа]	22.1	11.4	4.4	2.4	2.4	2.4	-1.3	-5.8	-5.4	4.3	12.6
D LEVELS	B71*	ISO 389-3 (ANSI S3.6)	[dB] re1µPa]	ı	67	58	48.5	42.5	36.5	31	30	35.5	40	40
NT THRESHOL	ER 3C	ISO 389-2 (ANSI S3.6)	[dB] re 20 μPa]	26	14	5.5	2	0	2	3	3.5	5.5	2	0
ICE EQUIVALE	DD45	ISO 389-1	[dB] re 20 μPa]	47.5	27	13	6.5	9	8	8	8	6	20.5	12
REFEREN	TDH 39	ISO 389-1 (ANSI S3.6)	[dB] re 20 µPa]	45	25.5	11.5	7.5 (8)	L	6.5	6	10	9.5	15.5	13
	Freq		[Hz]	125	250	500	750	1000	1500	2000	3000	4000	600	8000

(*) The calibration of the bone vibrator (B71) refers to positioning on the mastoid bone.

Reference standards: ISO 389-1, ISO 389-2, ISO 389-3, ISO 389-5, ISO 389-7, ISO 389-8 and ANSI S3.6

SPEECH AUDIOMETRY REFERENCE EQUIVALENT THRESHOLD VALUES IEC 60645-1							
TDH39	DD45	ER 3C	B71	Free field			
Coupler: ANSI S3 60318-	9.7 (NBS-9A) / IEC -3 (6cc)	IEC 60318-5	IEC 60318-6	Referred to 0 degrees incidence			
dB	dB	dB	dB	dB			
[re 20 µPa]	[re 20 µPa]	[re 20 µPa]	[re 1 µN]	[re 20 µPa]			
20.0	20.0	20.0	55.0	20.0			
Α	Applied to USB, MIC, SPEECH NOISE and WHITE NOISE signals						

SPEECH AUDIOMETRY REFERENCE EQUIVALENT THRESHOLD VALUES ANSI S3.6							
TDH39	DD45	ER 3C	B71	Free field			
Coupler: ANSI S3 60318	8.7 (NBS-9A) / IEC -3 (6cc)	IEC 60318-5	IEC 60318-6	Referred to 0 degrees incidence			
dB	dB dB		dB	dB			
[re 20 µPa]	[re 20 µPa]	[re 20 µPa]	[re 1 µN]	[re 20 µPa]			
19.5	18.5	12.5	55.0	14.5			
A	Applied to USB, MIC, SPEECH NOISE and WHITE NOISE signals						

STATIC CALIBRATION FORCE					
Transducer	Force				
Supra-aural headphones	$4.5~\mathrm{N}\pm0.5~\mathrm{N}$				
Bone vibrator	$5.4 \text{ N} \pm 0.5 \text{ N}$				

PATIENT - OPERATOR COMMUNICATION				
Talk-over by way of microphone connected to PC or to device	•			
Talk-back by way of microphone connected to PC or to device	•			
Patient response button	•			

Other specifications

CONTROL OF DEVICE

By way of Inventis Maestro software installed on PC.

•

POWER SUPPLY	
Main unit	15V, 2A d.c.
External power supply	Input 100-240Vac, 47-63Hz, 0.9-0.34A Output 15 Vdc, 2A max.

RADIO MODULES		
M	- 10 mW (Bluetooth module)	
Maximum emitted power	transmission module)	
	- 2402-2480 MHz (Bluetooth module)	
Operating frequency band	- 140 kHz (Wireless power transmission module)	

DURATION OF CALIBRATION	
Reference microphone	6-12 months
Transducers	12 months

WIRELESS REM PROBES		
Battery	Type designation: Varta CP 1654 A3, rechargeable Li-Ion cells, 3.7V, 120mAh Type number: 63165 Autonomy: minimum 3h of continuous use	
	Recharge time: max 2h	
Operation	Temperature: between +15°C and +35°C Relative humidity: between 30% and 90% (no condensation) Pressure: between 700 hPa and 1060 hPa	
Transport and storage	Temperature between 0°C and +45°C Relative humidity: 90% max, no condensation Pressure: between 500 hPa and 1060 hPa	

AMBIENT CONDITIONS		
	Temperature between +15°C and +35°C	
Operation	condensation)	
	Pressure: between 700 hPa and 1060 hPa	
	Temperature between -10°C and +50°C	
Transport and storage	Relative humidity: 90% max, no condensation	
	Pressure: between 500 hPa and 1060 hPa	
Warm-up time	1 minute	

	MECHANICAL
Dimensions (LxDxH)	15.5 x 10 x 24.5 cm / 6.1 x 3.9 x 9,7 in
Weight of instrument only	1.5 Kg / 3.5 lbs

BACK PANEL CONNECTORS		
Description	Type	Connector
Power supply	Input	DC socket 2.1/5.5 mm
AC transducers, L and R	Output	2 x mono audio jack, 1/4"
RECD insert transducers / L and R	Output	2 x mono audio jack, 1/4"
Bone vibrator	Output	Mono audio jack, 1/4"
Free field L and R (AUD)	Output	2 x Euro type spring clamp pluggable connector
SPK free field (SPARE, FOR FUTURE USES)	Output	Euro type spring clamp pluggable connector
Monitor earphones	Output	Stereo audio jack, 3.5 mm
Patient headphones	Output	Stereo audio jack, 3.5 mm
RECD external speaker	Input - Output	6 pin mini DIN
Patient response button	Input	Mono audio jack, 1/4"
External operator microphone (talk- forward)	Input	Mono audio jack, 3.5 mm
Patient microphone (talk-back)	Input	Mono audio jack, 3.5 mm
USB	Input - Output	USB type B
FRONT PANEL CONNECTORS		
REM probe	Input	8 pin DIN connector

HEARING SAFETY FEATURES		
	Alert situations: audible signal levels higher than 100 dB HL (IEC 60645-1)	
Audiometer	Safety measures:	
	 Activation of "higher dB" button required to raise intensity beyond 100 dB HL 	
	2) Warning on test screen	
	3) Activation of channel "Norm. ON" function inhibited	
	Alert situations: audible signal levels above UCL thresholds,	
	or above a customized value between 110 and 130 dB SPL,	
	registering in the auditory canal	
	Safety measures:	
DEM	1) Warning of level higher than UCL threshold and	
KEM	interruption of measurement (if threshold set at UCL)	
	2) Warning of level above customized threshold (between	
	110 and 130 dB SPL) in the auditory canal and	
	interruption of measurement (if threshold set at	
	customized value)	
	3) Warning if control is deactivated	

INPUTS: SPECIFICATIONS		
Input	Type of connector	Electrical specifications
Power supply	DC plug 2.5mm	Centre pin +15V, outer collar negative
Patient response button	Mono jack, 6.3mm	Switches 3V at logic input (switching current: 10mA) Impedance: 10KΩ
RECD speaker	6 pin DIN connector	-
Talk over	Mono jack, 3.5mm	Dynamic microphone, Electret or 200 Ω Impedance: 47K Ω
Talk back	Mono jack, 3.5mm	Frequency response: 100-12KHz +/- 3dB Electret bias: 2.2V with 2.2KΩ
REM probe	8 pin DIN connector	4 x omnidirectional microphone, Back Electret Condenser Impedance 2.2KΩ Frequency range: 20 – 16000 Hz Standard working voltage: 2 V Maximum working voltage: 10V

Current drawn 0.5 m A may		
Current draw: 0.3 IIIA IIIax		Current draw: 0.5 mA max

OUTPUTS: SPECIFICATIONS		
Output	Voltage available	Nominal impedance
TDH 39	16Vpp	10 Ω
DD45	16Vpp	10 Ω
ER 3C / RECD	16Vpp	10 Ω
B 71	16Vpp	10 Ω
Free field (SPKF),	16Vpp	4-8 Ω
Monitor	1Vrms wi	th 32 Ω load

SPECIFICATIONS OF SIGNALS AVAILABLE		
WN	Pure random noise, filtered Low cut-off frequency: 50 Hz High cut-off frequency: 20 kHz	
SN	As specified in IEC 60645-1	
PN	Pure random noise, filtered Low cut-off frequency: 50 Hz High cut-off frequency: 20 kHz	
Warble tone	Frequency of modulation signal: 10 Hz Modulation waveform: triangular Modulation coefficient: ±10%	
NBN	Band: ¹ / ₂ octave	
ISTS	International Speech Test Signal. See Holube, I., Fredelake, S., Vlaming, M. & Kollmeier, B. (2010). Development and analysis of an international speech test signal (ISTS). International Journal of Audiology, 49, 891-903 for further details.	
ICRA	International Collegium of Rehabilitative Audiology. See website "icra-audiology.org" for further details.	

APPLICABLE STANDARDS		
REM measurements	IEC 61669 / ANSI S3.46	
Audiometry	IEC 60645-1 / ANSI S3.6, Type 2 (pure tone), Class A (speech)	
Autothreshold	Amended by ISO 8253-1	
Calibration	<u>AC</u> : ISO 389-1 (TDH 39, DD 45), ISO 389-2 (ER 3C) <u>BC</u> : ISO 389-3 <u>FF</u> : ISO 389-7	
Electrical safety	EN 60601-1 Class I Type B	
EMC	EN 60601-1-2	
Radio	FCC CFR 47 Part 15 Subpart B and C and Part 1 Subpart I (§1.1310) ETSI standards for the compliance to the RED Directive (2014/53/EU)	

CE CERTIFICATE				
93/42 classification (DDM)	Class IIa			
Classification rules (Annex IX of 93/42)	10			
Notified body	TÜV SÜD Product Service GmbH Ridlerstrasse 65 D-80339 Műnchen			
Notified body number	0123			

Inventis will make available on request circuit diagrams, component part lists, descriptions, calibration instructions or other information that will help service personnel to repair those parts of the device that are designed by Inventis as repairable by service personnel.

APPENDIX B:

Specifications of speaker for free field tests

To produce the intensity levels in free field conditions specified in Appendix A, the speaker utilized must meet the specifications indicated below.

Speaker specifications:

Minimum input power	25 W
Nominal impedance	4-8 Ω
Sensitivity (at 1 m for 1 W, 1 kHz)	>93 dB SPL
Pass band	105 Hz – 20 kHz

APPENDIX C:

Troubleshooting

Problem	Possible cause	Solution
No signal from a transducer	Transducer not connected to correct output	Connect the transducer to the correct output
	Transducer damaged	Contact Inventis service department or dealer
No signal from patient response button when pressed	Wrong connection	Connect the patient response button to the correct socket
	Patient response button damaged	Contact Inventis service department or dealer
No signal from monitor earphones	Wrong connection	Check that the monitor headphones are plugged into the correct socket, according to the selected type of connection (instrument or PC)
	Headphone volume too low	Increase the audio volume by way of the Maestro control
Patient's voice cannot be heard	Problems with TALKBACK input connection	Check the connection to the TALKBACK input

	TALKBACK volume too low	Adjust the TALK BACK volume by way of the Maestro control	
Connection between PC and device cannot be established	USB port of the computer not working	Connect the instrument to a different USB port	
	USB cable damaged	Replace the USB cable (standard USB cable, type A/B)	
Results of test not credible	REM: Reference microphones not calibrated	Request calibration of reference microphones	
	REM: Equalization not valid	Run speaker equalization procedure	
	REM: Probe microphones not calibrated	Calibrate probe microphones	
	REM: Auditory canal obstructed by wax	Conduct an ear inspection and proceed to remove wax	
	AUD: Transducers not calibrated	Run transducer calibration procedure	
	AUD: Incorrect type of AC transducer selected (headphone or earphone)	Change the selected type of AC transducer in use by way of the Maestro control	
Speaker equalization procedure not successful	Probes not suitably positioned	Check that the probes are placed correctly on the patient, and that the patient is positioned facing the instrument, at a distance of around one metre	
	Background noise level too high	Move the instrument to a quieter room or take steps to reduce background noise as much as possible	
	REM reference microphone not functioning properly	If the problem persists after several attempts, contact customer service	
----------------------------------------------------------	--------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------	--
	Probe tube blocked or worn	Change the probe tube	
Probe tube calibration procedure not successful	Probe tube not inserted to the correct position for calibration purposes	See heading 2.2.4 <i>Probe tube calibration</i>	
	REM probes too far from speaker	Bring the REM probes nearer, to around 0.5 m	
	REM probes or reference microphone not functioning properly	If the problem persists after several attempts, contact customer service	
A test cannot be accessed	Optional test not activated	Contact the Inventis technical service department to obtain the license, supplying the serial number of the instrument	
FF (free field) output cannot be selected	Power adapter not connected	Connect the power adapter supplied with the instrument	
No sound from the required speaker in FF mode	FF output not correctly configured	Check that the correct configuration was selected during the calibration procedure	

Sound produced	Patient not positioned facing the device	Position the patient directly in front of the device, at a distance of around one metre	
distorted	Speaker faulty	Contact customer service	
	The probes are not paired correctly with the Trumpet	Perform the pairing procedure again (see Section 2.2.2)	
The wireless REM probes do not connect to the Trumpet (orange LED)	The Maestro REM module is not open or the set connection is wired	Access the Maestro REM module and set up the wireless connection with the probes (see the manual <i>Maestro – Hearing aid fitting</i> <i>and testing features – User</i> <i>Manual</i>)	
	The probes do not recognize the Trumpet	Turn off the probes by pressing the LED button for 10 seconds and disconnect the Trumpet from the USB socket. Reconnect the Trumpet and turn the probes on again by pressing the LED button. Perform the pairing procedure again if necessary.	
The wireless REM probes do not charge when in the charging	The probes are not centered correctly in the charging base	Remove the probes from the charging position, wait 5 seconds and then put them back in position, trying to center them correctly in the base	
position	The Trumpet is not connected to the power supply	Plug the Trumpet into the power supply	

	The Trumpet is not connected to a powered USB port (e.g. PC off)	Connect the Trumpet to a powered USB port or the USB port of an externally powered hub
The wireless REM probes do not turn on by pressing the LED button	The battery charge is too low	Connect the Trumpet to a powered USB socket and the power supply, then put the probes in the charging position and wait at least 20 minutes before removing them.



When using the instrument in conjunction with a soundproof booth, check that the connections both inside the booth and between the booth and the instrument are correct and secure.

APPENDIX D:

Electromagnetic emissions

The instrument has been tested and been found to comply with the limits for electromedical devices imposed by EN 60601-1-2 standard. These limits guarantee a reasonable level of protection against harmful interference in a typical medical installation.

The instrument generates, utilizes and radiates energy at radio frequency, and if not correctly installed and used as indicated in the instructions, can cause harmful interference to other devices located in the immediate vicinity. At all events, there can be no guarantee that interference will not occur in certain situations.



The Trumpet must not be used in close proximity to other equipment or stacked together with other equipment. If this cannot be avoided, check the operation of the instrument carefully to ensure that it is able to deliver normal levels of performance in the configuration adopted.

Should the instrument appear to interfere with other devices — this can be verified by switching off the instrument and then switching on again — it may be possible to remedy the problem by adopting one or more of the following solutions:

- change the orientation and/or the position of the device affected by interference;
- distance the devices one from another;
- consult the manufacturer or a service centre for assistance.

Cables, transducers and accessories

Inventis declares compliance with EN 60601-1-2 standard using the following cables, transducers and accessories:

- 1) Power cord (maximum length: 1.8 m)
- Medical grade power adapter, model n° MPU31-106-P01K002-B1F4-15V (L5)
- 3) REM probe (made by Inventis for the Trumpet)
- 4) TDH39 or DD-45 supra-aural headphones, 2 m shielded cable, dual signal (same as supplied with the device)

- 5) B-71 bone vibrator, 2 m braided cable, not shielded, dual signal (same as supplied with the device)
- 6) Patient headphones, manufacturer: Sennheiser, model: HD201
- 7) Monitor headphones with microphone, manufacturer: Sennheiser, model: PC3 (same as supplied with the device)
- 8) Clip-on talk-back microphone, 1 m shielded cable (same as supplied with the device)
- 9) Patient response button (made by Inventis), 2 m shielded cable
- 10) Insert earphones: ER-3C (manufacturer: Etymotic) (same as supplied with the device) or IP30 (manufacturer: RadioEar Corp.)
- 11) Speaker cable, not shielded, maximum length: 5 m
- 12) Shielded USB cable, maximum length: 2 m



The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by Inventis as replacement parts for internal components, can have the effect of increasing emissions from the device and reducing its immunity.

Anyone connecting other items of equipment must ensure that the overall system complies with EN 60601-1-2 standard.

Manufacturer's indications and declaration - electromagnetic emissions				
The Trumpet is designed for use in an electromagnetic environment as specified				
below. The user of the Trumpet must ensure that the instrument is operated in				
ambient conditions as specified.				
Emission tests Conformity		Electromagnetic environment - guidelines		
RF emissions CISPR11	Group 1	The Trumpet utilizes RF energy solely for internal functions of the instrument. Consequently, the RF emissions generated are minimal and unlikely to interfere with other equipment operating nearby.		
RF emissions CISPR11	Class B	The Trumpet is suitable for use in		
Harmonic emissions IEC 61000-3-2	Class A	any structure, including dwellings, and in environments connected		
Fluctuations in voltage/emissions (flicker) IEC 61000-3-3	Compliant	directly to the public low voltage grid supplying electrical power to residential buildings.		

Manufacturer's indications and declaration - electromagnetic immunity			
The Trumpet is designed for use in an electromagnetic environment as specified below. The user of the Trumpet must ensure that the instrument is operated in ambient conditions as specified.			
Immunity tests	Test level to EN 60601	Level of Compliance	Electromagnetic environment - guidelines
Electrostatic discharges (ESD) IEC 61000-4-2	± 6kV contact ± 8 kV air	± 6kV contact ± 8 kV air	The floor must be wood, concrete or ceramic tile. If floors are covered with materials that present synthetic coatings, relative humidity must be at least 30%.
Electrical fast transients / bursts IEC 61000-4-4	± 2 kV for electrical power supply lines ± 1 kV for input/output lines	± 2 kV for electrical power supply lines ± 1 kV for input/output lines	The quality of the electricity supply from the grid must be the same as for typical commercial or hospital premises.
Overcurrent IEC 61000-4-5	$\pm 1 \text{ kV from}$ line(s) to line(s) $\pm 2 \text{ kV from}$ line(s) to earth	$\pm 1 \text{ kV from}$ line(s) to line(s) $\pm 2 \text{ kV from}$ line(s) to earth	The quality of the electricity supply from the grid must be the same as for typical commercial or hospital premises.
Voltage dips, short interruptions and fluctuations affecting voltage on incoming electrical power supply lines IEC 61000-4-11		$< 5\% U_{T} (dip >95\% in U_{T}) for half cycle. 40% U_{T} (dip >60% in U_{T}) for 5 cycles. 70% U_{T} (dip >30% in U_{T}) for 25 cycles. <5% U_{T} (dip >95% in U_{T}) for 5 s.$	The quality of the electricity supply from the grid must be the same as for typical commercial or hospital premises. If the user of the Trumpet needs to rely on uninterrupted operation of the unit even in the absence of mains power, the instrument must be connected either to a UPS or to a battery.
Magnetic field at power frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields must correspond to the levels characteristic of sites located typically in standard commercial or hospital premises.

ŀ

Manufacturer's indications and declaration - electromagnetic immunity				
The Trumpet is designed for use in an electromagnetic environment as specified				
below. The user	below. The user of the Trumpet must ensure that the instrument is operated in			
ambient condition	ons as specified	d.		
Immunity tests	Test level to EN 60601	Level of Compliance	Electromagnetic environment - guidelines	
			In the event that there are portable and mobile RF communication appliances operating in the vicinity of the Trumpet, and the relative cables, be sure to maintain at least the recommended separation distance, calculated on the basis of the equation pertinent to the transmitter frequency. Recommended separation distance	
Conducted PE	3 Vrms			
	From150	3 Vrms	$d=1.2\sqrt{P}$	
IEC 61000-4-6	kHz to 80 Mhz			
	3 V/m			
Radiated RF	From 80	2 V/m	$d=1.2\sqrt{P}$ from 80 MHz to 800MHz	
IEC 61000-4-3	Mhz to 2.5 Ghz	3 V/m	$d=2.3\sqrt{P}$ from 800 MHz to 2.5 GHz	
			where <i>P</i> is the nominal maximum power output from the transmitter, expressed in watts (W) and indicated by the maker of the transmitter, and <i>d</i> is the recommended separation distance, expressed in metres (m).	
			The field strengths of fixed RF transmitters, established by way of an electromagnetic survey on site, <i>a</i>) must be lower than the level of compliance within each range of frequencies, <i>b</i>) may be associated with instances of interference near appliances bearing the symbol indicated below $(((\bullet)))$	
Note 1: At 80 MHz and at 800 MHz, the higher frequency range is applied.				
Note 2: These indications may not be valid in all situations. Electromagnetic propagation is				
affected by absorption and reflection on encountering structures, objects and persons.				

a) It is not possible to make a theoretical prediction with absolute certainty as to the field strengths of fixed transmitters, such as base radio stations (cell/cordless), terrestrial mobile telephones and radio, ham radio, AM and FM radio or TV transmissions. To evaluate the electromagnetic environment created by fixed RF transmitters, it may be necessary to

conduct an electromagnetic survey on site. If the field strength measured at the point where the Trumpet is in use exceeds the applicable RF level of compliance indicated above, check the operation of the instrument closely, to ensure that its performance levels meet the required standards. If abnormalities are observed, then additional measures may be necessary, such as changing the directional position or the location of the Trumpet.

b) In the frequency range from 150 kHz to 80 MHz, field strengths must be lower than 3 V/m.

Recommended separation distances between portable and mobile RF communication appliances and the device

The Trumpet is designed for use in electromagnetic environments where radiated RF interference phenomena are kept under control. The user of the Trumpet can help to prevent electromagnetic interference by ensuring that a minimum recommended distance is maintained between portable and mobile RF communication appliances (transmitters) and the Trumpet device, as indicated in the table below, according to the maximum power output from the appliances in question.

Nominal	Separation distance based on transmitter frequency (m)			
power output from	from 150 kHz to 80 MHz	from 80 MHz to 800 MHz	from 800 MHz to 2.5 GHz	
(W)	$d=1.2\sqrt{P}$	$d=1.2\sqrt{P}$	$d=2.3\sqrt{P}$	
0.01	0.12	0.12	0.23	
0.1	0.38	0.38	0.73	
1	1.2	1.2	2.3	
10	3.8	3.8	7.3	
100	12	12	23	

For transmitters with a nominal maximum power output not included in the table above, the recommended separation distance d in metres (m) can be estimated from an equation applicable to the transmitter frequency, where P is the nominal maximum power output from the transmitter in watts (W) specified by the manufacturer.

Note 1: At 80 MHz and at 800 MHz, the separation distance for the higher frequency range is applied.

Note 2: These indications may not be valid in all situations. Electromagnetic propagation is affected by absorption and reflection on encountering structures, objects and persons.

RF RADIATION EXPOSURE STATEMENT

This product complies with FCC and ISED radiation exposure limits set forth for an uncontrolled environment. The antenna should be installed and operated with minimum distance of 20 cm between the radiator and your body.

Cet appareil est conforme aux limites d'exposition aux rayonnements de l'ISED pour un environnement non contrôlé. L'antenne doit être installé de

façon à garder une distance minimale de 20 centimètres entre la source de rayonnements et votre corps.