

AUDIO LIBRARY HUB

PEHRATEK PRODUCTS



"Customize your Test playlist"

**REMOTE
CONTROL
OPERATION**



W-22



NU-6




SPONDEE




SPANISH

Combine all audio testing materials into one place -**AUDIO LIBRARY HUB** for quick easy access to specific test, title, or track for instant presentation, pausing or skipping to other stored files during testing. No more flipping CDs or searching tracks for the word list of choice as all files in the library are displayed with name, track, and elapsed time for easy identification. A built in bar graph displays track activity when playing. The wireless remote control makes for easy operation to select, skip tracks, or pause the presentation.




AUDIO LIBRARY HUB



Audiometer

- ◆ Desk top design offers small footprint and is easy to use and see
- ◆ Instant access to over 2000 audio files, test word lists, sound/noise files
- ◆ Display shows file name, title and track, elapsed time and bar graph
- ◆ Remote control over play, pause, file, right-left channel and volume
- ◆ Instant start up on selected file
- ◆ USB input, SD Memory card storage for program material
- ◆ FLAC processing is lossless with no compression to CD files

AUDIO LIBRARY HUB DISPLAY




“Easier to use than most PC based Audiometers”

FRONT PANEL



REMOTE CONTROL



800-470-3532