

PREMIUM NITRILE

AMMEX®

EXAM BLUE NITRILE GLOVES



Recommended Uses



General Description

Material	Nitrile	Type	Non-Sterile
Color	Blue	Usage	Single Use
Shape	Ambidextrous	Interior	Powder Free
Sizes	S - XL	Exterior	Textured/Micro-roughened
Cuff	Beaded	Process	Polymer Coated

Product Codes

Small	APFN42100
Medium	APFN44100
Large	APFN46100
X-Large	APFN48100

Packaging

Box	100 Gloves Per Box
Box Dims	219 x 111 x 60 mm
Case	10 Boxes Per Case
Case Dims	314 x 232 x 229 mm

Physical Properties

Tensile Strength	14 Mpa (min)
Tensile Strength (Aged)	14 Mpa (min)
Powder Content	< 2 mg
FDA Listing	510(k) Status
Elongation	500%
Elongation (Aged)	400%
Protein Content	N/A
AQL Level	1.5

Dimensions

Weight	4.2 ± 0.2 g	Length	240 ± 5 mm
Thickness		Width	
Cuff	3 ± 1 mil	Small	85 ± 5 mm
Palm	3 ± 1 mil	Medium	95 ± 5 mm
Finger	4 ± 1 mil	Large	106 ± 5 mm
		X-Large	116 ± 5 mm

Barcodes

 6 97383 10081 8 S Box	 6 97383 10082 5 M Box	 6 97383 10083 2 L Box	 6 97383 10084 9 XL Box
 6 97383 11081 7 S Case	 6 97383 11082 4 M Case	 6 97383 11083 1 L Case	 6 97383 11084 8 XL Case

