

Cerumen Microsuction With Confidence



- Gentle, safe suction
- Easy visualization
- Shorter and thinner suction tips allow fine control
- Thorough infection management

For the first time in the USA, we're offering our Cerumen Microsuction set-up featuring the Cliq[™] aspirator and Microsuction tips, to help doctors achieve excellent visualization of the ear canal.

Our unique Microsuction tips—with shorter working distance and narrow tubes—give every provider the confidence for quick and safe cerumen removal.


** Designed for use with Vorotek O Scope™ or a wall mount microscope

CIC ™ Vacuum Pump

- High-vacuum, low-flow suction
- Proprietary pump design for low-noise, highefficiency performance
- Noise reduction housing with collapsible handle
- 800 ml canister with overflow protection
- Microbial filter with intake protection
- Adjustable vacuum with display gauge







Cliq

78610 Cliq Portable Suction Pump

78615 Patient Tubing 6 ft each78616 Patient Tubing 12 ft each

78619 Canister + Lid Replacement

78622 Extension Tubing with Connector

89105 Replacement HEPA Filter 0.3 Micron



Microsuction Tips

78632 Ear Microscution Tips 14 ga **78637 Ear Microscution Tips** 17 ga

78624 Suction Controller

20/bx, Sterile, Single Use 20/bx, Sterile, Single Use 20/bx, Sterile, Disposable



Ear Loops

 78510
 Ear Loop
 2.0 mm

 78511
 Ear Loop
 2.5 mm

 78512
 Ear Loop
 3.0 mm

78513 Ear Loop 3.5 mm

78514 Ear Loop 4.0 mm

78515 Ear Loop 4.5 mm

78516 Small Hook78517 Large Hook90 degree hook90 degree hook

78518 Cotton Applicator



Ear Forceps

80108 Noyes Ear Forceps



