



Safe 'N' Clear®

CHANGING THE FACE OF HEALTHCARE



Industry:
Medical/Dental Consumables



Key Executive:
Dr. Anne McIntosh, President

SAFE'N'CLEAR & THE COMMUNICATOR™

A deaf-owned and woman-owned company whose mission is making sure medical and dental care are accessible and inclusive to patients and medical providers who may have special needs, such as pediatric patients, hearing loss, autism, learning disabilities, and geriatric patients. We are committed to doing our part to make communication more accessible and inclusive, which can lead to improved outcomes for patients and better profit margins for companies.

We distribute premium high-quality Made in USA ASTM Level 1 surgical facemask with a clear window that improves communication by reducing speech muffling and medical errors while boosting patient outcomes.

TECHNICAL SPEC SHEET

Description: FDA-approved F2100 ASTM Level 1 Surgical Face Mask with a Clear Window. Item #FM86000.

ASTM LEVEL 1

	ASTM F2100-11 LEVEL 1 STANDARD	THE COMMUNICATOR™ MASK
FLUID RESISTANCE WITH SYNTHETIC BLOOD	80 mm HG	80 mm HG
DELTA P DIFFERENTIAL PRESSURE	≤ 4.0 mm H ₂ O/cm ²	2.1 mm H ₂ O/cm ²
BACTERIAL FILTRATION EFFICIENCY (BFE)	≥ 95%	99.5%
SUBMICRON PARTICLE FILTRATION (PFE)	≥ 95%	99.33%
FLAMMABILITY	Class 1	Class 1