


5. EARIGATOR® Use and Cleaning

Before using, become familiar with the important parts of your EARIGATOR®.


5.1 Control Unit Parts

- A. Display panel
- B. Lid opening latch
- C. Hand piece cradle
- D. Hand piece holder
- E. Water reservoir
- F. Water level markings
- G. Thermistor
- H. Water level detector
- I. Drain port
- J. Thermal safety switch
- K. Heating element
- L. ON/OFF switch, fuse holder
- M. Power cord receptacle
- N. Water line connection
- O. Control signal connection
- P. Drain button



5.2 Hand Piece Parts (Patient Applied Part)

- A. Lens/Splash guard
- B. Nozzle
- C. LED bulbs
- D. Trigger
- E. Control signal cable connection
- F. Water tubing connection



5.3 Unit Markings



Patient applied part connection



Read operating manual before use