

# Wireless Patient Response Button

# **NEW! WIRELESS RESPONSE BUTTON**

- ✓ WIRELESS Bluetooth—Eliminates the cord!
- ✓ Easy-To-Use gentle touch activation button
- Provides LED light and vibration confirmation
- Easy-to-Clean—Wipe down after each use
- Rechargeable battery



- Small receiver gives selectable light and sound signal notification
- Can be plugged into audiometer or stand alone
- 3 Year Warranty



(ergo)

pehratek

www.pehratek.com

# **ERGO Wireless Patient Response Button**

## **WALL/TABLE-TOP WIRELESS RECEIVER**



- ✓ Wireless Bluetooth Response Button
- Easy for patient to grasp and use
- Durable—Easy-to-clean
- ✓ Wireless receiver gives light and sound indications
- Rechargeable operation—no batteries to change
- Operates weeks before recharging
- ✓ Hard wired output to plug into audiometer.

### **SYSTEM CONTENTS**

- Hand-held Response Button
- Wireless Receiver
- Power Adapter/Recharging Cord
- Audiometer Output Adapter Plug



# Gentle touch button — vibrates gently upon activation






Push button to select LED light color/sound combinations:

Hard-wired output jack to audiometer

# **RECEIVER (front)**



Desk top or wall mounting

- Light only
- Sound only
- Light & sound low
- Flashing light

### **Patient Benefits**

- Ergonomically shaped to grasp easily
- Provides light and vibrational feedback
- Rechargeable Lasts for weeks\*
- · Activated by gentle touch button
- Can be activated by finger, nose, toes
- Eliminates Perseveration response
- Easy-to-use for arthritic hands
- Extends use of portable audiometers

### **Examiner Benefit**

- High-Tech feel and look impresses patients
- More Accurate responses
- Sound and light indicators to flag attention
- Selectable color and sound signal combinations
- Can be positioned on a wall or desk top for best visibility
- Eliminates a cord
- · Low battery alarm
- · No switches to wear out

