



3-PLY SURGICAL FACE MASKS

COMFORTABLE LIGHTWEIGHT MATERIAL
MADE IN NEW YORK



Custom Packaging Available up to 20 Masks per pack

- ✓ Lightweight, non-woven polypropylene material
- ✓ Elastic loops fit comfortably around ears
- ✓ Nose clip molds mask securely to face
- ✓ Fluid resistant and disposable
- ✓ FDA Certified Facility
- ✓ Surgical Grade Mask
- ✓ 50 Per Box / 40 Boxes per Case / 20 Cases per Pallet
- ✓ SKU: 62212 (Level 1) / SKU: 62222 (Level 2)
- ✓ Made in NY, USA

ASTM F2100-11	Level 1 (62212)	Level 2 (62222)
Fluid Resistance, mmHg	80	120
BFE (Bacterial Filtration Efficiency)	>95%	>98%
PFE (Particulate Filtration Efficiency)	>95%	>98%
Differential Pressure, mm H ₂ O/cm ²	<5.0	<6.0
Flame Spread	Class 1	Class 1

*Tested and certified in an independent laboratory



Made with USA raw materials



100% made in the USA



Every mask gets inspected



Fast shipping



ALTOR SAFETY

855.378.5055 • altorsafety.com • 711 Executive Boulevard, Suite C, Valley Cottage NY 10989