


NEW PRODUCT BULLETIN

INTERCHANGEABLE TOY KIT(TOY-012)



The optional **INTERCHANGEABLE TOY KIT (TOY-012)** can be ordered for the **XM-1000** and **WR-500** and **WR-500 QUAD VRA** systems to permit the audiologist to swap out toy figures from time-to-time to present new and different toys to make the test interesting . The new Velcro pads allow for quick removal and replacement of each toy while firmly holding it in place throughout the test. Kit includes 2 Velcro mountable turn-styles and 6 alternative toy figures. (Of course, 2 toys would come installed in each VRA system shipped)



Each Toy figure includes a 2" Velcro mountable surface and quick change turn-style discs with mating Velcro pad .

Press-on Pull-off action

PEHRATEK PRODUCTS

PEHRATEK PRODUCTS
9285 Pierson Lake Drive
Chaska, MN 55318