Jodi®-vac Pro Models Compared











	Jodi [®] -XL with Air Pressure	Jodi [®] -XL	Jodi® Press/Vac	Jodi [®] Pro-Plus	Jodi [®] Pro
Oaktree Products Item #	JODIXL	JODIXL-NP	JODIPRESS	JODIPP	JODI
Dimensions (Includes Height of Tubing)	6 x 9 x 6 tall	6 x 9 x 6 tall	6.5 x 4.5 x 6 tall	6.5 x 4.5 x 6 tall	3 x 4 x 2.5 tall
My Needs Are Best Described As:	I deal with moisture issues I would like to clean vents with air pressure I want to remove wax from sound ports I may want to remove wax using liquids	I deal with moisture issues I want to remove wax from sound ports I have air pressure and/or clean vents with other tools I may want to remove wax using liquids	I would like to dean vents with air pressure I want to remove wax from sound ports I may want to remove wax using liquids	I want to remove wax from sound ports utilizing liquids I have air pressure and/or clean vents with other tools	I do service calls outside my office Finding electrical outlets are an issue
Key Features See Product Specs for Detailed Information About Each Model	Silent Running Desiccator Flexible Reinforced Air Pressure Needles to Clean Vents and BTE Tubing Highly Efficient Vacuum Needles Liquid Filter Option (you can add this at any time in your office)	Silent Running Desiccators Highly Efficient Vacuum Needles Liquid Filter Option (you can add this at any time in your office)	Flexible Reinforced Air Pressure Needles to Clean Vents and BTE Tubing Highly Efficient Vacuum Needles Liquid Filter Option (you can add this at any time in your office)	Highly Efficient Vacuum Needles Liquid Filter	Highly Efficient Vacuum Needles Optional Battery Pack Carry Bag Optional
Vacuum	22.5 inches of Mercury	22.5 inches of Mercury	22.5 inches of Mercury	22.5 inches of Mercury	22.5 inches of Mercury
Custom Vacuum Needles	Yes	Yes	Yes	Yes	Yes
Pressure	>38 PSI	-	>38 PSI	_	_
Custom Air Pressure Needles	Yes	Yes	Yes	_	_
Heated Desiccator with Silent Operation	Yes	Yes	_	_	_
Liquid Filter Option	Yes	Yes	Yes	Yes	_
Straight or Coiled Tubing (Field Replaceable)	Yes	Yes	Yes	Yes	_
Battery Pack	_	_	_	_	Yes
Carry Bag	_	_	_	_	Yes
2000-Hour, Serviceable Motor/Pump	Yes	Yes	Yes	Yes	Yes
Warranty	3 Years	3 Years	3 Years	3 Years	3 Years