Product Overview

The Mark V+ continues the legacy of the original Mark V. The Mark V+ offers the most powerful suction and quietest operation on the market.

With the Mark V+ you can suction wax from hearing aids and earmolds and blow air to remove dust and debris from devices and surfaces.

The large powerful drying chamber provides a powerful vacuum for completely drying out instruments quickly.

Warranty

The Mark V+ comes with a one year warranty on all parts and labor. The warranty does not cover malicious damage.

Technical Specifications

The Mark V+ is equipped with a 100-240V, 50/60Hz power supply so it can be used anywhere in the world with the appropriate plug adapter. The circuit breaker is 4 Amp.

Environment:

- a) Environment Temperature: 40°F 104°F
- b) Air Humidity: ≤80%RH (when temperature is 77°F)
- c) Atmospheric Pressure: 86kpa~106kpa
- d) Do not use in the presence of flammable or corrosive gases or conductive dust

Maximum Negative Pressure: ≤ -28 inHg

Vacuum Drying Chamber: 3.75"L x 3.75"H

Continuous Run Time: 15 minutes

Noise: ≤ 55 dB

External Dimension: 13.38"L x 7.5"W x 11.5"H

Weight: 9.64 lbs