

POCKETALKER ULTRA













(HED021) PKT-D1-EH

(EAR013) PKT-D1-EH

(EAR008) PKT-D1-E08

(EAR014) PKT-D1-E14

(NKL001) PKT-D1-N01

ITEM #: PKT-D1-XXX

DESCRIPTION: Non-Wireless Personal Listening System

WARRANTY: 5 Year

INCLUDES: Amplifier, Plug Mount Microphone, 12 ft TV Listening Cord, Belt Clip Case, 2 AAA

Batteries, 1 of 4 Earphone/Headphone Listening Options

WHAT IS THIS?

Williams Sound updated version of a non-wireless personal listening system with 4 different listening options for users.

HOW DOES IT WORK?

- Works using 2 AAA batteries (included)
- Plug 1 of 4 Listening Options (Headphone, Single Earbud, Dual Earbud, Wide Range Earphone, Neckloop) into Earphone/Headphone Jack
- Turn on unit using volume wheel, "ON" indicator on unit will illuminate
- Adjust volume using wheel
- Adjust tone using external tone control wheel and begin listening

WHAT IS SO COOL ABOUT THIS PRODUCT?

- Lightweight and easy to use
- New modern design with a sleek look
- External volume and tone control for user to adjust according to their personal preference
- Available in Neckloop Model for "T" coil equipped hearing aid users #PKT-D1-N01

ARE THERE ANY DOWNSIDES?

- Tone control may cause user confusion when adjusting
- AAA batteries may be difficult for older users or someone with dexterity problems