



PROFESSIONAL-GRADE
DRYING. **NOW AT HOME.**



A **GAME-CHANGER** FOR YOUR HEARING DEVICES.

Redux Home Dryer removes not just some, but **all** traces of harmful moisture from hearing devices so they can sound their best.

✓
COMPLETE MOISTURE REMOVAL

✓
SAFE FOR RECHARGEABLES


✓
NO DESICCANT REFILLS

HEAR THE DIFFERENCE.

Ever feel like your hearing devices lose effectiveness over time? With just one drying session, Redux improves clarity and increases volume capabilities back to peak functionality.

GET REDUX HOME DRYER TODAY.
VISIT WWW.REDUX.COM

REMOVING “A LITTLE MOISTURE” CAN MAKE A BIG DIFFERENCE.



DEVICES THOUGHT TO BE DEAD: **50%** REVIVED

DEVICES THAT SOUND WEAK: **82%** IMPROVED

DEVICES THAT SOUND GOOD: **80%** IMPROVED

89%

report subjective improvement in their devices

67%

of dead receivers (speakers) are recovered after Redux

24%

of significantly improved devices show improvement to microphone directionality¹

12%

of brand-new hearing aids significantly improve right out of the box

Reported findings from a top 10 audiology doctoral program study.

¹This is the first study that has looked at and identified moisture impact on microphones and directionality. Directionality determines how well the person will hear in noisy environments.

SIGNS YOUR HEARING DEVICES NEED MOISTURE REMOVAL:

- Not hearing as well as day one
- Sound is less clear or softer; it's hard to hear in noisy environments
- Battery issues, intermittency, or connectivity issues

BENEFITS OF REDUX HOME DRYER:

- Improved sound quality: louder, crisper, clearer hearing
- More reliable: less repairs, longer battery life, consistent performance
- Cost effective: proactive drying keeps hearing aids working like new for longer

THE GOLD STANDARD OF DRYING




	REDUX HOME	OTHER DRYERS
Drying Cycle Time	45 minutes (or until dry)	45 minutes (fixed timer)
Temperature Used to Dry Devices	No heat required	>100°F
H ₂ O Sensors for Endpoint Drying	✓	✗
Guaranteed 100% Moisture Removal	✓	✗
Recharges Desiccant	✓	✗
Bluetooth Connectivity	✓	✗
Moisture Removed Reporting	✓	✗
Machine Learning to Recommend Drying Frequency	✓	✗


HOW IT WORKS

1. Redux creates ultra-dry air that is “thirstier” and holds more moisture than ambient air. Using a regenerative desiccant, the environment in the Redux dryer attracts moisture out of the hearing device.
2. Once the surrounding air is saturated with moisture from the device, Redux cycles it through the desiccant again until the environment is returned to ultra-dry.
3. Steps 1 and 2 repeat until the humidity sensors detect no more moisture inside the hearing device.
4. Once dry, the desiccant is regenerated with heat.
5. Download the Redux app to see the amount of moisture removed from your hearing device.

THE TECHNOLOGY BEHIND REDUX HOME DRYER

Redux Home Dryer is a complete moisture removal system, backed by over 70 patents and trusted by thousands of hearing care professionals worldwide.





“

I'm having great success with my home dryer. My hearing aids sound much better. Such a knowledgeable staff to work with. Thank you.

DAVID H.

Our Redux Home dryer works great! This is the same brand dryer that our audiologist has in her office. We have tried all different types of home hearing aid dryers for our active 8-year-old, and this is hands down the best we've ever had. It takes about 2 hours to complete the drying cycle, so we use it overnight. It gets all the moisture out! Highly recommend! Amazing customer service!

RYAN B.

Unlike other options on the market, the Redux Home Dryer doesn't require desiccant refills, additional equipment, fees, or subscriptions. It has made a significant difference in the maintenance of my devices, and I couldn't be more satisfied with my purchase.

JASON

The user-friendly design of the dryer is also worthy of praise. It's simple to use and has clear instructions, making the entire process hassle-free. The compact size means it doesn't take up much space. My grandfather can easily take it when he travels.

ANGELICA

”

OUR PARTNERS

Redux is proud to work with hearing care providers across North America.



GIVING BACK TO OUR COMMUNITY.

HEARING THE CALL

We work with Hearing the Call, a nonprofit that provides access to hearing health care in developing regions. The Redux team has attended several mission trips in Mexico, fitting over 500 patients with hearing aids. Through Hearing the Call, we also donate Redux Pro systems to low-income areas in humid climates around the world like Guatemala, Mozambique, and Tasmania.

GIFT OF WATER

We support Gift of Water, a nonprofit providing filtered, clean, and available drinking water to developing countries. We have donated \$10,000 to their cause, helping 200 families in Haiti.



GET REDUX HOME DRYER TODAY.

Try it risk-free with our 30-day money-back guarantee.

www.redux.com



- @Get-Redux
- @GetRedux
- @GetRedux

Questions? Contact us at support@redux.com.