

SCREENING PORTABLE TOUCH SCREEN AUDIOMETER

PRODUCT SPECIFICATIONS

GENERAL SPECIFICATIONS

DIMENSIONS AND WEIGHT

- L x W x H: 180 x 135 x 41 mm
- Net weight: 450 g.
- Body material: self-extinguishing ABS

TEST TYPES

Pure Tone test, Speech Test (performed with on-board recorded speech material)

DISPLAY

- 7" TFT Resistive touch screen display
- 800 by 480 resolution at 133 pixels per inch (ppi) 262K colors

USER INTERFACE

- Intuitive icon system
- Internal patient database

REPORTS

- Data transfer to PC with Resonance Management Data
 Suite installed
- Dedicated report template

DATA TRANSFER TO PC

- Bluetooth 2.1 • USB cable
- · ODD cable

USB PORT

- Nr.1 USB host type A
 (pen drive connection for software upgrade)
- Nr.1 mini USB slave (for patient response connection)
- WINDOWS® COMPATIBLE SOFTWARE
- Resonance MDS Management Data Suite

MAIN SYSTEM COMPONENTS

- Operative system: Microsoft[®] Windows[®] Embedded Compact
- CPU: TEGRA® 2 dual core @ 1Ghz

POWER

POWER SUPPLY

- Built-in 3.6V/3.1 Ah rechargeable lithium-ion battery
- External Medical power supply/charger: Input 100 - 240 Vac / 0.3 A - 50/60 Hz; Output +5 Vdc / 2 A

CONSUMPTION (WITH RECHARGEABLE BATTERY)

- Up to 3 hours of testing working on speech test modality.
- Charge time with external medical power supply: up to 80% in 2 hours

ENVIRONMENTAL

OPERATING ENVIRONMENT

- Storage: -20° C up to +50° C
- Operating: +15° C up to +35° C
- Humidity: up to 90%, (non-condensing)
- Ambient pressure: from 700 hPa up to 1060 hPa

Start

r17a

AUDIOMETRY OPERATING SPECIFICATIONS

RANGE

- Frequency range (AC):
- 125 8000 Hz • Range stimuli level -10 up to 110 dB HL

ACCURACY

- Frequency < 0.5%
- Distortion < 1%
- Crosstalk >70dB
- Attenuator linearity 1 dB per 5 dB step, max 3 dB whole range

TYPE OF SIGNALS

- Pure tone: sine wave 125 to 8 Khz signal
 Warble: +/- 5% frequency sine wave modulated,
- modulation: sine wave 5 Hz • Narrow band noise: 24 dB/oct filtered noise
- Speech noise: 1 khz 12 dB/oct filtered noise
- White noise
- Speech material recorded on internal SD-card
 On/Off rise fall time: 40 msec

OUTPUT TRANSDUCERS

• ACR, ACL: 10 ohm DD45 matched pair earphone. IP30 (optional)

STIMULUS PRESENTATION MODALITY

- Presentation: Normal
- Modality: Steady, Pulsed (rate 2 Hz)

• dB step variation: 1 dB, 2 dB, 5 dB.

QUALITY SYSTEM

R

Manufactured, designed, developed and marketed under an ISO 13485, ISO 9001 certified quality system. Medical CE marks and FDA approval.

0

resonance

COMPLIANCE/REGULATORY STANDARDS

Designed, tested and manufactured to meet the European and International Standards:

- MDD 93/42/EEC and its revised versions: Class IIa (as referred to in Annex IX, rule 10 of said MDD 93/42 EEC) • Safety: IEC 60601-1, 3rd edition, Class II Type BF
- EMC: IEC 60601-1-2
- Audiometer:
- IEC 60645-1and ANSI S3.6: Type 4 (R17A-AC) IEC 60645-2 (only for the functions regarding the speech test; the device is not an A or B model)
- RoHS II compliant

STANDARD ACCESSORIES

- DD45 headphones
- Built in talk over microphone for patient communication
- Multilingual Quick user's handbook
- Resonance[®] MDS software
- Medical Power supply (100 Vac 240 Vac)
- Carrying bag
- Charging Cradle
- Pen Drive (with Multilingual Quick user's handbook)
- Touch screen pen

OPTIONAL

- External USB Keyboard
- IP30 Insert headphones
- Pediatric Headset
- RE-Play Audiometry for children
- Wallmount kit
- RE-VO USB Video otoscope
- Wi-Fi thermal printer
- Autothreshold and patient response push button
- Amplifier and free field speakers



SKU: 4300104000

The information, pictures and specifications found in this datasheet are intended as a general guideline for customers seeking information about supplies for Resonance equipment. Resonance makes no warranty, nor assumes any legal liability or responsibility for the accuracy, typing errors or mistakes, correctness or completeness of any information in this datasheet. The information in this datasheet was correct to the best of our knowledge at the time of printing.