



CLINICAL MULTIFREQUENCY MIDDLE EAR ANALYZER

GENERAL SPECIFICATIONS

Tympanometry, Acoustic Reflex, Reflex Decay, nr. 3 Quick Tests, Acoustic Reflex Latency Test (ARLT), ETF

(Intact, Perforated and Patulous), Special tests

• External fast thermal printer available as optional

• .pdf report created directly from the device and stored

on USB Pen drive with possibility to add patient data

• Data transfer to PC using Resonance Management

· To help keep the child involved while running

a test, a series of animated images appears on the

WINDOWS® COMPATIBLE SOFTWARE

Resonance MDS Management Data Suite

and tests comments via the USB Keyboard (optional)

Non-acoustic),

DIMENSIONS AND WEIGHT

• L x W x H: 370 x 290 x180 mm • Net weight: 3.5 kg

TEST TYPES

(Growth-DLI and

• 7" TFT Color display

USER INTERFACE

"CHILDREN" FEATURE

DATA TRANSFER TO PC

Via cable through USB port

• Nr.1 USB host type A

• Nr.1 USB slave type B

POWFR 1

POWER SUPPLY

CONSUMPTION

• Max current: 0.15 A

• Fuses: 2 x T 1 A L 250 V

COMMUNICATION PORT

Tympanometry

DISPLAY

Multilingual

PRINTER

color display



PRODUCT SPECIFICATIONS

TYMPANOMETRY OPERATING **SPECIFICATIONS**

PROBE TONE

- 226 Hz for Admittance (Y) curve tympanometry
- 1000 Hz for Admittance (Y) curve tympanometry with added Susceptance (B) and Conductance (G) curves

INTENSITY

- 226 Hz: 85 dB SPL ± 2dB
- 678, 800 and 1000 Hz: 75 dB SPL ± 2 dB
- Frequency Accuracy: ± 0.5%
- Harmonic distortion: Less than 1%

ADMITTANCE MEASUREMENTS

- Compliance range at 226 Hz: 0.05 up to 7 ml
- · Sensitivity scale: Autoscale to appropriate range; available scale at 226 Hz: 1.5, 2, 5 or 7 ml
- \bullet Compliance range at 678, 800 and 1000 Hz: 0 to +25
- Sensitivity scale HF: Autoscale to appropriate range, available scales at 678, 800 and 1000 Hz: 5, 10, 15, 20, 25 mmho

AIR PRESSURE

- Control: Automatic and Manual
- Range: from +400 up to -600 daPa adjustable in 50 daPa steps
- Pressure accuracy: ±10 daPa or ±10%
- Sweep rate: 50, 100, 200, 300 daPa/sec and automatic
- Indicator: Measured value is displayed
- Safety limitations: -800 up to +600 daPa

EUSTACHIAN TUBE FUNCTION

- ETF test for use with both intact and perforated eardrums
- · Available also test for patulous tympanic membrane

ACOUSTIC REFLEX TESTS

- Reflex test method: Threshold, Automatic, Manual
- Stimulus duration: 0.5, 1 or 2 sec.
- User selectable protocols for all test methods
- Ipsi or Contralateral stimulation for all reflex test
- · Automated peak search functions available for all
- Manual Reflex: Pump manual control of all stimuli
- Reflex Decay: Threshold, Automatic, Manual
- Stimulus duration: 10 or 20 sec.
- · ARLT: Threshold, Automatic, Manual
- Stimulus duration: 1 sec. fixed
- Non-Acoustic: 10/20 sec.

FREOUENCIES AND INTENSITY RANGES

Ipsilateral:

- Pure tone level range (dB HL) from 50 to 110
- Noise level range (dB SPL) from 50 to 100
- Frequency: 500, 1000, 2000, 3000, 4000 Hz
- Noise: BBN, HP or LP
- Frequency accuracy: ±1%
- Harmonic distortion (THD): less than 3%

Contralateral: Pure tone level range (dB HL) from 50 to 120

• Noise level range (dB SPL) from 50 to 115

SKU: 4300102830

- Frequency: 250, 500, 1000, 2000, 3000, 4000, 6000, 8000 Hz
- Noise: BBN, HP or LP
- Frequency accuracy: ±1%
- Calibration accuracy ±3 dB
- Level steps: 1, 2, 5 or 10 dB • On/Off ratio: 70 dB minimum

QUALITY SYSTEM

Manufactured, designed, developed and marketed under an ISO 13485, ISO 9001 certified quality system. Medical CE mark and FDA approval.

COMPLIANCE/REGULATORY

Designed, tested and manufactured to meet the European and International Standards:

- MDD 93/42/EEC and its revised versions: Class IIa (as referred to in Annex IX, rule 10 of said MDD 93/42 EEC)
- Safety: IEC 60601-1, 3rd edition, Class 1 Type B
- EMC: IEC 60601-1-2
- Impedance: IEC 60645-5/ANSI S.3.39 Type 1

STANDARD ACCESSORIES

- Probe HF Multifrequency
- · Headband and probe handle
- · Probe cleaning tool kit
- DD45 for Contra
- Insert for Contra
- Kit of assorted ear tips
- Calibration cavity with probe holder • Device dust cover
- Multilingual Quick user's handbook
- · Pen Drive
- USB cable device to PC
- Resonance® MDS software with NOAH® module included
- Power supply cable (110 220 V)
- Spare fuse

OPTIONALS

- External USB Keyboard
- · External fast thermal printer
- · Carrying bag

Resonance makes no warranty, nor assumes any legal liability or responsibility for the accuracy, typing errors or mistakes, correctness or completeness of any information in this datasheet. The information in this datasheet was correct to the best of our knowledge at the time of printing.

FNVIRONMENTAL OPERATING ENVIRONMENT

• 110 - 240 V AC 50/60 Hz 40 VA

• Storage: -20° C up to +50° C

• Power consumption: 40 VA

- Operating: +15° C up to +35° C
- Ambient pressure: from 700 hPa up to 1060 hPa